WELLSAWIPES

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REFI WCW407

REFI WCW409

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REFI WCW409

Product description

WELLSAWIPES GETRÄNKTE DESINFEKTIONSTÜCHER are ready to use alcohol-soaked disinfection wipes in a refill bag. They are applied with a dispenser bucket which is available separately. The wipes are soaked with WELLSACLEAN S+ SCHNELLDESINFEKTION (containing alcohol and quarternary ammonium compounds) and can be used directly after insertion into the dispenser bucket. WELLSAWIPES GETRÄNKTE DESINFEKTIONSTÜCHER provide an universal and time-saving application.

Areas of application

For rapid disinfection of alcohol-resistant surfaces of non-invasive, non-critical medical devices and medical equipment.

Suitable for:

- Surfaces close to patients, dental chairs, examination couches, operating tables
- Areas in medical and dental practices which are used for medical purposes
- Rapid disinfection between surgeries and treatments
- Clothes like OR shoes, OR aprons, X-ray aprons
- Alcohol-resistant surfaces of medical equipment

Not tested for food and veterinary areas. Not suitable for disinfection of skin and hands.

Not suitable for disinfection of invasive medical devices. Do not use for disinfection of semi-critical and critical medical devices.

Special Advice

Do not reuse the product. If the disinfectant wipe is reused, the indicated effectiveness is no longer ensured.

Do not use the product if it is damaged or soiled.

Application only by qualified personnel.

Report serious incidents related to this product to the manufacturer and the responsible authority.

Application

Put on protective gloves, face masks, and eye protection prior to disinfection! Provide adequate ventilation.

1. In case of contaminated surfaces: Remove visible contaminations with a disinfection wipe and discard it.
2. Disinfect surface with a fresh disinfection wipe. Wipe the surface applying mild pressure (wet wiping). Keep the surface completely moist during contact time.
3. After contact time has elapsed: Remove residues with a clean, dry, disposable cloth or allow to dry.

Shelf life after opening

Shelf life of opened disinfection wipes in dispenser bucket: 28 days from day of opening if the lid is closed properly and tightly every time a wipe is taken out.

The indicated shelf life after opening can only be guaranteed when using the dispenser bucket tested by the manufacturer.

Summarized efficacy spectrum

Efficacy under dirty conditions (product ready to use, at room temperature)	Contact time
Bactericidal/Yeasticidal (VAH standard methods)	5 min
Bactericidal/Yeasticidal (EN13727, EN 13624, EN 16615)	5 min
Mycobactericidal/TB (EN 14348)	
Limited spectrum of virucidal activity (EN 14476, EN 16777)	30 sec
Virucidal activity against enveloped viruses (EN 14476)	

See page 2 for a detailed efficacy spectrum, material compatibility, and advice on filling of the dispenser bucket.

WELLSAWIPES

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REFI WCW407

REFI WCW409

Detailed spectrum of efficacy

Efficacy (product ready to use, at room temperature)	Standard	Contact time
BACTERIA		
S. aureus, E. hirae, P. aeruginosa, E. coli, MRSA → dirty conditions	EN 13727	30 sec
S. aureus, E. hirae, P. aeruginosa → dirty conditions	EN 16615	5 min
S. aureus, E. hirae, P. aeruginosa → dirty conditions	Standard methods VAH*	5 min
YEASTS		
C. albicans → dirty conditions	EN 13624	30 sec
C. albicans → dirty conditions	EN 16615	5 min
C. albicans → dirty conditions	Standard methods VAH*	5 min
MYCOBACTERIA/TB		
M. terrae, M. avium → dirty conditions	EN 14348	1 min
VIRUSES		
Enveloped viruses		
MVA Modified Vacciniavirus Ankara → dirty conditions	EN 14476	30 sec
BVDV → dirty conditions		
Non-enveloped viruses		
Adenovirus type 5 → dirty conditions	EN 14476	5 min
MNV murine Norovirus → dirty conditions		
Human Rotavirus → clean conditions		
MNV murine Norovirus → dirty conditions	EN 16777	
Adenovirus type 5 → dirty conditions		

*standard methods of the VAH (Association of Applied Hygiene) for surface disinfection by wiping (quantitative carrier test with mechanical action, practical conditions).

Material compatibility

Metals	aluminium, stainless steel, copper, brass
Plastics	PSU, PVC, silicone
Rubber	EPDM, latex
Advice	Upon compliance to the intended use and the indicated contact time, no material damage is to be expected with the above mentioned materials. Avoid contact of longer than 48 hours with linoleum flooring. Avoid contact with acrylates. May cause stress corrosion on acrylic glass. Test material compatibility for polyacrylates, polycarbonates as well as not mentioned materials before use.

Filling of the dispenser bucket

Comply with the dispenser bucket manufacturer's recommendations for reprocessing as well as national requirements concerning the reprocessing of dispenser buckets. Put on protective gloves, face mask, and eye protection. Provide adequate ventilation.

1. Record product data (provided on the refill bag) on the additional label and attach it to the dispenser bucket.
2. Put the refill bag into the prepared dispenser bucket.
3. Cut open the refill bag at the top, insert the wipe through the opening of the lid and close the bucket thoroughly. The product is now ready for use.

Advice

Avoid contact with other disinfectants. Do not mix with other disinfectants. Do not dilute. Apply on smooth surfaces only. On coarse surfaces lint formation cannot be excluded. Do not dispose of WELLSAWIPES GETRÄNKTE DESINFEKTIONSTÜCHER down the toilet.

Storage advice

Store the product between 5 and 25 °C. Protect from heat and sunlight. Store in an upright position and tightly closed in a well-ventilated place. Keep cool. Keep away from sources of ignition – No smoking. Take precautionary measures against static discharge. Keep away from food, drink, and animal feedingstuffs. Keep locked up and out of the reach of children.

Further information on cleaning and disinfection of surfaces and medical devices:

Recommendation of the Commission for Hospital Hygiene and Infection Prevention at the Robert-Koch Institute (RKI) entitled „Hygiene Requirements for Cleaning and Disinfecting Surfaces“. Available online.

Composition

100 g of soaking solution contain:

33.66 g ethanol, 15.84 g 2-propanol, 0.05 g didecyldimethylammonium chloride 70 %, perfume*2 (not applicable for perfume-free)

Chemical/physical parameters

appearance: clear and colourless
 density (20 °C): 0.900 — 0.920
 flash point (EN 22719): 22.75 + 0.5 °C

*2 Please note that perfume may cause allergies.

WELLSAWIPES

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REF WCW407

REF WCW409

Expert opinions & further information

WELLSAWIPES GETRÄNKTE DESINFEKTIONSTÜCHER have been assessed by independent laboratories. The shelf life after opening of 28 days was confirmed in an expert opinion.

The listings, reports, and notes in this product information refer to the product PLIWA Quick & Clean S-Wipes, which is identical in formulation.

To exclude adverse effects on the skin upon unintentional contact with the product, a dermatologically survey has been carried out with the soaking solution by the company Dermatatest.

For further details concerning expert opinions, information on products, and answers to specific questions contact us by e-mail (info@wellsamed.com).

The manufacturer has a manufacturing authorisation according as per §13 German medicines act as well as a GMP- and an EN 13485 certificate.

Labelling according to Regulation (EC) No 1272/2008 (CLP)

For information on labelling according to Regulation (EC) No 1272/2008 (CLP) see safety datasheet and label.

Identification, certificates

CE 0123 MD IIb

The PLIWA Quick & Clean S-Wipes, which has the same formulation as the WELLSAWIPES GETRÄNKTE DESINFEKTIONSTÜCHER, is VAH-listed.

PLIWA Hygiene GmbH,
Pliwastraße 2, 34323 Malsfeld, Germany

Distributed by:
wellsamed GmbH
Brahestraße 16, 04347 Leipzig, Germany

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