WELLSATRAY LÖFFELMATERIAL

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Patientenzielgruppe

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en INSTRUCTIONS FOR USE

Patient target group

All persons for whom individual trays are required for treatment.

Intended user

Dentist

Composition

Qualitative data

Oligomeric urethane acrylate, vinyl ester resin, 2,2-dimethoxy-1,2-diphenylethane-1-one or by 2,2-dimethoxy-2-phenylacetophenone, (+)-2,3-dioxo-1,7,7-trimethyl-bicyclo-(2,2,1)-heptane di-bornane-2,3 dione, diphenyl(2,4,6-trimethylbenzoyl) phosphinoxide, 2-hydroxy-2-methyl-1-phenyl-propane-1-one, 2-ethylhexyl-4-(dimethylamino)benzoate, 2,2-diethoxyacetophenone, soda-lime glass, amorphous silica, dyes (Pigment Red 146 = CI 12485 / Pigment Blue 29 = CI 77007 / D&C Red 30 = CI 73360 / Pigment Blue CI 74160 / Pigment Violet 19 = C.I. 73900), flavourings (Optamint).

Quantitative data

In 100 g there are 19 g acrylates/methacrylates and 80 g glass.

Intended purpose of the product

WELLSATRAY Löffelmateriale is a light-curing tray material in the form of prefabricated, baseplates to make rigid, custom impression trays and bite-blocks, as well as for other indications. WELLSATRAY Löffelmateriale can be cured in all ultraviolet and Halogen light-curing units. The baseplates are supplied in a re-usable light-proof box that can be used for transporting or storing models, etc.

Indication

Custom impression tray

Contraindications

The material contains methacrylates. Do not use in patients with an allergy to methacrylates or other substances contained in the product.

Performance characteristics of the product

WELLSATRAY Löffelmateriale polymerises with light to form a solid and compact polymer.

Specific precautions

Avoid any contact between the eyes and skin and the uncured material, i.e. also with the inhibited layer.

Always use a dust extractor when trimming the tray!

Preparation/Manufacture

Mark the outline of the future tray on the model. Block out the undercuts on the model. To create a spacer for the impression material, block out the entire model (e.g. with a sheet of wax) and cut out stops. If the plaster is very dry, soak the model briefly in water or apply dye isolation or an alginate isolation. Remove the baseplate from the box and close it. Once the model has been blocked out, adapt the baseplate to it and cut to shape using, e.g. a scalpel. Use the remaining material to form a handle and attach it to the tray. Place the model and tray in a light-curing unit and cure without vacuum.

Remove the tray from the model and polymerize once again from the underside. If the tray has been cured without a vacuum. Use alcohol to remove the smear layer.

Times

Working at daylight: 20 minutes

Polymerisation in an ultraviolet or Halogen light-curing unit, depending on the unit in use:

1. The top of the tray, on the model: 3-5 minutes
2. The underside of the tray, off the model: 3-5 minutes

Trimming

Trim the custom tray with tungsten carbide burs and smooth it with sandpaper.

Hygiene

Immerse impressions taken with WELLSATRAY Löffelmateriale customs trays in a disinfectant solution based on glutar aldehyde. The time depends on the disinfectant manufacturer's directions.

The finished medical product is used by the dentist, e.g. for individual trays with impression material.

Caution

The product has been developed for use in the dental sector and must be used according to the directions. The manufacturer accepts no liability for damage that results from use other than the above. Furthermore the user is obliged under his own responsibility to examine the material for its suitability and its usage possibilities for the intended applications, especially if the latter are not detailed in the directions for use.

Side effects

Adverse side effects are not known.

Interactions

Interaction of the product with other medical devices, medicinal products or other substances is not known.

Storage

The Product should be stored between 0 °C - 22 °C.

Technical data

Flexural strength min. 130 N/mm²

Minimum shelf life

Do not use WELLSATRAY Löffelmateriale after the expiry date.

Variants

WELLSATRAY Löffelmateriale Upper jaw template:

WIT151 Natur/natural

WIT156 Blau/blue

WIT158 Rosa/pink

Disposal

Uncured material can be polymerised with light (sunlight or artificial light) and then disposed of as plastic waste.

All serious incidents related to the device shall be reported to the manufacturer and to the competent authority of the Member State where the user and/or the patient is established.

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