





Dear dental professionals,

It is with great pleasure that we present our latest catalogue 2018 – finally also available in English for our numerous international customers. At this very moment, you are taking part in a very important step for our brand. Broadening our scope of trading internationally was the next logical step in our company's development. Having received extensive demand from outside Germany, we are very happy to now act on a more global basis while at the same time preserving the importance of our local values.

With wellsamed being an owner-operated business, we bring a lot of personal expertise to the table deriving from yearlong work and experience in the international dental industry. The personal involvement furthermore results in a very personal responsibility for the quality of our dental products. It is always our most important goal to offer the highest standards possible. That is why we produce locally, always seeking the merits of world famous German excellence.

At the same time, we keep looking towards the future. These are exciting times for dental professionals and we are doing our best to support our customers in being the best in their field. State-of-the-art, top-notch products are what differentiates modern dentists from their peers. Our customers are the ones who inspire us to always look further and anticipate dental trends. In this catalogue you will find many new wellsamed products we have recently developed. To name one of many, our newly introduced wellsamed Supersoft 3 toothbrush is the result of many years of research and development.

Talking about innovations and future – what could be more important than caring for our children? wellsamed has dedicated many efforts to its kids product line wellsakids. We believe that oral health needs to be taken care of right from the start. And we also understand that it needs to be fun! wellsakids combines both dental purpose and cheery motivation so that children are ready and willing to learn how to take care of their oral hygiene. This learning process is supported by colourful designs, child-appropriate incentives and easy application. We are constantly working on improving our prophylaxis products as our clients' health is always the focus of our aspirations of continuous improvement.

For many years now, the combination of tradition and innovation has been the basis of our brand philosophy. We work hard to make quality products and services that improve our customers' practice. Now that you know what drives us forward every day, be inspired and benefit from wellsamed as your reliable, competent and long-standing partner.

Yours faithfully,

A handwritten signature in blue ink, appearing to be 'Denny Weidner', written over a light blue circular stamp or watermark.


Denny Weidner



IMPRESSION page 6



HYGIENE & DISPOSABLES page 12



PROPHYLAXIS page 25



KIDS PROPHYLAXIS page 38



REFRIGERANTS page 44



ENDODONTICS page 46



DENTAL JEWELLERY page 54




DENTAL RESTORATION page 59




INFECTION CONTROL page 69



DEVICES & PRACTICE EQUIPMENT page 77



LABORATORY page 88



PATIENT BROCHURES page 90



INSTRUMENTS page 91



IMPRESSION

TOP
seller

wellsachrom Alginate

- dental impression material with 4-phases-colour indicator
- dust-free, homogeneous
- colour change: violet to pink, pink to light blue, light blue to white
- mixing time: 35 sec.
- processing time: 45 sec.
- in-mouth setting time: 40 sec.
- setting time: 2 min.
- impression can be stored up to 100 hours if sealed hermetically
- DIN EN 21563/ISO 1563
- flavour: peppermint
- packaging: 450 g

REF	packaging
WIA 110	450 g



WIA 110



WIA 100

wellsaprint Basic Alginate

- fast setting alginate for complete and partial impressions
- smooth, creamy consistency with excellent detail reproduction
- dust-free, high tear strength
- great compatibility with dental stone
- free of lead and cadmium, high biocompatibility
- can still be poured after as much as 5 days time
- exceeds standards of ISO 1563/ADA18
- in-mouth setting time: 40 sec.
- setting time: 2 min.
- flavour: tutti frutti
- packaging: 500 g incl. measuring equipment



WIA 120

wellsaprint Ortho Alginate

- fast setting alginate specifically for orthodontic requirements
- pleasant for children because of short in-mouth setting time
- homogeneous, steadfast consistency to create necessary counter-pressure during impression taking
- dust-free, high tear strength
- great compatibility with dental stone
- free of lead and cadmium, high biocompatibility
- can still be poured after as much as 5 days time
- exceeds standards of ISO 1563/ADA18
- in-mouth setting time: 35 sec.
- setting time: 1 min. 45 sec.
- flavour: tutti frutti
- packaging: 450 g incl. measuring equipment

REF	colour	packaging
WIA 100	orange	500 g

REF	colour	packaging
WIA 120	yellow	450 g



WIT 151

wellsatray Tray Material

- light-curing hybrid composite resin
- curing time: both sides 3 - 5 min.
- universal plates for upper and lower jaw
- packaging: 50 pcs.

REF	colour	packaging
WIT 151	transparent	50 pcs.
WIT 156	blue	50 pcs.
WIT 158	rose	50 pcs.

Disposable Mixing Pad, PE-paper

- all-purpose paper
- also available with non-slip base (WMB 237)
- packaging: 1 pad of 50 sheets



WMB 207



WMB 212



WMB 237

WMB 215

Disposable Mixing Pad, parchment

- all-purpose paper
- packaging: 1 pad of 50 sheets (WMB 256) or 75 sheets (WMB 254, WMB 257)



WMB 256



WMB 254



WMB 257

REF	material	dimensions (W x L)	packaging
WMB 254	parchment paper	4 x 4 cm	75 sheets
WMB 256		6 x 6 cm	50 sheets
WMB 257		7 x 10 cm	75 sheets
WMB 207	PE-paper	7 x 9.5 cm	50 sheets
WMB 212		12.5 x 19 cm	50 sheets
WMB 215		15 x 24 cm	50 sheets
WMB 237	PE-paper anti-slip	7 x 9.5 cm	50 sheets

wellsakon A-Silicone Putty

- addition curing putty material with new filling materials ensuring a careful and tack-free mixing procedure
- suitable for both double mixing technique (single-step impression) and corrective impression (two-step impression technique)
- resilience: $\geq 99.5\%$
- processing time: variable up to 2 min.
- in-mouth setting time: 2 min.
- final hardness: Shore A 70
- mint-flavoured
- packaging: 300 ml pot base, 300 ml pot catalyst



REF	packaging
WIS 210	2 x 300 ml

WIS 210



WMS 101

Alginate Mixing Spatula

- flexible
- handy, curved handle with loop
- packaging: 1 spatula



REF	material	colour	size
WMS 101	plastic	assorted colours	19.5 cm

Mixing Spatula RS, non-flexible, curved

- for gypsums and alginates
- non-flexible
- packaging: 1 spatula



E730

REF	material	size
E730	steel, stainless	20 cm

WMB 104



WMB 103



WMB 101

Alginate Mixing Bowl

- for mixing alginates, gypsums and potting compounds
- smooth and self-cleaning surface
- silicone, flexible
- packaging: 1 bowl

REF	colour	dimensions
WMB 101	assorted colours	Ø 85 x 50 mm
WMB 102	assorted colours	Ø 115 x 70 mm
WMB 103	assorted colours	Ø 125 x 85 mm
WMB 104	assorted colours	Ø 140 x 110 mm

WMB 102



WIS 230

wellsakon A-Silicone putty soft 5:1

- addition curing putty material with new filling materials ensuring a careful and tack-free mixing procedure
- suitable for both double mixing technique (single-step impression) and corrective impression (two-step impression technique)
- 5:1 full cartridge for automatic mixing systems
- resilience: $\geq 99.5\%$
- processing time: variable up to 2 min.
- in-mouth setting time: 2 min.
- final hardness: Shore A 60
- mint-flavoured
- packaging: 2 x 380 ml cartridges 5:1, 1 fix cap

REF	colour	packaging
WIS 230	blue	2 x 380 ml, 1 Fix cap



WIS 251

wellsakon A-Silicone Correction Material, light

- light bodied, addition curing silicone
- extremely fine flowing material
- hydrophilic with optimal thixotropic properties for high precision
- excellent elastic properties prevent tearing of material during removal of impression from the mouth
- suitable for both double mixing technique (single-step impression) and corrective impression (two-step impression technique)
- processing time: variable up to 2 min.
- in-mouth setting time: 2 min.
- mint-flavoured
- packaging: 2 x 50 ml cartridges, 12 mixing-tips yellow (WIS 920)

REF	colour	packaging
WIS 251	pink	2 x 50 ml, 12 Mixing Tips



WIS 271

wellsakon A-Silicone Correction Material, medium

- medium bodied, addition curing silicone
- a pronounced hydrophilic nature and high thixotropy enable an optimal flow and exact detail reproduction in moist mouth condition
- suitable for both double mixing technique (single-step impression) and corrective impression (two-step impression technique)
- processing time: variable up to 2 min.
- in-mouth setting time: 2 min.
- mint-flavoured
- packaging: 2 x 50 ml cartridges, 12 mixing-tips green (WIS 922)

REF	colour	packaging
WIS 271	lilac	2 x 50 ml, 12 Mixing Tips



WBS 620

wellsabite Bite Registration, Hard A85

- silicone-based bite registration material in safety cartridge system
- excellent thixotropic, flows only under pressure
- excellent milling and cutting features due to high final hardness (Shore A85)
- high dimension stability
- long lasting, multiple-use bite registration material
- fixing of support and face-bow registration
- processing time: 30 sec.
- in-mouth setting time: 30 sec.
- mint-flavoured
- packaging: 2 x 50 ml cartridges, 12 mixing-tips green (WIS 922)

REF	colour	packaging
WBS 620	yellow	2 x 50 ml, 12 Mixing Tips



WBS 640

wellsabite Bite Registration, Heavy D45

- extra hard with a final hardness of: Shore D45
- extra fast-setting
- more comfortable for patients, avoids warped bite registrations
- high thixotropy and dimension stability
- simple handling thanks to milling and breaking
- long lasting, multiple-use bite registration material
- securing needle point registrations and bite ratios with a facebow
- processing time: 30 sec.
- in-mouth setting time: 30 sec.
- mint-flavoured
- packaging: 2 x 50 ml cartridges, 12 mixing-tips green (WIS 922)

REF	colour	packaging
WBS 640	turquoise	2 x 50 ml, 12 Mixing Tips

wellsamed Mixing Gun 1:1/2:1

- i.a. suitable for wellsakon a-silicone correction material light and medium, wellsabite bite registration
- packaging: 1 mixing gun



WIS 910

REF	packaging
WIS 910	1 piece

wellsamed Intra Oral Tips

- nose for mixing-tips
- packaging: 100 pcs.



WIS 951



WIS 952



WIS 953

REF	material	colour	size	packaging
WIS 951	plastic	transparent	Small	100 pcs.
WIS 952	plastic	transparent	Large	100 pcs.
WIS 953	plastic	transparent	X-Large	100 pcs.

Disposable Dispensing Syringe

- flexible, bendable tip
- tip can be shortened by knife or scissors according to individual preferences
- suitable for all impression materials
- packaging: 50 pcs.



WMS 201



REF	material	colour	size	packaging
WMS 201	plastic	transparent/white	5 ml	50 pcs.



wellsamed Mixing Tips

wellsamed Mixing Tips are compatible with a multitude of materials by original manufacturers and all wellsakon and wellsabite products (p. 6-8).

Packaging: 1 bag of 50 pcs. respectively 40 pcs. (WIS 941)

yellow

mixing ratio: 1:1/2:1
compatible with:

- Affinis light body and regular body (Coltène Whaledent)
- Aquasil LV/ULV (Dentsply)
- Honigum Automix Light and Light Fast, LuxaCore (DMG)

- Dimension Garant L and L Quick, Imprint II Garant Light Body/Regular Body, Imprint II Quick Step Light Body – Regular Body (3M ESPE)
- Flexitime Correct Flow, Provil Novo (Heraeus Kulzer)

pink

mixing ratio: 1:1/2:1
combinable with:

- Detaseal Quick and Light, Detaseal bite, One-time perfect (Detax)
- Examix NDS, Exafast NDS (GC)

- Betasil Vario (Müller Dental)
- Softreliner (Tokuyama)
- Registrado X-tra, Ufi Gel hard C, Ufi Gel SC (VOCO)

green

mixing ratio: 1:1/2:1
combinable with:

- Affinis heavy body (Coltène Whaledent)
- Greenbite Apple (Detax)
- Honigum Automix Mono (DMG)
- Imprint Bite, Imprint II Garant Monophase, Imprint II Quick Step Heavy, Imprint SBR (3M ESPE)

- Flexitime Monophase, Flexitime Heavy Tray, Memosil 2, Memoreg 2, Provil Novo (Heraeus Kulzer)
- Futar D occlusion, Futur System 50, Futar D System 50, Monopren transfer, Mucopren Soft, Mucopren transfer, Panasil Contact two in one (Kettenbach)

violet

mixing ratio: 1:1/2:1
combinable with:

- Aquasil Monophase and Heavy Body F/S (Dentsply)
- Detaseal Mono (Detax)

white

mixing ratio: 1:1/2:1
combinable with:

- Impregum Garant L Duosoft, Permadyne Garant 2:1 (3M ESPE)

light blue-orange

mixing ratio: 4:1/10:1
combinable with:

- Luxatemp (DMG)
- Predemp 3 Garant (3M ESPE)

- Prevision CB (Heraeus Kulzer)
- Omnitemp KB (Omnident)
- System c & b (Ivoclar Vivadent clinical)

yellow-blue

mixing ratio: 4:1/10:1
combinable with:

- Silasoft Direct/Express (Detax)
- Xantopren Comfort (Heraeus Kulzer)

REF	mixing ratio	material	colour	dimensions
WIS 920	1:1/2:1	plastic	yellow	Ø 4.2 x 70 mm
WIS 921	1:1/2:1	plastic	pink	Ø 5.4 x 83 mm
WIS 922	1:1/2:1	plastic	green	Ø 6.5 x 88 mm
WIS 923	1:1/2:1	plastic	violet	Ø 7.5 mm
WIS 924	1:1/2:1	plastic	white	-/-
WIS 931	4:1/10:1	plastic	light blue-orange	-/-
WIS 932	4:1/10:1	plastic	yellow-blue	-/-
WIS 941	5:1	plastic	yellow	-/-

WIS 920

WIS 921

WIS 922

WIS 923

WIS 924

WIS 931

WIS 932

WIS 941

wellsacup Paper Cups

- polymer coated paper cups
- capacity 180 ml/7 oz.
- various designs available
- packaging: carton of 1,000 cups



HYGIENE & DISPOSABLES

wellsamed Flower Set

package deal includes:

- 1 box of wellsacup paper cups "Flower" (WPC 007)
- 4 boxes of wellsamask Face Mask with Earloop "Flowers" (WFM 100)



WPC 007



WFM 100

REF	description	packaging
WFM 100S	Flower Set with WFM. 100 + WPC 007	1 set

wellsacup Paper Cups

- polymer coated paper cups
- capacity 180 ml/7 oz.
- various designs available
- packaging: box of 1,000 cups

**NEW
DESIGN**



WPC 006



WPC 001

WPC 007

WPC 008

WPC 005

WPC 004

Plastic Cups

- plastic
- robust, steady
- top quality
- suitable for all standardised cup dispensers
- bottom diameter: 4,5 cm
- top diameter: 6,3 cm
- capacity: 180 ml
- packaging: box of 3,000 cups

*Also
available as
colour set,
see page 22*



**NEW
COLOURS**

*Add to your
order!*

NEW



WCD 010

WCD 020

Cup Dispenser Stainless Steel

- classy dispenser, practical and sanitary
- made of stainless steel
- elegant und stable
- easy and fast removal of cups
- suitable for wall mounting
- a maximum of 70 cups, 180 ml and 150 ml
- holds dimensions (W x D x H): 85 x 85 x 360 mm
- packaging: 1 dispenser

REF	material	colour	packaging
WPC 110	plastic	white	3,000 pcs.
WPC 120		yellow	3,000 pcs.
WPC 135		light green	3,000 pcs.
WPC 150		light blue	3,000 pcs.
WPC 160		dark blue	3,000 pcs.
WPC 170		purple	3,000 pcs.
WPC 180		red	3,000 pcs.
WPC 185		pink	3,000 pcs.
WPC 190		orange	3,000 pcs.
WPC 999		assorted colours	3,000 pcs.
REF	material	colour	packaging
WPC 001	paper	Bubble	1,000 pcs.
WPC 004		Bunny	1,000 pcs.
WPC 005		Animal Train	1,000 pcs.
WPC 006		wellsakids	1,000 pcs.
WPC 007		Flowers	1,000 pcs.
WPC 008		Spots	1,000 pcs.
REF	material	colour	dimensions (W x D x H)
WCD 010	stainless steel	white	85 x 85 x 360 mm
WCD 020		silver	85 x 85 x 360 mm

wellsamask Tie-On Face Mask

- three-ply filter medium: middle ply with high filtration efficiency against bacteria
- integrated nose bridge
- latex free, hypoallergenic
- soft, non-irritating ties
- high wearing comfort
- Filtration Efficiency: BFE > 99.5 %
- EN 14683: Typ IIR
- packaging: box of 50 masks




WFM 360



wellsamask Face Mask Eco

- with earloop
- made of paper, single-layered
- packaging: box of 100 masks



WFM 001

NEW

*BFE > 99.5%
high filtration
efficiency*



WFM 090



WFM 105



WFM 100



wellsamask Face Mask with Earloop

- free of fiberglass and latex
- integrated nose bridge
- smooth and non-irritating non-woven material
- soft and skin-friendly elastic straps to put behind the ears
- high wearing comfort
- Filtration Efficiency: BFE > 99.5 %
- EN 14683: Typ IIR
- available in a variety of colours
- packaging: box of 50 masks



WFM 090

REF	description	colour	packaging
WFM 310	Tie-On Face Mask, non-woven, 3-ply	white	50 pcs.
WFM 320		yellow	50 pcs.
WFM 330		green	50 pcs.
WFM 360		blue	50 pcs.
WFM 380		pink	50 pcs.
WFM 090	Face Mask with Earloop, non-woven, 3-ply	wellsakids	50 pcs.
WFM 100		flowers	50 pcs.
WFM 105		cats	50 pcs.
WFM 110		white	50 pcs.
WFM 120		yellow	50 pcs.
WFM 130		green	50 pcs.
WFM 160		blue	50 pcs.
WFM 180		pink	50 pcs.
WFM 001		white	100 pcs.



Headrest Cover Standard

- 1-ply polyethylene film, 1-ply cellulose
- size (w x l): 33 x 25 cm
- packaging: box of 500 pcs.

REF	description	colour	packaging
WHC 301	Headrest Cover, standard	white	500 pcs.
WHC 302		yellow	500 pcs.
WHC 303		green	500 pcs.
WHC 305		aqua	500 pcs.
WHC 306		blue	500 pcs.
WHC 307		lavender	500 pcs.
WHC 308		pink	500 pcs.
WHC 309		peach	500 pcs.



WHC 210

Headrest Cover

- to protect headrest against contamination caused by sweat and haircare products
- 100 % pure cellulose and polyethylene film
- easy to take out of dispenser
- dimensions (w x l): 30 x 27 cm
- material: non-woven fabric
- packaging: box of 100 pcs.

REF	description	colour	packaging
WHC 210	Headrest Cover	white	100 pcs.



WHC 111



WHC 118



WHC 110

Headrest Cover

- disposable headrest covers for optimal hygiene
- soft, absorbent non-woven material, breathable
- protection against contamination caused by e.g. haircare products or sweat
- suitable for all common treatment units
- size: universal
- packaging: box of 100 pcs.

REF	description	colour	packaging
WHC 110	Headrest Cover	white	100 pcs.
WHC 111		pink	100 pcs.
WHC 118		blue	100 pcs.



**NEW
DESIGNS**



wellsabib Premium Patient Towels in the Dispenser

- for optimum protection and comfort of patients
- 2-ply cellulose, 1-ply waterproof polyethylene film
- high absorbency
- size (w x l): 33 x 45 cm
- packaging: 500 pcs. (4 dispenser of 125 pcs.)

REF	description	colour	packaging
WPB 500	Premium Patient Towels	black	500 pcs.
WPB 507		flower	500 pcs.
WPB 509		wellsakids	500 pcs.
WPB 510		white	500 pcs.
WPB 520		yellow	500 pcs.
WPB 530		green	500 pcs.
WPB 535		lime green	500 pcs.
WPB 550		light blue	500 pcs.
WPB 560		blue	500 pcs.
WPB 570		purple	500 pcs.
WPB 585		pink	500 pcs.
WPB 588		bordeaux	500 pcs.
WPB 590		orange	500 pcs.



Tray Cover

- to use on trays and in drawers
- highly absorbent, impedes flow of moisture
- chlorine-free, dioxin-free
- material: cellulose
- size (w x l): 28 x 36 cm
- packaging: box of 250 sheets

**NEW
COLOURS**

REF	description	colour	packaging
WCP 010	Tray Cover	white	250 sheets
WCP 020		yellow	250 sheets
WCP 030		green	250 sheets
WCP 040		purple	250 sheets
WCP 050		pink	250 sheets
WCP 060		blue	250 sheets



Tray Paper

- paper inlays for standard trays
- highly absorbent, high-wet-strength
- material: cellulose
- size (w x l): 18 x 28 cm
- packaging: box of 250 sheets

**NEW
COLOURS**

REF	description	colour	packaging
WTP 010	Tray Paper	white	250 sheets
WTP 020		yellow	250 sheets
WTP 030		green	250 sheets
WTP 040		purple	250 sheets
WTP 050		pink	250 sheets
WTP 060		blue	250 sheets




WBH 001



WBH 006

Plastic Bib Clip

- secure fixation of bib/napkin around the patient's neck during treatment
- two clips for attaching clothes protection
- packing: 1 bib clip




WBH 020

Coloured Plastic Bib Clip


- secure fixation of bib/napkin during treatment
- coloured chain in spiral shape
- elastic and stable plastic
- two metal clips for attaching clothes protection
- packaging: 1 bib clip



WBH 111



WBH 101



WBH 102

Metal Bib Clip, not adjustable

- secure fixation of bib/napkin during treatment
- two clamps keep the clothes guard in place
- length: WBH 101 = 37.5 cm/WBH 102 = 30 cm/WBH 111 = 45 cm
- WBH 102 with assorted coloured plastic ends in red, yellow, green or blue

REF	description	colour	length
WBH 001	Plastic Bib Clip	yellow	49 cm
WBH 006		blue	49 cm
WBH 101	Metal Bib Clip not adjustable	silver	37.5 cm
WBH 102		assorted colours	30 cm
WBH 111	Metal Bib Clip adjustable	silver	45 cm
REF	description	colour	
WBH 020	Coloured Plastic Bib Clip	multi-colour	



WNK 010

Wet Crepe Napkins

- protect patients against contamination and soginess
- enhanced grammage (about 26 g/m²)
- high wet strength, tearproof
- strong and soft
- size (w x l): 38 x 40 cm
- packaging: carton of 1,000 pcs.

REF	description	colour	dimensions (w x l)
WNK 010	paper, ca. 26 g/m ²	white	34 x 38 cm

Wellsamed Napkins

- 2- or 3-ply napkins made of cellulose
- 1/4 fold
- size (w x l): 38 x 40 cm
- packaging: 20 packs of 50 pcs./pack



WPB 412

WPB 462

WPB 413

REF	description	colour	dimensions (w x l)	packaging
WPB 412	cellulose, 2-ply	white	38 x 40 cm	1,000 pcs.
WPB 462		blue	38 x 40 cm	1,000 pcs.
WPB 413	cellulose, 3-ply	white	38 x 40 cm	1,000 pcs.




wellscape Patient Aprons, 2-ply

- optimum protection against contamination
- 100 % pure cellulose and polyethylene film
- with neckline
- with ties/laces to tie in the neck area
- perforated, on roll
- impervious to liquids
- dimensions (w x l): 53 x 60 cm
- packaging: 1 roll of 80 pcs./roll

wellsachild Bibs »Boy&Girl« or »wellsakids« or »Teeth«

- for young patients
- 2-ply
- 100 % pure cellulose and polyethylene film
- perforated, on roll
- with child-friendly motifs
- dimensions (w x l): 35 x 40 cm
- packaging: 1 roll of 100 pcs./roll



WCB 160



WCB 290




WCB 370

REF	description	colour	dimensions (w x l)	packaging
WPB 210	Patient Apron paper, 2-ply	white	54 x 60 cm	80 pcs. per roll
WPB 230		green	54 x 60 cm	
WPB 260		blue	54 x 60 cm	
REF	description	colour	dimensions (w x l)	packaging
WCB 160	wellsachild Bibs Teeth	blue	35 x 40 cm	100 pcs. per roll
WCB 190		apricot	35 x 40 cm	
WCB 260	wellsachild Bibs Boy&Girl	blue	35 x 40 cm	
WCB 290		apricot	35 x 40 cm	
WCB 370	wellsachild Bibs wellsakids	magenta	35 x 40 cm	

WPA 110

WPA 130

WPA 160



wellsplast Patient Aprons

- optimum protection against contamination
- 1-ply polyethylene film
- with neckline
- with ties/laces to tie in the neck area
- perforated, on roll
- impervious to liquids
- dimensions (w x l): 54 x 80 cm
- packaging: 1 roll of 200 pcs./roll



REF	description	colour	dimensions (w x l)	packaging
WPA 110	Patient Apron	white	54 x 80 cm	200 pcs. per roll
WPA 130		green	54 x 80 cm	
WPA 160		blue	54 x 80 cm	


Apron Dispenser

- to be filled with wellsamed Patient Aprons
- impact-resistant plastic
- easy to refill
- packaging: 1 dispenser



WAD 100

REF	material	colour	dimensions (W x D x H)
WAD 100	plastic	white	27.5 x 13 x 13.5 cm



WPT 110



WPT 250

WPT 230

WPT 210

wellsamed Paper Towels, 1-ply

- hand towels for daily use at the dental practice
- 100 % recycled material
- Z-fold
- 1-ply paper, embossed
- dimensions (w x l): 25 x 23 cm
- packaging: carton of 5,000 towels

wellsamed Paper Towels, 2-ply

- hand towels for daily use at the dental practice
- 100 % recycled material
- Z-fold
- double-layered paper, embossed
- dimensions (w x l): 25 x 23 cm
- packaging: carton of 3,200 towels

REF	description	colour	packaging
WPT 110	Paper Towels, 1-ply	natural white	5,000 pcs.
WPT 210	Paper Towels, 2-ply	bright white	3,200 pcs.
WPT 230		green	3,200 pcs.
WPT 250		blue	3,200 pcs.

Paper Towel Dispenser plastic

- ABS-plastic
- lockable
- capacity: 600 towels
- packaging: 1 dispenser

Paper Towel Dispenser

- stainless steel (silver) or sheet steel (white)
- lockable
- capacity: 600 towels
- packaging: 1 dispenser



A108360

A108362

A108361

REF	material	colour	dimensions (W x D x H)
A108360	plastic (ABS)	white	275 x 126 x 365 mm
A108362	sheet steel	white	250 x 125 x 330 mm
A108361	stainless steel polished	silver	250 x 125 x 330 mm



WCS 101

well-sazell Cellulose Pads

- bleached cellulose, DAB
- soft, abrasion-resistant, portioned
- highly absorbent
- packaging: 2 packs of 500 pads

REF	description	colour	dimensions (w x l)
WCS 102	cellulose, 8-ply	white	4 x 5 cm
WCS 101	cellulose, 12-ply	white	4 x 5 cm



WFP 004



Foam Pellets

- used for vitality test and application of liquids or pastes
- colour: white
- packaging: pack of 500, 1,000 or 3,000 pcs.

REF	material	size	packaging
WFP 004	foam	Ø 4 mm, round	3,000 pcs.
WFP 005		Ø 5 mm, round	1,000 pcs.
WFP 008		Ø 8 mm, round	500 pcs.
WFP 086		8 x 6 mm, conical	500 pcs.

6 packs
of 10 pcs.



WPT006/001

well-samed Facial Tissues

- supersoft and highly absorbent
- tissues made of 4-ply cellulose
- easy removal
- dimensions (w x l): 20 x 20 cm
- packaging: 6 packs of 10 pcs.

REF	description	colour	dimensions (w x l)
WPT 006/001	cellulose, 4-ply	white	20 x 20 cm



WCD 210

Cotton Roll Dispenser

- practical storage und protection against contamination
- apply slight pressure on lid to dispense the cotton rolls separately
- skidproof due to anti-slip base
- remaining cotton rolls are protected and uncontaminated
- suitable for sizes 1 and 2
- packaging: 1 dispenser

REF	material	colour	dimensions (W x D x H)
WCD 210	plastic	white	10 x 48 x 134 mm



WCD 110

Cotton Roll Dispenser Eco

- practical storage und protection against contamination
- skidproof due to anti-slip base
- removal of cotton rolls through pull-out drawer
- suitable for sizes 1, 2, and 3
- packaging: 1 dispenser



WCD 120



WCD 130



WCD 160



WCD 180

REF	material	colour	dimensions (W x D x H)
WCD 110	plastic	white/grey	100 x 50 x 100 mm
WCD 120		yellow	100 x 50 x 100 mm
WCD 130		green	100 x 50 x 100 mm
WCD 160		blue	100 x 50 x 100 mm
WCD 180		red	100 x 50 x 100 mm

wellsaroll Dental Cotton Rolls

- 100 % pure long-fibre cotton
- chlorine-free bleached
- lint-free
- colour: bright white
- high absorption capacity
- packaging: 1 bag of 300 g



WCR 001

WCR 002

WCR 003

REF	material	size	packaging
WCR 001	100% cotton	Nr. 1	300 g
WCR 002		Nr. 2	300 g
WCR 003		Nr. 3	300 g



**NEW
COLOURS**



SHOW YOUR COLOUR

Colour Set disposables

content:

- 500 saliva ejectors
- 500 patient towels
- 500 tray paper sheets
- 3,000 plastic cups
- 50 face masks with earloop
- 100 disposable toothbrushes
- 10 aspirator tips standard
- 2 cans of ice spray

REF	colour	content
VDA121	blue	over 4,500 items
VDA122	yellow	over 4,500 items
VDA123	green	over 4,500 items
VDA124	pink	over 4,500 items
VDA125	purple	over 4,500 items



WELLSA GLOVE

wellsaglove Examination Gloves Premium Nitril

- high-quality nitrile material
- latex-free, powder-free, non-sterile
- AQL 1.5
- sizes: x-small to x-large
- packaging: box of 100 gloves

wellsaglove Examination Gloves Premium Grip Latex

- high-quality examination gloves
- latex, powder-free, non-sterile
- roughened fingertips to ensure better handling (normal)
- AQL 1.5
- sizes: x-small to x-large
- packaging: box of 100 gloves



WEG 422

REF	description	colour	size
WEG 421	nitrile, powder-free	blue	XS
WEG 422			S
WEG 423			M
WEG 424			L
WEG 425			XL



WEG 823

REF	description	colour	size
WEG 821	latex, powder-free	white	XS
WEG 822			S
WEG 823			M
WEG 824			L
WEG 825			XL



WIC 278

WIC 274



WIC 306

WIC 302



wellsaject Disposable Injection Needles 27 G

- suitable for all cylindrical ampoule syringes
- sterile packaging
- easy handling and disposal
- each needle with sterility seal
- three-edge lancet grinding for smooth entry
- siliconised cannula
- flexible
- packaging: 100 pcs.



wellsaject Disposable Injection Needle 30 G

- suitable for all cylindrical ampoule syringes
- sterile packaging
- easy handling and disposal
- each needle with sterility seal
- three-edge lancet grinding for smooth entry
- siliconised cannula
- flexible
- packaging: 100 pcs.

REF	additional specification	dimensions	colour
WIC 274	27G, short	0.41 x 21 mm	yellow
WIC 275	27G, medium	0.41 x 25 mm	brown
WIC 276	27G, long	0.41 x 30 mm	orange
WIC 277	27G, x-long	0.41 x 35 mm	red
WIC 278	27G, xx-long	0.41 x 38 mm	purple
WIC 302	30G, xx-short	0.31 x 12 mm	blue
WIC 303	30G, x-short	0.31 x 16 mm	light blue
WIC 304	30G, short	0.31 x 21 mm	green
WIC 306	30G, long	0.31 x 25 mm	violet

wellsaliva Saliva Ejectors

- integrated metal wire
- wire is completely coated
- bendable and inherently stable
- transparent, removable tip
- length: 12.5 cm
- colours blue (WSE 165) and green (WSE 135) available in 15 cm length
- packaging: bag of 100 pcs.


REF	colour	length	packaging
WSE 100	transparent	12.5 cm	100 pcs.
WSE 110	white	12.5 cm	100 pcs.
WSE 120	yellow	12.5 cm	100 pcs.
WSE 130	green	12.5 cm	100 pcs.
WSE 160	blue	12.5 cm	100 pcs.
WSE 165	blue	15.0 cm (special length)	100 pcs.
WSE 135	green	15.0 cm (special length)	100 pcs.
WSE 180	pink	12.5 cm	100 pcs.
WSE 190	orange	12.5 cm	100 pcs.

wellsamed Aspirator Tip Cleaning Brush

- special brush to clean aspirator tips
- flexible wire
- loop at the end of the brush for sanitary storage
- packaging: pack of 10 pcs. respectively 6 pcs. (WAC 200)

NEW

Add to your order!



REF	dimensions	packaging
WAC 200	Ø 17 mm/197 mm length	6 pcs.
WAC 301	Ø 1.5 mm/120 mm length	10 pcs.
WAC 303	Ø 4.8 mm/160 mm length	10 pcs.
WAC 305	Ø 5 mm/230 mm length	10 pcs.

wellsamed Aspirator Tips

- autoclavable up to 134 °C
- low noise
- ergonomic design
- for adult treatment
- diameter 16 mm/length 125 mm
- packaging: pack of 10 pcs.



REF	colour	dimensions	packaging
WAC 110	white	Ø16 mm/125 mm length	10 pcs.
WAC 120	yellow	Ø16 mm/125 mm length	10 pcs.
WAC 130	green	Ø16 mm/125 mm length	10 pcs.
WAC 135	fresh-green	Ø16 mm/125 mm length	10 pcs.
WAC 160	blue	Ø16 mm/125 mm length	10 pcs.
WAC 170	purple	Ø16 mm/125 mm length	10 pcs.
WAC 180	red	Ø16 mm/125 mm length	10 pcs.
WAC 185	pink	Ø16 mm/125 mm length	10 pcs.
WAC 190	orange	Ø16 mm/125 mm length	10 pcs.

wellsamed Aspirator Tips Kids

- autoclavable up to 134 °C
- low noise
- ergonomic design
- for treatment of children
- diameter 16 mm/length 108 mm
- packaging: pack of 10 pcs.



REF	colour	dimensions	packaging
WAC 210	white	Ø16 mm/108 mm length	10 pcs.
WAC 220	yellow	Ø16 mm/108 mm length	10 pcs.
WAC 230	green	Ø16 mm/108 mm length	10 pcs.
WAC 235	fresh-green	Ø16 mm/108 mm length	10 pcs.
WAC 260	blue	Ø16 mm/108 mm length	10 pcs.
WAC 270	purple	Ø16 mm/108 mm length	10 pcs.
WAC 280	red	Ø16 mm/108 mm length	10 pcs.
WAC 285	pink	Ø16 mm/108 mm length	10 pcs.
WAC 290	orange	Ø16 mm/108 mm length	10 pcs.




**TOP
seller**

PROPHYLAXIS


*With protective cover
in convenient travel
box*

welsamed Travelbrush

- special long-wire brush to clean interdental spaces and hard-to-reach areas around bridges, implants, orthodontics and dentures
- long handle with thumb-rest
- bendable brush head
- each brush comes with protection cap in convenient travel box
- packaging: box of 6 pcs or 7 pcs incl. cap in colour assorted set (WTS 020)



REF	size	Ø brush head	Ø wire	colour	packaging
WTB 003	4S	1.7 mm	0.4 mm	white	6 pcs.
WTB 004	3S	1.9 mm	0.5 mm	orange	6 pcs.
WTB 005	2S	2.0 mm	0.55 mm	red	6 pcs.
WTB 006	S	3.3 mm	0.7 mm	yellow	6 pcs.
WTB 007	M	5 mm	0.9 mm	green	6 pcs.
WTB 011	L	2-3.8 mm	0.6 mm	blue	6 pcs.
WTB 013	2L	3-5.5 mm	0.8 mm	purple	6 pcs.
WTS 020	assorted				7 pcs.



REF	size	Ø brush head	Ø wire	colour	packaging
WMP 003	4S	1.9 mm	0.4 mm	grey	10 pcs.
WMP 004	3S	2.2 mm	0.5 mm	pink	10 pcs.
WMP 005	2S	2.4 mm	0.6 mm	yellow	10 pcs.
WMP 006	S	2.6 mm	0.7 mm	red	10 pcs.
WMP 007	M	2.5 – 3.0 mm	0.8 mm	blue	10 pcs.
WMP 008	L	3.0 – 5.0 mm	1.0 mm	green	10 pcs.
WMP 009	2L	4.0 – 6.0 mm	1.2 mm	black	10 pcs.
WMP 001	assorted sizes (4S - 2L)			assorted	10 pcs.



SPECIAL OFFER

incl. display,
see page 36



welsamed MiniPik Interdental Brush


- effective and gentle
- for cleaning of interdental spaces and posterior teeth, removal of plaque
- strong wire, twistable handle
- packaging: bag of 10 pcs.



welsamed Smilebrush Long-Wire

- special long-wire brush to clean interproximal spaces and hard-to-reach areas around bridges, implants, retainers, and dentures
- with long handle
- coated wire
- nylon bristles
- packaging: pack of 10 pcs. or bulk pack of 500 pcs.

REF	size	Ø brush head	Ø wire	colour	packaging
WSL 007	x-fine	2.7 mm	0.7 mm	green	10 pcs.
WSL 008	fine	4.0 mm	0.8 mm	blue	10 pcs.
WSL 009	medium	5.0 mm	0.85 mm	purple	10 pcs.
WSL 010	large	8.0 mm	0.9 mm	red	10 pcs.
REF	size	Ø brush head	Ø wire	colour	packaging
WSL 507	x-fine	2.7 mm	0.7 mm	green	500 pcs.
WSL 508	fine	4.0 mm	0.8 mm	blue	500 pcs.
WSL 509	medium	5.0 mm	0.85 mm	purple	500 pcs.
WSL 5010	large	8.0 mm	0.9 mm	red	500 pcs.



wellsabrush Interdental Starter Kit

- sizes: assorted (4S-2L)
- packaging: box of 10 assorted brush sizes, 1 white handle

wellsabrush Replacement Brushes

- refill for wellsamed Interdental Starter Kit (WSK 001)
- packaging: box of 10 pcs.

REF	size	Ø brush head	Ø wire	colour	packaging
WIB 003	4S	1.9 mm	0,4 mm	grey	10 pcs.
WIB 004	3S	2.2 mm	0,5 mm	pink	10 pcs.
WIB 005	2S	2.4 mm	0,6 mm	yellow	10 pcs.
WIB 006	S	2.6 mm	0,7 mm	red	10 pcs.
WIB 007	M	2.5 – 3.0 mm	0,8 mm	blue	10 pcs.
WIB 008	L	3.0 – 5.0 mm	1,0 mm	green	10 pcs.
WIB 009	2L	4.0 – 6.0 mm	1,2 mm	black	10 pcs.

REF	packaging
WSK 001	1 Kit

wellsamed Interdental Brush Holder »mono«

- single-ended handle
- slightly angled head
- screwtop to lock brush into position
- material: aluminium
- packaging: 1 holder



wellsamed Interdental Brush Holder »duo«

- double-ended handle
- slightly angled head
- screwtop to lock brush into position
- material: aluminium
- packaging: 1 holder



REF	material	colour	wrapping	packaging
WMH 414	aluminium	silver	Cello	1 piece
WMH 412		pink	Cello	1 piece
WMH 410		green	Cello	1 piece
WMH 413		gold	Cello	1 piece
WMH 411		blue	Cello	1 piece
WMH 424	aluminium	silver	Cello	1 piece
WMH 422		pink	Cello	1 piece
WMH 420		green	Cello	1 piece
WMH 423		gold	Cello	1 piece
WMH 421		blue	Cello	1 piece



WFU 12N



WTU 21M



WMW 51M



WFW 51M



WTU 12M

well-samed Floss

- effectively cleans interdental spaces
- frees teeth from plaque and food particles
- available in different versions

well-samed PTFE Tape

- effectively cleans interdental spaces
- frees teeth from plaque and food particles
- available in different versions

well-samed MagicFloss

- effectively cleans interdental spaces
- frees teeth from plaque and food particles
- expands when exposed to moisture
- packaging: 1 blister or 1 dispenser



well-samed Dental Floss Holder »MegaFloss«

- with cutter to cut the excess floss
- clean teeth by moving the holder gently up and down
- incl. 10 meters of PTFE Tape
- packaging: 1 holder

REF	description	wrapping	dimensions
WFU 10N	Floss, unwaxed, unflavoured	1 dispenser	10 m
WFU 12N		1 blister	2 x 10 m
WFU 51N		1 blister	50 m
WFW 200N	Floss, waxed, mint	1 dispenser	200 m
WFW 10M		1 dispenser	10 m
WFW 12M		1 blister	2 x 10 m
WFW 51M		1 blister	50 m
WTU 100N		1 dispenser	100 m
WTU 05M	PTFE Tape, mint	1 dispenser	5 m
WTU 25M		1 dispenser	25 m
WTU 21M		1 blister	25 m
WTU 12M		1 blister	2 x 5 m
REF	description	wrapping	dimensions
WMW 50M	Magicfloss, mint	1 dispenser	50 m
WMW 51M		1 blister	50 m
REF	description	wrapping	dimensions
WMF 01	Dental Floss Holder incl. PTFE Tape	blister	

wellsamed Floss Threaders

- helps threading floss under bridges, implants, orthodontics, and crowns
- packaging: 1 pack of 10 pcs.



WFT 010

REF	packaging
WFT 010	10 pcs.

wellsabox

- threepart box with sieve
- denture box, also perfectly suitable for orthodontic appliances
- material: plastic
- packaging: 1 box, threepart



WDC 211

REF	material	packaging
WDC 211	plastic	1 piece

wellsamed Retainer Box, Standard



- impact-resistant
- with eyelet to tie ribbon to box
- convenient click-closure

wellsamed Retainer Box, Mini



- dimensions Mini: 65 x 70 x 25 mm
- dimensions Standard: 70 x 75 x 35 mm
- packaging: 1 box

wellsamed Orthodontic Box

- with vent holes
- with eyelet
- packaging: 1 box



REF	material	colour	dimensions (W x D x H)	packaging
WBX 001	plastic	white	78 x 78 x 25 mm/Mini	1 piece
WBX 020		assorted colours	65 x 70 x 25 mm/Mini	1 piece
WBX 021		blue	65 x 70 x 25 mm/Mini	1 piece
WBX 022		pink	65 x 70 x 25 mm/Mini	1 piece
WBX 023		yellow	65 x 70 x 25 mm/Mini	1 piece
WBX 024		purple	65 x 70 x 25 mm/Mini	1 piece
WBX 025		assorted colours	70 x 75 x 35 mm/Standard	1 piece
WBX 026		blue	70 x 75 x 35 mm/Standard	1 piece
WBX 027		pink	70 x 75 x 35 mm/Standard	1 piece
WBX 028		yellow	70 x 75 x 35 mm/Standard	1 piece
WBX 029		purple	70 x 75 x 35 mm/Standard	1 piece
REF	material	colour	dimensions (W x D x H)	packaging
WBX 041	plastic	yellow	10 x 6.4 x 4.5 cm	1 piece
WBX 042		pink	10 x 6.4 x 4.5 cm	1 piece
WBX 043		light green	10 x 6.4 x 4.5 cm	1 piece
WBX 044		orange	10 x 6.4 x 4.5 cm	1 piece
WBX 045		dark blue	10 x 6.4 x 4.5 cm	1 piece
WBX 046		black	10 x 6.4 x 4.5 cm	1 piece

wellsamed Retainer/Denture Box, Mini or Standard



- impact-resistant
- convenient click-closure
- material: plastic
- dimensions Mini: 70 x 95 x 35 mm
- dimensions Standard: 70 x 95 x 55 mm
- packaging: 1 box

REF	colour	dimensions (W x D x H)	packaging
WBX 030	assorted colours	70 x 95 x 35 mm/Mini	1 piece
WBX 031	blue	70 x 95 x 35 mm/Mini	1 piece
WBX 032	pink	70 x 95 x 35 mm/Mini	1 piece
WBX 033	yellow	70 x 95 x 35 mm/Mini	1 piece
WBX 034	purple	70 x 95 x 35 mm/Mini	1 piece
WBX 035	assorted colours	70 x 95 x 55 mm/Standard	1 piece
WBX 036	blue	70 x 95 x 55 mm/Standard	1 piece
WBX 037	pink	70 x 95 x 55 mm/Standard	1 piece
WBX 038	yellow	70 x 95 x 55 mm/Standard	1 piece
WBX 039	purple	70 x 95 x 55 mm/Standard	1 piece

wellsamed Tooth-brush Multi-Head

- non-slip material
- single tuft brush at the tip for interproximal spaces and hard-to-reach areas
- longer bristles on the side to massage gums
- short bristles in the middle for cleaning tooth surfaces
- packaging: 1 piece



REF	colour	packaging
WTB 031	assorted colours	1 piece

wellsamed Brushing Timer with rectangular case

- 2-minute timer
- stand-up
- dimensions WBC 010:
20 x 25 x 66.5 mm
- assorted colours
- packaging: 1 timer



wellsamed Brushing Timer without case

- 2-minute timer
- for setting up
- dimensions WBC 013:
Ø 20 x 90 mm
- assorted colours
- packaging: 1 timer



REF	description	colour	dimensions	packaging
WBC 010	Brushing Timer with rectangular case	assorted colours	20 x 25 x 66.5 mm	1 piece
WBC 013	Brushing Timer without case	assorted colours	Ø 20 x 90 mm	1 piece

wellsamed Toothbrush cup

- available in white, pink, blue
- Ø 7 x 9.5 cm
- material: plastic
- see-through
- packaging: 1 cup



REF	description	colour	packaging
WPC 210	Toothbrush Cup	white	1 piece
WPC 260		blue	1 piece
WPC 285		pink	1 piece



WTB 081

wellsamed Toothbrush »Soft«

- microfine bristles
- narrow handle
- small brush head
- particularly suitable for sensitive gums
- assorted colours
- packaging: 1 piece

REF	colour	packaging
WTB 081	assorted colours	1 piece



WTB 095



WTB 096

wellsamed Toothbrush CleanX

- cross-angled bristles
- extra-long, coloured bristles
- tuft at the tip of the brush head
- ergonomic handle
- assorted colours
- packaging: 1 piece

REF	wrapping	packaging
WTB 096	cello	1 piece
WTB 095	blister	1 piece



WWS 001

wellsastrips Teeth Whitening without Peroxide

- 30 minutes express
- 14 applications
- also suitable for sensitive teeth
- packaging: 1 set of 28 strips (14 x upper jaw, 14 x lower jaw)

REF	description	packaging
WWS 001	wellsastrips Teeth Whitening	1 piece
WTB 097	Black'n White Charcoal Toothbrush	1 piece

wellsamed Black'n White Charcoal Toothbrush

- charcoal infused bristles to remove plaque and discolourations (e.g. tea, coffee, red wine) gently and naturally
- single-component toothbrush, design "Swing" with straight cut
- carefully end-rounded bristles
- packaging: 1 piece, blister



NEW



WTB 097

SUPERSOFT 3

wellsamed SuperSoft 3

- with over 2,000 soft Tynex®-filaments
- particularly gentle on gums
- compact brush head
- plane bristles
- gummed and ergonomic handle
- packaging: 1 piece

compact head

0.15 mm tynex filaments

ergonomic handle



WTB 310

REF	colour	packaging
WTB 310	assorted colours	1 piece

NEW



WTB 351

wellsamed Denture Brush

- brush head with two sides for different purposes
- dense, high-/low-cut bristles for cleaning denture surfaces
- single-tufted head for cleaning interproximal spaces and denture ridges
- assorted colours
- packaging: 1 piece

REF	colour	packaging
WTB 351	assorted colours	1 piece



WTB 105

wellsamed Tooth-brush Megasoft

- extrasoft, rounded bristles
- specifically suitable for sensitive gums
- assorted colours
- packaging: 1 piece

REF	colour	packaging
WTB 105	assorted colours	1 piece

SPECIAL OFFER

incl. display,
see page 37

Demonstration Model

- demonstrative model to teach the correct method of brushing and flossing
- set of 28 teeth, with metal strap articulator
- possibility to demonstrate overbite and underbite
- packaging: 1 model with demonstration toothbrush



REF	dimensions (W x D x H)	packaging
WDM 010	10 x 12 x 10 cm	1 Set
REF	colour	packaging
WTB 311	assorted colours	1 piece

wellsamed Travel Toothbrush

- 2-component travel toothbrush with cap
- cap transforms into into handle
- dual-coloured bristles
- assorted colours
- packaging: 1 piece



wellsabrush Disposable Toothbrushes

- single use
- coated with mint tooth-paste
- available in various colours
- ideal for travelling, dental practices, hotels, care facilities
- packaging: box of 100 pcs.



REF	colour	packaging
WDB 020	yellow-transparent	100 pcs.
WDB 025	yellow green-transparent	100 pcs.
WDB 030	green-transparent	100 pcs.
WDB 060	blue-transparent	100 pcs.
WDB 070	purple-transparent	100 pcs.
WDB 080	red-transparent	100 pcs.
WDB 085	pink-transparent	100 pcs.
WDB 095	orange-transparent	100 pcs.
WDB 010	assortment	100 pcs.

wellsamed Universal Tongue Cleaner

- triple row, with lamellae
- handy design with curved handle
- material: plastic
- packaging: 1 piece



wellsamed Tongue Cleaner Clean 2 in 1

- double function: tongue cleaner and tongue scraper
- ergonomically designed
- flexible head
- packaging: 1 piece



wellsamed Tongue Cleaner Expert

- packaging: 1 piece



wellsamed Stain Remover

- safe removal of tea and coffee stains on teeth
- rubber tip
- various colours
- assorted colours
- packaging: 1 piece, blister



REF	colour	packaging
WTC 015	turquoise	1 piece
WTC 011	assorted colours	1 piece
WTC 020	yellow	1 piece

REF	colour	packaging
WSR 019	assorted colours	1 piece

Prophy Brushes

- for contra-angles
- fine bristles
- gentle and in-depth cleaning
- packaging: box of 50 pcs.



REF	size	material	shape	packaging
WZR 062	normal, flat	wool	SP Cup	50 pcs.
WZR 061		nylon, soft	SP Cup	50 pcs.
WZR 052		natural	SP Cup	50 pcs.
WZR 051		nylon	SP Cup	50 pcs.
WZR 053	small, flat	nylon	SP Brush	50 pcs.
WZR 055	small, tapered	nylon	SP Brush	50 pcs.
WZR 054	small, flat	natural	SP Brush	50 pcs.
WZR 110	normal, flat	nylon, megasoft	CC Cup	50 pcs.
WZR 111		nylon, normal	CC Cup	50 pcs.
WZR 112		nylon, soft	CC Cup	50 pcs.
WZR 113		nylon, supersoft	CC Cup	50 pcs.
WZR 120	small, flat	nylon, megasoft	CC Brush	50 pcs.
WZR 122		nylon, soft	CC Brush	50 pcs.
WZR 121		nylon, normal	CC Brush	50 pcs.
WZR 123		nylon, supersoft	CC Brush	50 pcs.
WZR 131		silicone carbide	SK Brush	50 pcs.
WZR 132	normal, tapered	silicone carbide	SK Brush	50 pcs.
WZR 130	normal, flat	silicone carbide	SK Cup	50 pcs.



Order your free sample

- SP Cup	WZR 051 Muster
- CC Cup	WZR 111 Muster
- SK Cup	WZR 130 Muster

Prophy Cups

- webbed and ribbed (WPC 302), soft
- attach to contra-angle
- prevents contamination of handpcs. and contra-angles
- not suitable for sterilisation
- packaging: box of 144 pcs.

REF	colour	shape	packaging
WPC 301	grey	-/-	144 pcs.
WPC 302	grey	ribbed inside	144 pcs.

WPP 004

WPP 025

wellsaprophy Prophylactic Paste

- for final polishing of natural teeth and fillings
- effective protection against cavity and periodontal diseases
- gentle removal of discolourations and plaque from teeth surfaces
- mint flavour
- flouride 0.2 %
- free of gluten and phthalate
- packaging: 1 tube of 95 ml

REF	grit	colour	packaging
WPP 004	x-fine (RDA 40)	orange	95 g
WPP 025	coarse (RDA 250)	blue	95 g

wellsamed Prophy Angles

- disposable contra-angle with 105° angle
- fits most standard handpcs.
- 1.200 - 3.000 rpm
- packaging: pack of 100 pcs.

Order your free sample

- Prophy Angles regular A119986 Muster
- Prophy Angles soft A119987 Muster

REF	description	colour	packaging
A119986	Prophy Angles	regular, lavender	100 pcs.
A119987		soft, green	100 pcs.



800 Piks



wellsamed MiniPik incl. Display

- counter display for dental practice with MiniPik Interdental Brushes
- dimensions (W x D x H): 286 x 125 x 590 mm
- packaging: 1 display, compatible hooks, 8 x 10 packs of MiniPik Interdental Brushes (size: 4S-2L)



REF	description	dimensions (W x D x H)
WDP 101	800 MiniPiks incl. Display	286 x 125 x 580 mm

SPECIAL OFFER

incl. Display Free of Charge

96 toothbrushes

96 wellsamed Megasoft incl. Display

- counter display for dental practice with wellsamed toothbrushes Megasoft
- dimensions (W x D x H): 286 x 125 x 580 mm
- packaging: 1 display, compatible hooks, 96 blister pack Megasoft



REF	description	dimensions (W x D x H)
WDP 102	96 Megasoft Toothbrushes incl. Display	286 x 125 x 580 mm



**NEW
DESIGN**

KIDS PROPHYLAXIS

wellsacup paper cups

- polymer coated paper cups
- capacity 180 ml/7 oz.
- various designs
- packaging: carton of 1,000 cups



WPC 008

WPC 004



WPC 006

REF	material	colour	packaging
WPC 004	paper	Bunny	1,000 pcs.
WPC 008		polka dots	1,000 pcs.
WPC 006		wellsakids	1,000 pcs.

Milk Tooth Box

- tooth box for safekeeping
- tooth-shaped, lockable
- with soft ribbon to be worn as a necklace
- assorted colours
- dimensions: 16 x 18 mm
- packaging: 1 box



WBX 010

ages
3 - 8

4 - 36
months

ages
2 - 4

ages
3 - 8

wellsamed kids Toothbrush Fantasy

- short head and super soft bristles
- ergonomic handle for child-friendly handling
- for children aged 4 - 36 months (WKB 091) or ages 3 - 8 (WKB 092)
- packaging: 1 piece



WKB 092

WKB 091

wellsamed kids Toothbrush Karli

- short head with multituft cleaning section and flat cut
- soft Tynexbristles
- for children ages 2-4
- packaging: 1 piece



WKB 012



WKB 010



wellsamed kids Toothbrush Eco

- ergonomic handle with thumb-rest for child-friendly handling
- for children ages 3 - 8
- packaging: 1 piece

ribbon

- with eyelet
- packaging: 1 piece



LMB 10

NEW

wellsamed kids Toothbrush

- soft Tynexbristles
- for children ages 0 - 6 years
- available in 2 different colours and designs
- packaging: 1 piece



WKB 013

ages
0 - 6



Milk Tooth Box

- tooth container to safekeep lost baby teeth
- see-through
- packaging: 1 box

REF	description	dimensions (W x H)	packaging
WBX 010	Milk Tooth Box	16 x 18 mm	1 piece
REF	description		packaging
WKB 010	kids Toothbrush	ages 3 - 8	1 piece
WKB 012		ages 2 - 4	1 piece
WKB 013		ages 0 - 6	1 piece
WKB 091	kids Toothbrush Fantasy	4 - 36 months	1 piece
WKB 092		ages 3 - 8	1 piece
REF	description	colour	packaging
LMB 02	Milk Tooth Box	transparent	1 piece
LMB 03		orange	1 piece
LMB 04		red	1 piece
LMB 05		blue	1 piece
LMB 06		mintgreen	1 piece
LMB 10	ribbon	blue	1 piece

wellsamed Infant Brush PREMIUM

- teething ring with bumps on one end
- soft bristles on the other end
- assorted colours
- packaging: 1 piece



WKB 004

NEW**wellsamed Turtle Infant Brush**

- turtle-shaped with dual-sided, short, and small bristles
- assorted colours: blue or pink
- packaging: 1 piece



WKB 003

NEW**wellsamed Infant Brush ECO**

- 28 tuft head on a 2 finger size oval handle
- with soft 0.005 nylon
- facilitates transition to conventional children toothbrushes
- assorted colours
- packaging: 1 piece



WKB 005

NEW**wellsamed Grip Reinforcement for Toothbrushes**

- soft sponge rubber handle, washable, multipurpose
- suitable for people with disabilities and fine motor skills issues
- packaging: 1 piece



WTB 100

NEW

REF	description	packaging
WKB 001	wellsakids Finger Toothbrush	1 piece
WKB 003	wellsamed TURTLE Infant Brush	1 piece
WKB 004	wellsamed Infant Brush PREMIUM	1 piece
WKB 005	wellsamed Infant Brush ECO	1 piece
WTB 100	wellsamed Grip Reinforcement for Toothbrushes	1 piece
WTE 106	wellsamed FunnyClean for Kids Electric Toothbrush, butterfly	1 piece
WTE 108	wellsamed FunnyClean for Kids Electric Toothbrush, dog	1 piece
WTE 110	wellsamed FunnyClean Brush Heads	2 pcs.

wellsamed FunnyClean for Kids Electric Toothbrush

- battery-operated toothbrush with 2 soft, replaceable brush heads
- specially designed for children from 5 years
- 8,800 oscillating movements per minute
- motif: dog, blue (WTE 108)
- motif: butterfly, pink (WTE 106)
- packaging: 1 electric toothbrush (batteries not included)

wellsakids Finger Toothbrush

- gentle cleaning of baby's gums and first teeth
- washable, hygienic
- 100 % silicone, without plasticisers
- assorted colours
- packaging: 1 piece, blister

NEW

WTE 106

WTE 108

wellsamed FunnyClean Brush Heads

- colour: white
- packaging: 2 brush heads, blister

WTE 110



WKB 001





age 7
and above

wellsakids Children Toothbrush

- improved plaque removal
- small brush head, child-friendly handle
- soft bristles
- children ages 3-7 or age 7 and above
- packaging: 1 piece, blister

ages
3-7



WKB 023

WKB 022

wellsakids Toothbrush cup

- assorted motifs (girls, boys)
- available in white, pink, and blue
- Ø 7 x 9.5 cm
- packaging: 1 cup



WPC 360



WPC 385

WPC 310

REF	description	colour	packaging
WPC 310	Toothbrush cup, girls/boys	white	1 piece
WPC 385		pink	1 piece
WPC 360		blue	1 piece
WPC 310A	Toothbrush cup Set	white	1 Set
WPC 385B		pink	1 Set
WPC 460C		blue	1 Set
REF	description	packaging	
WKB 022	Children Toothbrush	ages 3-7	1 piece
WKB 023		+7 years	1 piece

wellsakids Toothbrush cup Set Putzi

- dimensions: Ø 7 x 9.5 cm
- material: plastic
- packaging: 1 cup + 1 toothbrush (WKB 022) + 1 PUTZI toothpaste



WPC 310A

wellsakids Toothbrush cup Set STAR WARS

- dimensions: Ø 7 x 9.5 cm
- material: plastic
- packaging: 1 cup + 1 toothbrush (WKB 022) + 1 Stages STAR WARS toothpaste



WPC 460C

NEW

wellsakids Toothbrush cup Set FROZEN

- dimensions: Ø 7 x 9.5 cm
- material: plastic
- packaging: 1 cup + 1 toothbrush (WKB 022) + 1 Stages FROZEN toothpaste



WPC 385B

wellsakids Bibs

- dimensions (w x l): 35 x 40 cm
- packaging: 1 roll of 100 pcs./roll

NEW

wellsamask Face Mask with Earloop

- smooth and non-irritating non-woven material
- soft and skin-friendly elastic straps to put behind the ears
- packaging: box of 50 masks

wellsacup paper cups

- polymer coated paper cups
- capacity 180 ml/7 oz.
- various designs
- packaging: carton of 1,000 cups

WCB 370

WFM 090

WPC 006

**Waiting at the dentists
can be great fun**

wellsakids Coloured Pencils

- wooden coloured pencils in 6 different colours
- length: 8.7 cm
- packaging: 6 pencils incl. sharpener

wellsakids Colouring Book

- adventure story with dental care advices
- 12 pages to colour
- packaging: 1 book



WCP 806



WKP 010



REF	description	colour	packaging
WCP 806	Coloured Pencils	6 different colours	6 pcs.
WCB 370	wellsakids Bibs	magenta	100 pcs.
WFM 090	wellsamask	wellsakids	50 pcs.
WPC 006	wellsacup	wellsakids	1,000 pcs.
REF	description	packaging	
WKP 010	Colouring Book	1 piece	

wellsakids Brushing Timer

- 2-minute-timer
- round and rectangular
- packaging: 1 timer



WBC 011



WBC 012



WBC 014

wellsamed Denture Flapper

- funny game for children
- packaging: 12 pcs.



WKM 212

wellsakids Certificate

- certificate of bravery for children
- blank space for practice stamps and the young patient's name
- packaging: 50 certificates



WKU 005

wellsakids Bouncing Ball

- bouncing balls in 4 different colours
- dimensions: Ø 26 mm
- packaging: 1 pack of 400 pcs.



WKM 106

REF	description	shape	dimensions
WBC 011	Brushing Timer	round	Ø 85 x 13 mm
WBC 012		rectangular	80 x 110 x 13 mm
WBC 014		rectangular	80 x 110 x 13 mm
REF	description	packaging	
WKU 005	Certificate	50 pcs.	
WKM 212	wellsamed Denture Flapper	12 pcs.	
REF	description	colour	packaging
WKM 106	Bouncing Ball	assorted colours	400 pcs.



NEW

REFRIGERANTS

wellsafresh Oshibori

- refreshing towel made of premium terry cloth
- for hot and cold application
- without alcohol
- smells fresh as dew
- packaging: 1 refreshing towel

REF	description	packaging
WOT 001	wellsafresh Oshibori	1 piece



WOT 001

**TOP
seller**

wellsacool Cooling Compress

- handy compress with gel filling
- use in case of inflammations, swellings and toothache etc.
- cooling effect lasts about 30 minutes
- packaging: 1 pack



WCP 901

wellsacool Ice Spray

- contains hydrocarbon mixture, CFC-free
- release temperature: -45 °C
- exact application due to fine spray head
- no expiration date
- in accordance with directive 93/42/EWG, class 1
- mint-flavoured
- packaging: 1 bottle of 200 ml



WIS 101

REF	description	colour	packaging
WCP 901	wellsacool Cooling Compress	blue	1 piece
REF	description		packaging
WIS 101	wellsacool Ice Spray		200 ml



ENDODONTICS

wellsastand Foam Inserts

- for interim storage of endodontic instruments
- round foam inserts
- packaging WFI 001: 1 pack of 50 pcs.
- packaging WFI 002: 1 pack of 25 pcs.
- packaging WFI 003: 1 pack of 25 pcs.



REF	compatibility	dimensions	packaging
WFI 001	Endo Cleaner, CleanStand, CleanGrip, InterimStand	H 0,5 cm, Ø 5 cm	50 pcs.
WFI 002	CleanStand, CleanGrip, InterimStand	H 1 cm, Ø 5 cm	25 pcs.
WFI 003	Twin-Bloc	H 1 cm, Ø 6,5 cm	25 pcs.

wellsamed Files

- considerably enhanced flexibility compared to steel instruments
- high resilience
- ISO colour-coded plastic handle
- pre-assembled silicone-stops
- packaging: 1 box of 6 pcs.

wellsamed NiTi K-File



ISO	ISO 8	ISO 10	ISO 15	ISO 20	ISO 25	ISO 30	ISO 35	ISO 15-40
	grey	purple	white	yellow	red	blue	green	assorted
21 mm	WNF 1108	WNF 1110	WNF 1115	WNF 1120	WNF 1125	WNF 1130	WNF 1135	WNF 1101
25 mm	WNF 1208	WNF 1210	WNF 1215	WNF 1220	WNF 1225	WNF 1230	WNF 1235	WNF 1201

ISO	ISO 40	ISO 45	ISO 50	ISO 55	ISO 60	ISO 70	ISO 80	ISO 45-80
	black	white	yellow	red	blue	green	black	assorted
21 mm	WNF 1140	WNF 1145	WNF 1150	WNF 1155	WNF 1160	WNF 1170	WNF 1180	WNF 1102
25 mm	WNF 1240	WNF 1245	WNF 1250	WNF 1255	WNF 1260	WNF 1270	WNF 1280	WNF 1201


wellsamed NiTi H-File



ISO	ISO 8	ISO 10	ISO 15	ISO 20	ISO 25	ISO 30	ISO 35	ISO 15-40
	grey	purple	white	yellow	red	blue	green	assorted
21 mm	WNF 2108	WNF 2110	WNF 2115	WNF 2120	WNF 2125	WNF 2130	WNF 2135	WNF 2101
25 mm	WNF 2208	WNF 2210	WNF 2215	WNF 2220	WNF 2225	WNF 2230	WNF 2235	WNF 2201

ISO	ISO 40	ISO 45	ISO 50	ISO 55	ISO 60	ISO 70	ISO 80	ISO 45-80
	black	white	yellow	red	blue	green	black	assorted
21 mm	WNF 2140	WNF 2145	WNF 2150	WNF 2155	WNF 2160	WNF 2170	WNF 2180	WNF 2102
25 mm	WNF 2240	WNF 2245	WNF 2250	WNF 2255	WNF 2260	WNF 2270	WNF 2280	WNF 2202

wellsamed NiTi Reamer



ISO	ISO 15	ISO 20	ISO 25	ISO 30	ISO 35	ISO 40	ISO 15-40
	white	yellow	red	blue	green	black	assorted
21 mm	WNF 3115	WNF 3120	WNF 3125	WNF 3130	WNF 3135	WNF 3140	WNF 3101
25 mm	WNF 3215	WNF 3220	WNF 3225	WNF 3230	WNF 3235	WNF 3240	WNF 3201

wellsamed NiTi Spreader

- considerably enhanced flexibility compared to steel instruments
- high resilience
- ISO colour-coded plastic handle
- smooth conical shape with pointed tip
- designed for lateral condensation of gutta percha and other filling materials
- pre-assembled silicone-stops
- packaging: 1 box of 6 pcs.



ISO	ISO 15	ISO 20	ISO 25	ISO 30	ISO 35	ISO 40	ISO 15-40
	white	yellow	red	blue	green	black	assorted
21 mm	WNF 4115	WNF 4120	WNF 4125	WNF 4130	WNF 4135	WNF 4140	WNF 4101
25 mm	WNF 4215	WNF 4220	WNF 4225	WNF 4230	WNF 4235	WNF 4240	WNF 4201

wellsamed NiTi Plugger

- considerably enhanced flexibility compared to steel instruments
- high resilience
- smooth conical shape with flat tip
- designed for both vertical and lateral compacting of gutta percha and other filling materials
- ISO colour-coded plastic handle
- pre-assembled silicone-stops
- packaging: 1 box of 6 pcs.



ISO	ISO 15	ISO 20	ISO 25	ISO 30	ISO 35	ISO 40	ISO 15-40
	white	yellow	red	blue	green	black	assorted
21 mm	WNF 5115	WNF 5120	WNF 5125	WNF 5130	WNF 5135	WNF 5140	WNF 5101
25 mm	WNF 5215	WNF 5220	WNF 5225	WNF 5230	WNF 5235	WNF 5240	WNF 5201

wellsamed NiTi Large Taper

- considerably enhanced flexibility compared to steel instruments
- high resilience
- ISO colour-coded plastic handle
- pre-assembled silicone-stops
- packaging: 1 box of 6 pcs.



ISO	D1	D2	DX	M1	M2	M3	D1-M3
	purple	white	orange	light green	red	blue	assorted
21 mm	WNF 6103	WNF 6104	WNF 6105	WNF 6106	WNF 6107	WNF 6108	WNF 6101
25 mm	WNF 6203	WNF 6204	WNF 6205	WNF 6206	WNF 6207	WNF 6208	WNF 6201

wellsamed Paper Points Color 2%

- composition: 100 % paper
- highly absorbent
- for removing moisture from the root canal or applying medicaments
- guaranteed consistency, smooth surface
- fibre-free
- ISO colour-coded and sterilised
- packaging: 1 box of 200 pcs.



ISO 15	ISO 20	ISO 25	ISO 30	ISO 35	ISO 40	ISO 15-40
white	yellow	red	blue	green	black	assorted
WPP 315	WPP 320	WPP 325	WPP 330	WPP 335	WPP 340	WPP 301
ISO 45	ISO 50	ISO 55	ISO 60	ISO 70	ISO 80	ISO 45-80
white	yellow	red	blue	green	black	assorted
WPP 345	WPP 350	WPP 355	WPP 360	WPP 370	WPP 380	WPP 302

wellsamed Gutta Percha Points Color 2%

- composition: gutta percha 25 %, zinc oxide 62.5 %, barium sulfate 12 %, coloring agent 0.5 %
- for filling a prepared root canal prior to tooth restoration
- ISO standard radiopaque
- ISO colour-coded
- packaging: 1 box of 120 pcs.



ISO 15	ISO 20	ISO 25	ISO 30	ISO 35	ISO 40	ISO 15-40
white	yellow	red	blue	green	black	assorted
WGP 315	WGP 320	WGP 325	WGP 330	WGP 335	WGP 340	WGP 301
ISO 45	ISO 50	ISO 55	ISO 60	ISO 70	ISO 80	ISO 45-80
white	yellow	red	blue	green	black	assorted
WGP 345	WGP 350	WGP 355	WGP 360	WGP 370	WGP 380	WGP 302



wellamed Paper Points Isotaper 4 %

- composition: 100 % paper
- use in root canal to remove moisture or to introduce medicaments
- ISO colour-coded
- do not use in autoclave
- packaging: 1 box of 100 pcs.

ISO 15	ISO 20	ISO 25	ISO 30	ISO 35
white	yellow	red	blue	green
WPP 415	WPP 420	WPP 425	WPP 430	WPP 435
ISO 40	ISO 45	ISO 20-45		
black	white	assorted		
WPP 440	WPP 445	WPP 401		



wellamed Gutta Percha Points Isotaper 4 %

- composition: gutta percha 25 %, zinc oxide 62.5 %, barium sulfate 12 %, coloring agent 0.5 %
- for filling a prepared root canal prior to tooth restoration
- ISO standard radiopaque
- ISO colour-coded
- not autoclavable
- packaging: 1 box of 60 pcs.

ISO 15	ISO 20	ISO 25	ISO 30	ISO 35
white	yellow	red	blue	green
WGP 415	WGP 420	WGP 425	WGP 430	WGP 435
ISO 40	ISO 45	ISO 20-45		
black	white	assorted		
WGP 440	WGP 445	WGP 401		



wellamed Paper Points Isotaper 6 %

- composition: 100 % paper
- ISO colour-coded
- not autoclavable
- packaging: 1 box of 100 pcs.

ISO 15	ISO 20	ISO 25	ISO 30	ISO 35
white	yellow	red	blue	green
WPP 615	WPP 620	WPP 625	WPP 630	WPP 635
ISO 40	ISO 45	ISO 20-45		
black	white	assorted		
WPP 640	WPP 645	WPP 601		



wellamed Gutta Percha Points Isotaper 6 %

- composition: gutta percha 25 %, zinc oxide 62.5 %, barium sulfate 12 %, coloring agent 0.5 %
- use in root canal
- ISO colour-coded
- not autoclavable
- packaging: 1 box of 60 pcs.

ISO 15	ISO 20	ISO 25	ISO 30	ISO 35
white	yellow	red	blue	green
WGP 615	WGP 620	WGP 625	WGP 630	WGP 635
ISO 40	ISO 45	ISO 20-45		
black	white	assorted		
WGP 640	WGP 645	WGP 601		

Metapaste

- premixed, water soluble calcium hydroxide paste with barium sulfate
- formation of a hard tissue barrier (Apexification)
- temporary root canal filling
- internal and external root resorption
- excellent antibacterial effect and radiopacity
- packaging: 2.2 g paste in a syringe, 20 x disposable tips



ADSEAL

- epoxy resin based root canal sealer
- easy to mix (paste-paste)
- hermetic sealing ability
- non-staining to the teeth
- insoluble in tissue fluids
- good radiopacity and apical tolerance
- packaging: base 9 g, catalyst 4.5 g, spatula, mixing plate

Z.O.B. Seal

- zinc oxide based sealer
- permanent root canal sealing
- good adhesion to tooth material
- good mixing consistency
- excellent radiopacity and anti-bacterial effect
- packaging: 15 g zinc oxide powder, 7 g liquid (Eugenol), spoon, mixing pad



MD-ChelCream

- 19 % EDTA cream
- root canal preparation ahead of filling with gutta percha and cement
- for root canal cleaning and smear layer removal
- removal of inorganic substances by chelation (EDTA-CA^{2+})
- packaging: 2 x 7 g syringe MD-ChelCream 19 % EDTA, 10 x disposable tips




MD-Cleanser

- 17 % EDTA solution
- for root canal cleaning and smear layer removal
- removal of inorganic substances by chelation (EDTA-CA^{2+})
- easier instrumentation by softening dentin canal walls
- packaging: 100 ml MD-cleanser solution, 17 % EDTA, 1 cc syringe, 3 x disposable tips

Endo Box, round

- metal box for burs or scaler tips
- dimensions: Ø 4.5 x H 5.5 cm
- sterilisable and autoclavable
- packaging: 1 box
- do not use thermodisinfectant



REF	material	colour	version
WSB 101	anodized metal	blue	for 9 burs, FG-shaft
WSB 102		red	for 9 burs, FG-shaft
WSB 103		blue	for 6 burs, RA-shaft
WSB 104		red	for 6 burs, RA-shaft
WSB 105		blue	for 6 scaler tips
WSB 106		red	for 6 scaler tips
REF	material	colour	dimensions
WSB 109	anodized metal	blue	Ø 4.5 x H 5.5 cm
WSB 110		red	Ø 4.5 x H 5.5 cm
REF	material	colour	dimensions
WSB 111	anodized metal	blue	Ø 4.5 x H 5.5 cm
WSB 112		red	Ø 4.5 x H 5.5 cm

Endo Box, round for 20 endodontic instruments

- metal box for 20 endodontic instruments
- do not use thermodisinfectant
- sterilisable and autoclavable
- packaging: 1 box



Endo Box, round with 6 colour-coded plastic cups

- metal box with colour-coded plastic cups
- sterilisable and autoclavable
- do not use thermodisinfectant
- packaging: 1 box



Endo Box, small

- rectangular metal box for burs and scaler tips
- sterilisable and autoclavable
- do not use thermodisinfectant
- dimensions: W 6.3 x D 5 x H 5.5 cm
- packaging: 1 box



REF	material	colour	version
WSB 201	anodized metal	blue	for 14 burs, FG-shaft
WSB 202		red	for 14 burs, FG-shaft
WSB 203		blue	for 6 burs, RA-shaft
WSB 204		red	for 6 burs, RA-shaft
WSB 205		blue	for 6 scaler tips
WSB 206		red	for 6 scaler tips

Endo Box, small for 30 endodontic instruments

- rectangular metal box for 30 endodontic instruments
- sterilisable and autoclavable
- do not use thermoisinfectant
- dimensions: W 6.3 x D 5 x H 5.5 cm
- packaging: 1 box



Endo Box, small with 6 colour-coded plastic cups

- rectangular metal box with colour-coded plastic cups
- sterilisable and autoclavable
- do not use thermoisinfectant
- packaging: 1 box

REF	material	colour	dimensions (W x D x H)
WSB 209	anodized metal	blue	6.3 x 5 x 5.5 cm
WSB 210		red	6.3 x 5 x 5.5 cm
REF	material	colour	dimensions (W x D x H)
WSB 211	anodized metal	blue	6.3 x 5 x 5.5 cm
WSB 212		red	6.3 x 5 x 5.5 cm

Endo Box, medium

- rectangular metal box for burs
- sterilisable and autoclavable
- do not use thermoisinfectant
- dimensions: W 10 x D 5.5 x H 5.5 cm
- packaging: 1 box



WSB 301
WSB 302

WSB 303
WSB 304

Endo Box, medium

- rectangular metal box for burs respectively scaler tips
- sterilisable and autoclavable
- do not use thermoisinfectant
- dimensions: W 10 x D 5.5 x H 5.5 cm
- packaging: 1 box



WSB 305
WSB 306

WSB 307
WSB 308

REF	material	colour	version
WSB 301	anodized metal	blue	for 36 burs, FG-shaft
WSB 302		red	for 36 burs, FG-shaft
WSB 303		blue	for 16 burs, RA-shaft
WSB 304		red	for 16 burs, RA-shaft
WSB 305	anodized metal	blue	for 16 scaler tips
WSB 306		red	for 16 scaler tips
WSB 307		blue	for 10 scaler tips + 18 burs, FG-shaft
WSB 308		red	for 10 scaler tips + 18 burs, FG-shaft



Endo Box, medium for 60 endodontic instruments

- rectangular metal box for 60 endodontic instruments
- do not use thermoisinfector
- sterilisable and autoclavable
- packaging: 1 box

REF	material	colour	dimensions (W x D x H)
WSB 309	anodized metal	blue	10 x 5.5 x 5.5 cm
WSB 310		red	10 x 5.5 x 5.5 cm
WSB 311	anodized metal	blue	10 x 5.5 x 5.5 cm
WSB 312		red	10 x 5.5 x 5.5 cm



Endo Box, medium with 18 colour-coded plastic cups

- rectangular metal box with colour-coded plastic cups
- do not use thermoisinfector
- sterilisable and autoclavable
- packaging: 1 box




WSB 314

Tip Stand Endo Box for 6 Scaler Tips

- stainless steel
- rotatable cover can be used as a stand
- incl. 6 silicone inserts

REF	material	colour	dimensions (W x D x H)
WSB 314	stainless steel	silver	8.3 x 2.5 x 4.5 cm



WSB 313

Plastic Bur Box

- for 16 scaler tips
- autoclavable
- packaging: 1 box

REF	material	colour	dimensions (W x D x H)
WSB 313	plastic	transparent/blue	10 x 5.5 x 5.5 cm



WEC 110



WEC 120

Endo-Cleaner

- useful for cleaning, sterilising, and staging endodontic instruments during treatment
- autoclavable up to 134 °C
- do not use thermoisinfector
- suitable foam inserts: for WEC 110: WFI 001 for WEC 120: WFI 002
- packaging: 1 piece

REF	material	colour	dimensions
WEC 110	aluminium	blue	Ø 40 mm, H: 40 mm
WEC 120		silver	Ø 40 mm, H: 30 mm



WSB 010

Bur Stand »Silver Bell«

- for 14 FG-burs and 7 angle-burs
- packaging: 1 stand

REF	material	colour
WSB 010	stainless steel	silver



DENTAL JEWELLERY



Brilident™
Zahnschmucksteine

Brilident Tooth Gem

• packaging: 1 gem

REF	colour	dimensions
B 101	crystal clear	Ø 1.8 mm
B 301		Ø 2.0 mm
B 401		Ø 2.6 mm
B 102	aquamarine	Ø 1.8 mm
B 302		Ø 2.0 mm
B 103		Ø 1.8 mm
B 303	light blue	Ø 2.0 mm
B 104		Ø 1.8 mm
B 304		Ø 2.0 mm
B 107	light green	Ø 1.8 mm
B 307		Ø 2.0 mm
B 109		Ø 1.8 mm
B 309	green	Ø 2.0 mm
B 111		Ø 1.8 mm
B 311		Ø 2.0 mm
B 112	yellow	Ø 1.8 mm
B 312		Ø 2.0 mm
B 114		Ø 1.8 mm
B 314	pink	Ø 2.0 mm
B 116		Ø 1.8 mm
B 316		Ø 2.0 mm
B 118	red	Ø 1.8 mm
B 318		Ø 2.0 mm
B 120		Ø 1.8 mm
B 320	pearl	Ø 2.0 mm
B 420		Ø 2.6 mm
B 121		Ø 1.8 mm
B 122	light-rose	Ø 1.8 mm
B 123	ocean-blue	Ø 1.8 mm
B 124	sandy-cool	Ø 1.8 mm
B 125	sunny-yellow	Ø 1.8 mm
B 125	fire-red	Ø 1.8 mm





Brilident DIN A2 Poster Butterfly

- DIN A 2 poster for dental practices
- beautiful brand advertisement and incentive for patients
- packaging: 1 poster


REF	material	dimensions	packaging
B 503	paper	DIN A2	1 poster



Brilident Set assorted

- all 12 colours
- packaging: 12 gems


REF	colour	dimensions	packaging
B 151	assorted colours	Ø 1.8 mm	12 pcs.
B 351	assorted colours	Ø 2.0 mm	12 pcs.



Brilident Set assorted

- all 12 colours + 10 x clear gem
- packaging: 22 gems

REF	colour	dimensions	packaging
B 152	assorted + clear	Ø 1.8 mm	22 pcs.
B 352	assorted + clear	Ø 2.0 mm	22 pcs.



Brilident Set assorted

- packaging: 6 gems

REF	colour	dimensions	packaging
B 161	reddish colours	Ø 1.8 mm	6 pcs.
B 361	reddish colours	Ø 2.0 mm	6 pcs.
B 162	bluish colours	Ø 1.8 mm	6 pcs.
B 362	bluish colours	Ø 2.0 mm	6 pcs.

Brilident™
Zahnschmucksteine

Brilident™
Zahnschmucksteine

Brilident™
Zahnschmucksteine



izi
jewellery



iziTooth Gem Star Sign

- tooth jewellery for natural teeth
- 22 ct yellow gold or 18 ct white gold
- incl. bonding guide
- packaging: 1 gem




REF	description	colour
22 YZ 12	izi Pisces	gold
18 WZ 12		white gold
22 YZ 06	izi Virgo	gold
18 WZ 06		white gold
22 YZ 04	izi Cancer	gold
18 WZ 04		white gold
22 YZ 05	izi Leo	gold
18 WZ 05		white gold
22 YZ 09	izi Sagittarius	gold
18 WZ 09		white gold
22 YZ 08	izi Scorpio	gold
18 WZ 08		white gold
22 YZ 10	izi Capricorn	gold
18 WZ 10		white gold
22 YZ 02	izi Taurus	gold
18 WZ 02		white gold
22 YZ 07	izi Libra	gold
18 WZ 07		white gold
22 YZ 11	izi Aquarius	gold
18 WZ 11		white gold
22 YZ 01	izi Aries	gold
18 WZ 01		white gold
22 YZ 03	izi Gemini	gold
18 WZ 03		white gold

izi
jewellery



iziTooth Gem

- tooth jewellery for natural teeth
- 22 ct yellow gold or 18 ct white gold
- incl. bonding guide
- packaging: 1 gem



REF	description	colour
22 Y 640	izi round/diamond 0.01 ct	gold
18 W 640		white gold
22 Y 645	izi round/black diamond 0.01 ct	gold
18 W 645		white gold
22 Y 646	izi round/black brilliant 0.02 ct	gold
18 W 646		white gold
22 Y 644	izi round/white brilliant 0.02 ct	gold
18 W 644		white gold
22 Y 642	izi round with ruby	gold
18 W 642		white gold
22 Y 643	izi round with sapphire	gold
18 W 643		white gold
22 Y 641	izi round with zirconia	gold
18 W 641		white gold
22 Y 506	izi heart, large	gold
18 W 506		white gold
22 Y 505	izi heart, small	gold
18 W 505		white gold
22 Y 620	izi heart/diamond 0.01 ct	gold
18 W 620		white gold
22 Y 625	izi heart/black diamond 0.01 ct	gold
18 W 625		white gold
22 Y 622	izi heart with ruby	gold
18 W 622		white gold
22 Y 623	izi heart with sapphire	gold
18 W 623		white gold
22 Y 621	izi heart with zirconia	gold
18 W 621		white gold
22 Y 507	izi moon	gold
18 W 507		white gold
22 Y 630	izi moon/diamond 0.01 ct	gold
18 W 630		white gold
22 Y 635	izi moon/black diamond 0.01 ct	gold
18 W 635		white gold
22 Y 632	izi moon with ruby	gold
18 W 632		white gold
22 Y 633	izi moon with sapphire	gold
18 W 633		white gold
22 Y 631	izi moon with zirconia	gold
18 W 631		white gold




izi Tooth Gem

- tooth jewellery for natural teeth
- 22 ct yellow gold or 18 ct white gold
- incl. bonding guide
- packaging: 1 gem



REF	description	colour
22 Y 502	izi star, large	gold
18 W 502		white gold
22 Y 501	izi star, small	gold
18 W 501		white gold
22 Y 600	izi star/diamond 0.01 ct	gold
18 W 600		white gold
22 Y 605	izi star/black diamond 0.01 ct	gold
18 W 605		white gold
22 Y 602	izi star with ruby	gold
18 W 602		white gold
22 Y 603	izi star with sapphire	gold
18 W 603		white gold
22 Y 601	izi star with zirconia	gold
18 W 601		white gold
22 Y 504	izi droplet, large	gold
18 W 504		white gold
22 Y 503	izi droplet, small	gold
18 W 503		white gold
22 Y 610	izi droplet/diamond 0.01 ct	gold
18 W 610		white gold
22 Y 615	izi droplet/black diamond 0.01 ct	gold
18 W 615		white gold
22 Y 613	izi droplet with ruby	gold
18 W 613		white gold
22 Y 612	izi droplet with sapphire	gold
18 W 612		white gold
22 Y 611	izi droplet with zirconia	gold
18 W 611		white gold
22 YG 12	izi skull	gold
18 WG 12		white gold
22 YG 04	izi tribal	gold
18 WG 04		white gold
22 YG 05	izi female	gold
18 WG 05		white gold
22 YG 15	izi Superman	gold
18 WG 15		white gold
22 YG 22	izi unicorn	gold
18 WG 22		white gold
22 Y 517	izi yin yan	gold
18 W 517		white gold
22 Y 519	izi cherries	gold
18 W 519		white gold
22 YG 10	izi rose	gold
18 WG 10		white gold
22 YG 24	izi springbok	gold
18 WG 24		white gold
22 Y 511	izi club	gold
18 W 511		white gold
22 Y 512	izi spade	gold
18 W 512		white gold
22 Y 520	izi @-sign	gold
18 W 520		white gold
22 Y 510	izi diamonds	gold
18 W 510		white gold






izi Tooth Gem

- tooth jewellery for natural teeth
- 22 ct yellow gold or 18 ct white gold
- incl. bonding guide
- packaging: 1 gem

REF	description	colour
22 Y 521	izi moon & star	gold
18 W 521		white gold
22 Y 515	izi alien	gold
18 W 515		white gold
22 YG 09	izi peace	gold
18 WG 09		white gold
22 Y 518	izi lips	gold
18 W 518		white gold
22 YG 19	izi maple leaf	gold
18 WG 19		white gold
22 YG 01	izi bunny	gold
18 WG 01		white gold
22 Y 509	izi dolphin	gold
18 W 509		white gold
22 Y 508	izi dollar sign	gold
18 W 508		white gold
22 YG 21	izi falcon	gold
18 WG 21		white gold
22 Y 516	izi soccerball	gold
18 W 516		white gold
22 YG 17	izi gecko	gold
18 WG 17		white gold
22 Y 522	izi tick	gold
18 W 522		white gold
22 YG 16	izi hemp leaf	gold
18 WG 16		white gold
22 YG 20	izi Harley	gold
18 WG 20		white gold
22 YG 06	izi horseshoe	gold
18 WG 06		white gold
22 YG 02	izi cat	gold
18 WG 02		white gold
22 YG 13	izi clover	gold
18 WG 13		white gold
22 Y 523	izi music note enamel	gold
18 W 523		white gold
22 YG 07	izi Mickey Mouse	gold
18 WG 07		white gold
22 YG 08	izi motorbike	gold
18 WG 08		white gold
22 YG 18	izi music note	gold
18 WG 18		white gold
22 YG 03	izi treble clef	gold
18 WG 03		white gold
22 YG 23	izi scissors	gold
18 WG 23		white gold
22 Y 513	izi butterfly	gold
18 W 513		white gold
22 YG 11	izi scorpion	gold
18 WG 11		white gold
22 Y 514	izi smiley	gold
18 W 514		white gold
22 YG 14	izi sun	gold
18 WG 14		white gold



wellsamed Etching Gel

- 36 % phosphoric acid etchant
- to etch tooth enamel and dentin
- blue colour is easily distinguishable
- packaging: 1 x 6.5 g (5 ml) syringe Etch Gel, 2 x straight disposable etchant tips, 2 x angled disposable etchant tips, 2 adapter



DENTAL RESTORATION

wellsamed Etching Gel Maxi

- 36 % phosphoric acid etchant
- to etch tooth enamel and dentin
- blue colour is easily detectable
- packaging: 1 x 62 g (50 ml) syringe Etchant Gel, 10 x 1.9 g delivery syringes, 50 x angled disposable etchant tips, 2 x adapters



REF	description	packaging
WET 005	Etching Gel	1 x 6.5 g (5 ml) syringe Etch Gel, 2 x straight disposable etchant tips, 2 x angled disposable etchant tips, 2 adapter
WET 015	Etching Gel Maxi	1 x 62 g (50 ml) syringe Etchant Gel, 10 x 1.9 g delivery syringe, 50 x angled disposable etchant tips, 2 x adapter



WMA 020



WAH 251



WAH 252



WAH 253



WAH 254



WMA 118 WMA 113 WMA 112 WMA 116



WMA 210 WMA 218 WMA 213 WMA 212 WMA 216

Wellsatip Starter Set

- assortment of disposable micro applicators
- convenient set for first-time users
- made of plastic, flocked
- assorted colours
- packaging: 1 box of 400 pcs., 2 holders

Instantly ready
for use:
all Applitip colours
incl. 2 handles

REF	material	colour	size	packaging
WMA 020	plastic, flocked	assorted colours	short	400 pcs. + 2 holders

Wellsamed Applicator Brush Handle

- plastic, double ended
- for short, disposable micro applicators
- can be assembled on both sides
- ergonomically designed handle for precise application
- packaging: 1 bag of 2 pcs.

REF	material	colour	packaging
WAH 251	plastic	orange	2 pcs.
WAH 252		blue	2 pcs.
WAH 253		white	2 pcs.
WAH 254		red	2 pcs.

Wellsamed Applitip Short

- disposable micro applicators, short
- made of plastic, flocked
- precise application of small, controlled amounts of material in limited access areas
- packaging: bag of 100 pcs.

REF	material	colour	packaging
WMA 118	plastic, flocked, short	red	100 pcs.
WMA 113		green	100 pcs.
WMA 112		yellow	100 pcs.
WMA 116		blue	100 pcs.

Wellsamed Applitip Long

- disposable micro applicators, long
- made of plastic, flocked
- improved control and accurate placement of solution in hard-to-reach areas
- packaging: 1 can of 100 pcs.



REF	description	colour	packaging
WMA 210	Applip Long	black	100 pcs.
WMA 218		red	100 pcs.
WMA 213		green	100 pcs.
WMA 212		yellow	100 pcs.
WMA 216		blue	100 pcs.

wellsamed Disposable Micro Applicators

- standard handles
- precise and consistent application
- length: 100 mm
- packaging: 1 can of 100 pcs.



REF	material	colour	size	packaging
WMA 233	plastic, flocked, long	blue	medium	100 pcs.
WMA 232		green	fine	100 pcs.
WMA 231		purple	ultrafine	100 pcs.
WMA 241		white	cylinder	100 pcs.



wellsamed Disposable Composite Brush Premium

- nylon hair
- material: plastic
- packaging: 100 pcs.

**TOP
seller**

REF	material	colour	dimensions	packaging
WAB 020	plastic	black	medium (short)	100 pcs.
WAB 022		black		100 pcs.
WAB 010		black	x-fine (short)	100 pcs.
WAB 011		white		100 pcs.



**Instantly ready
for use:**
all Brushtip
Applicator Brushes
incl. 1 handle

wellsamed Brushtip Applicator Brush

- disposable composite brushes
- brush tips with flat or tapered trim and nylon bristles
- transparent/red (WAB 182) shaft
- for precise and easy application
- material: plastic
- packaging: 50 or 100 pcs. + 1 handle

REF	material	colour	dimensions	packaging
WAB 102	plastic	black	x-fine (short, tapered)	100 pcs.
WAB 100		transparent		100 pcs.
WAB 152		black	medium (short, tapered)	100 pcs.
WAB 150		transparent		100 pcs.
WAB 182		black	medium (short, flat)	100 pcs.
WAB 157		black	medium (short, flat)	100 pcs.
WAB 352		black	medium (long, tapered)	50 pcs.
WAB 357		black	medium (long, flat)	50 pcs.



wellsamed Brush Holder

- for composite brushes
- single-ended plus spatula
- material: plastic
- packaging: 1 holder




WAH 164

WAH 165

wellsamed Brushtip Holder

- single-ended or double-ended holder for brushes or applicators
- material: plastic
- packaging: 1 holder



WAH 171

WAH 173

WAH 182


wellsamed Brushtip Holder Eco

- long holder for brushtip Applicator Brushes
- material: plastic
- packaging: 1 bag of 10 pcs.



WAH 110

REF	description	colour	packaging
WAH 165	Brush Holder plus spatula	white	1 piece
WAH 164		black	1 piece
WAH 171	Brushtip Holder, double-ended	orange	1 piece
WAH 173		white	1 piece
WAH 182	Brushtip Holder, single-end	blue	1 piece
WAH 110	Brushtip Holder Eco	white	10 pcs.



WAB 218

WAB 212

WAB 213

WAB 216

wellsamed Brushtip Long

- disposable composite brushes
- long and easily bendable
- for precise application
- length: 103 mm
- material: plastic/nylon hair
- packaging: 1 can of 100 pcs.



REF	material	colour	dimensions	packaging
WAB 218	plastic/nylon hair	red	103 mm (long)	100 pcs.
WAB 212		yellow	103 mm (long)	100 pcs.
WAB 213		green	103 mm (long)	100 pcs.
WAB 216		blue	103 mm (long)	100 pcs.

wellsamed Disposable Applicator Tips

- bent needle tips for use in placing viscous to thin materials, e.g. etchants
- non-sterile
- disposable
- luer lock type
- fits all standard injection plungers
- packaging: 100 pcs.

REF	description	colour	dimensions	packaging
WAC 352	18G	pink	1.2 mm x 15 mm	100 pcs.
WAC 353	20G	yellow	0.9 mm x 15 mm	100 pcs.
WAC 355	22G	black	0.7 mm x 15 mm	100 pcs.
WAC 354	20G	black opaque	0.9 mm x 15 mm	100 pcs.
WAC 357	25G	blue	0.5 mm x 15 mm	100 pcs.




WAC 352

WAC 353

WAC 355

WAC 354

WAC 357



Amalgam Gun Straight and with 45° Tip

- autoclavable up to 134 °C
- used to transfer the mixed amalgam into the prepared cavity
- screw-on tips with Ø 2.5 mm
- packaging: 1 gun incl. 1 straight/45° tip

REF	description	colour	packaging
WAP 001	straight	white	1 piece
WAP 002	with 45° Tip	white	1 piece



Disposable Mixing Tray

- made of plastic
- 12 wells per tray
- packaging: 50 pcs.

REF	material	colour	packaging
WMB 301	plastic	white	50 pcs.



Disposable Spatulas

- made of plastic
- very flexible
- composite materials do not stick to the spatula
- packaging: 100 pcs.

REF	material	colour	packaging
WMS 111	plastic, flexible	white	100 pcs.

Disposable Dappen Dishes

- to mix or supply liquids, medicaments or pastes
- material: plastic
- include finger ring
- not sterilisable, single use only
- packaging: 1 pack of 50 pcs.

REF	material	colour	packaging
WDD 100	plastic	black	50 pcs.
WDD 110		white	50 pcs.
WDD 120		yellow	50 pcs.
WDD 160		blue	50 pcs.
WDD 180		red	50 pcs.





WDG 100



WDG 135



WDG 195



WDG 190



WDG 130



WDG 150



WDG 160


Glass Dappen Dish

- to mix or supply liquids, medicaments or pastes
- coloured glass
- suitable for thermal disinfection or sterilisation
- packaging: 1 dish

REF	material	colour	packaging
WDG 100	glass	clear	1 piece
WDG 130		light green	1 piece
WDG 135		dark green	1 piece
WDG 150		light blue	1 piece
WDG 160		dark blue	1 piece
WDG 190		brown	1 piece
WDG 195		dark brown	1 piece

Disposable Pipettes

- 2 compartments
- drop-wise application or irrigation
- packaging: 1 pack of 100 pcs.



WMD 010

REF	material	packaging
WMD 010	plastic	100 pcs.

Silicone Mixing Bowls

- used for mixing dental acrylics, resins, and cements
- suction cup on bottom prevents sliding
- clear and durable
- packaging: 1 bowl

WMB 120



WMB 121



WMB 122



WMB 123

REF	material	colour	dimensions	packaging
WMB 120	silicone	transparent	XS/Ø 2.5 x H 2.5 cm	1 bowl
WMB 121		transparent	S/Ø 3.5 x H 3 cm	1 bowl
WMB 122		transparent	M/Ø 4.5 x H 4.5 cm	1 bowl
WMB 123		transparent	L/Ø 6.5 x H 5.5 cm	1 bowl



Biner LC

- light-curing, radiopaque cavity liner
- contains hydroxy calcium phosphate
- base material for use with adhesives and composites and with conventional restorative materials
- protects dental pulp
- lining under all filling materials
- packaging: 2 g syringe, 10 disposable tips

NETC

- non-eugenol, temporary cementation of crowns and bridges
- easy to mix by paste-paste type
- 1.5 minutes working time and 6 minutes setting time
- packaging: Bulk Pack (135293): 60 g base (zinc oxide), 20 g catalyst, mixing pad
Standard Pack (135294): 40 g base (zinc oxide), 14 g catalyst, mixing pad

REF	packaging
135294	Standard Pack
135293	Bulk Pack



Metacem

- high strength, permanent dual cure resin cement
- crowns and bridges (ceramic, composite, porcelain)
- inlay/onlay cementation
- bonding for porcelain veneers
- endodontic posts cementation
- core build-up
- high resistance to compressive force
- easy to use for cosmetic, indirect restorations
- packaging: Metacem Set: 2 x 9 g syringe (A3, translucent white), 5 g adhesive, 3 g etchant, 20 x mixing tips, 20 x etchant disposable tips
Metacem Refill: 1 x 9 g syringe

REF	colour	packaging
135295	A3, translucent white	Set
135296	A1	Refill syringe
135297	A3	
135298	B2	
135299	translucent white	
135300	opaque	

MD-Temp white/pink

- ready-for-use, self-curing, temporary cavity filling based on zinc oxide-polyvinyl acetate
- easy application and removal
- hermetic sealing of the cavity
- excellent stability when set and proper crush strength
- packaging: 40 g white or pink paste in a jar

REF	colour	packaging
135304	white	40 g
135305	pink	40 g



Meta Etchant

- 37 % phosphoric acid etchant
- conditioning of surface for dentine and enamel
- applicable with composite and inlay/onlay
- easy dispensing and application, no residue
- thixotropic characteristics (not runny)
- blue colouring for better recognisability
- packaging:
Kit: 3 x 3g syringe, 20 x disposable tips
Refill: 3 g syringe

REF	packaging
135301	Kit
135302	Refill Syringe



Meta P & Bond

- 5th generation, one step for priming and bonding
- bonding for direct composite, amalgam, veneers and posts
- for indirect bonding procedures with Metacem
- superior bonding strength
- ethanol base adhesive
- packaging: 5 g bottle, 10 x micro brushes





135306



135307



135310

Nexcomp

- aesthetic, light-curing nano-hybrid composite for anterior and posterior restorations
- composite and porcelain repairs, core build-up
- for diastema closures and direct veneers
- low 1.6 % shrinkage
- excellent polishability and sculptability
- packaging:
Kit: 7 x 4 g syringe (A1, A2, A3, A3.5, B2, A2E, B2E), 3 g Meta Etchant, 5 g Meta B&Bond, 20 x etchant disposable tips
Refill: 4 g syringe

REF	colour	type	packaging
135306	A1, A2, A3, A3.5, B2, A2E, B2E	Kit	7 x 4 g syringe, 3 g Meta Etchant, 20 x disposable tips
135307	A1	Refill syringe	4 g
135308	A2		4 g
135309	A3		4 g
135310	A3.5		4 g
135311	A4		4 g
135316	B1		4 g
135317	B2		4 g
135318	B3		4 g
135320	A2E		4 g
135321	B2E		4 g
135322	C2		4 g
135323	C3		4 g
135324	D2		4 g
135325	D3		4 g
135326	TL		4 g
135327	WT		4 g
135328	OP		4 g
135329	OA2		4 g
135330	OA3		4 g

Nexcomp

- aesthetic, light-cured nano-hybrid composite for anterior and posterior restorations
- composite and porcelain repairs, core build-up
- for diastema closures and direct veneers
- low 1.6 % shrinkage
- excellent polishability and sculptability
- packaging: 20 x 0.25 g



135333

REF	colour	type	packaging
135331	A1	Refill Unidose	20 x 0.25 g
135332	A2		20 x 0.25 g
135333	A3		20 x 0.25 g
135334	A3.5		20 x 0.25 g
135335	A4		20 x 0.25 g
135340	B1		20 x 0.25 g
135341	B2		20 x 0.25 g
135342	B3		20 x 0.25 g
135344	A2E		20 x 0.25 g
135345	B2E		20 x 0.25 g
135346	C2		20 x 0.25 g
135347	C3		20 x 0.25 g
135348	D2		20 x 0.25 g
135349	D3		20 x 0.25 g
135350	TL		20 x 0.25 g
135351	WT		20 x 0.25 g
135352	OP		20 x 0.25 g
135353	OA2		20 x 0.25 g
135354	OA3		20 x 0.25 g



Nexcomp Flow

- aesthetic, flowable nano-hybrid composite for minimally invasive cavity preparation
- pits and fissure sealing
- low viscosity
- light-cured, radiopaque
- excellent polishability
- packaging:
 Kit (135355): 5 x 2 g syringe (A1, A2, A3.5, B2, B3), 20 x disposable tips
 Refill (135356-135369): 2 x 2 g syringe, 10 x disposable tips



REF	colour	type	packaging
135355	A1, A2, A3.5, B2, B3	Kit	5 x 2 g syringe, 20 x disposable tips
135356	A1	Refill syringe	2 x 2 g syringe + 10 disposable tips
135357	A2		2 x 2 g syringe + 10 disposable tips
135358	A3		2 x 2 g syringe + 10 disposable tips
135359	A3.5		2 x 2 g syringe + 10 disposable tips
135360	A4		2 x 2 g syringe + 10 disposable tips
135361	B1		2 x 2 g syringe + 10 disposable tips
135362	B2		2 x 2 g syringe + 10 disposable tips
135363	B3		2 x 2 g syringe + 10 disposable tips
135364	OA2		2 x 2 g syringe + 10 disposable tips
135365	OA3		2 x 2 g syringe + 10 disposable tips
135366	C2		2 x 2 g syringe + 10 disposable tips
135367	D2		2 x 2 g syringe + 10 disposable tips
135368	WT		2 x 2 g syringe + 10 disposable tips
135369	TL		2 x 2 g syringe + 10 disposable tips



INFECTION CONTROL

wellsaclean Wet Wipes Dispenser Starter Set

- round dispenser bucket for a max. of 2.5 litres
- wipes: grammage 60 g/m²
- unpunched
- wipes: non-woven material, bucket: 100 % PET
- packaging: 1 dispenser bucket, 6 rolls of 30 wipes



wellsaclean Wet Wipes Dispenser Wall Bracket Stainless Steel

- suitable for WCW 106

REF	material	dimensions	packaging
WCW 106	100% PET	max. 2.5 litres	1 bucket, 6 rolls of 30 wipes
WCW 206		30 x 30 cm	6 rolls of 30 wipes

REF	material	packaging
WCW 106-1	Stainless steel	1 piece

**TOP
seller**



WCC 113

WCC 106



WCW 105



NEW



WCW 305

NEW



WCW 302

**Only while
stock lasts!**



WCW 301/001

wellsaclean S Rapid Disinfectant

- ready-for-use rapid disinfectant for wiping and spraying
- aldehyde-free, on ethanol basis
- for rapid disinfection and cleaning of alcohol-resistant, non-invasive medical inventory and surfaces
- packaging: canister of 10 litres or bottle of 1 litre

Contact Time:

wipe disinfection bactericidal incl. MRSA, yeasticidal (minimal contamination): 1 min., spray disinfection bactericidal incl. MRSA, yeasticidal (high contamination): 2 min., tuberculocidal (M. terrae): 2 min., inactivation of enveloped viruses incl. Vaccinie, HBV, HIV, HCV, BVDV: 1 min.

REF	scent	packaging
WCC 106	Neutral	10 litres
WCC 116	Lemon	10 litres
WCC 126	Peach	10 litres
WCC 113	Lemon	1 litre

wellsaclean Refillable Dispenser Bucket

- 1 bucket for max. 5 litres respectively 90 wipes, rectangular
- resealable dispenser bucket
- dry wipes to soak with disinfectant or cleaning agent
- suitable for all standard wipes
- packaging: 1 bucket

REF	material	dimensions	packaging
WCW 105	plastic	max. 5 litres	1 piece

wellsaclean PET Wipes

- 100 % PET, dry, smooth, 50 g/m², to be soaked with wellsaclean S Rapid Disinfectant
- sheet dimension: 32 x 30 cm, roll height: ca. 16 cm
- recommended amount of impregnating liquids for 90 wipes a minimum 2.5 litres
- non-woven
- packaging: 3 rolls of 90 wipes

wellsaclean PET Wipes

- 100 % PET, dry, smooth, 50 g/m², to be soaked with wellsaclean S Rapid Disinfectant
- sheet dimension: 25 x 25 cm, roll height: ca. 16 cm
- recommended amount of impregnating liquids for 70 wipes a minimum 1.5 litres
- non-woven
- packaging: 3 rolls of 70 wipes

wellsaclean Non-Woven Wipes

- for disinfection-set
- dry wipes, 30 g/m²
- to be soaked with spray disinfectant for buckets from 3 to 5 litres
- wipe dimension: 25 x 25 cm
- non-woven
- packaging: 3 rolls of 150 wipes

REF	packaging
WCW 302	3 rolls of 70 wipes
WCW 305	3 rolls of 90 wipes
WCW 301/001	3 rolls of 150 wipes

TOP seller



WCC 204



wellsaclean Disinfection Wipes Refill

- ready-to-use disinfection wipes, aldehyde-free
- ethanol-based
- dermatologically tested
- impregnating liquids VAH-tested
- packaging: 1 refill package of 120 pcs.

Application and contact times:

Wipe surfaces and objects thoroughly and allow disinfectant to act. Do not use for cleaning of alcohol-sensitive surfaces. Not suitable for acrylic glass. Bactericidal incl. MRSA, tuberculocidal, yeasticidal acc. to VAH: conc. 2 min., limited virucidal incl. HBV, HIV, HCV: conc. 1 min.

Formulation
VAH-
tested



wellsaclean Disinfection Wipes Empty Dispenser

- packaging: 1 empty dispenser for wellsaclean Disinfection Wipes Refill

REF	scent	packaging
WCC 204	Neutral	120 wipes
WCC 214	Lemon	120 wipes
WCC 224	Peach	120 wipes

REF	scent	packaging
WCC 201	Neutral	1 piece
WCC 211	Lemon	1 piece
WCC 221	Peach	1 piece

Ready-to-use:

Disinfection Wipes, impregnated with wellsaclean S Rapid Disinfectant



WCW 407

NEW

Formulation
VAH-
tested

wellsaclean Disinfection Wipes

- wet wipes made of PET
- contains disinfectant wellsaclean S neutral
- aldehyde-free
- fast acting
- open bag and put the wet wipes in the according dispenser bucket (WCW 105, pg. 70)
- ready-to-use

REF	description		packaging
WCW 418	Disinfection Wipes	152 x 300 mm	180 wipes
WCW 409		305 x 320 mm	90 wipes
WCW 407		250 x 250 mm	70 wipes

wellsaclean Orange Clean

- cleansing spray for the thorough removal of band-aid residue, residuals of adhesive (stickers, labels), impression materials and cements from laquer, plastic-coated surfaces and stainless-steel instruments
- removes zinc plaster, residuals of plaster casts and iodine stains

- easy and targeted application
- colourless and residue free
- good material compatibility
- fresh orange scent with natural orange oils
- packaging: aerosol can of 500 ml



WOS 201

REF	packaging
WOS 201	500 ml



WCC 303

WCC 306

wellsclean F Surface Disinfectant

- for both disinfection and cleaning of non-invasive medical devices
- alcohol-free and unscented
- aldehyde-free and phenol-free
- suitable for acrylic glass
- packaging: bottle of 1 litre or canister of 10 litres

Concentration and contact times:

bactericidal incl. MRSA, yeasticidal acc. to VAH: 7.5% (75 g/l) - 5 min., 2.5% (25 g/l) - 15 min., 2% (20 g/l) - 30 min., 1% (10 g/l) - 60 min., limited virucidal incl. HBV, HIV, HCV: 2% (20 g/l) - 15 min., 2.5% (25 g/l) - 5 min.

REF	packaging
WCC 303	1 litre
WCC 306	10 litres



WCC 503

WCC 506

wellsclean I Instrument Disinfectant

- highly effective, strong cleaning power
- suitable for ultrasonic baths
- aldehyde-free, phenol-free
- packaging: bottle of 1 litre or canister of 10 litres

Concentration and contact times:

bactericidal incl. MRSA, yeasticidal, tuberculocidal, limited virucidal (incl. HBV, HIV, HCV, BVDV, Vaccinie): 2% - 60 min., 4% - 30 min.

REF	packaging
WCC 503	1 litre
WCC 506	10 litres



WCC 705

wellsclean A Impression Disinfectant

- highly effective, ready-to-use solution disinfection of dental impression material
- formaldehyde-free
- contact time: 5 min. (incl. HBV/HCV/HIV)
- packaging: canister of 5 litres

REF	packaging
WCC 705	5 litres



WCC 605

wellsclean B Bur Bath

- ready-to-use universal solution for disinfecting and cleaning in ultrasonic baths and dental bur baths
- extensive microbiological effectiveness
- good corrosion inhibition
- packaging: canister of 5 litres

Contact times:

instrument disinfection acc. to VAH: conc. - 5 min., HBV/HIV/HCV inactivation: conc. - 5 min.

REF	packaging
WCC 605	5 litres



WCC 402

WCC 403



WCC 412

WCC 413



WCC 802

WCC 803



Dispenser wellisaclean

- suitable for liquid soaps, disinfectants etc.
- easy installation
- lockable
- transparent closure plate
- optionally plastic or stainless steel
- capacity of 500 ml (liquid soaps, disinfectants etc.)
- dimensions (W x D x H): 18,5 x 9,5 x 26 cm
- packaging: 1 empty bottle, wall plate, key, lockable closure plate



A108363

A108364

wellisaclean H Hand Disinfectant

- alcoholic solution for hygienic hand disinfection
- free of perfume
- packaging: bottle of 500 ml or 1 litre
- Ready-to-use solution to rub into dry

hands. During entire application time hands need to be kept moist by undiluted disinfectant. Acc. to EN 1500 hygienic hand disinfection: 3 ml - 30 seconds. Acc. to EN 12791: 2 x 5 ml - 1.5 minutes

REF	packaging
WCC 402	500 ml
WCC 403	1 litre

wellisaclean H^{plus} Hand Disinfectant

- alcoholic solution for hygienic hand disinfection
- free of perfume
- with refatting substances
- packaging: bottle of 500 ml or 1 litre

- Ready-to-use solution to rub into dry hands. During entire application time hands need to be kept moist by undiluted disinfectant. Acc. to EN 1500 hygienic hand disinfection: 3 ml - 30 seconds. Acc. to EN 12791: 2 x 5 ml - 1.5 minutes

REF	packaging
WCC 412	500 ml
WCC 413	1 litre

wellisacare W^{plus} Wash Lotion

- mild cleansing, skincare wash lotion for frequent use
- also suitable for sensitive skin

- with phanthenol
- packaging: bottle of 500 ml or 1 litre

REF	packaging
WCC 802	500 ml
WCC 803	1 litre

wellisamed Dispenser Euro

- specifically for medical practices, autoclavable
- for alcohol-containing hand disinfectants, liquid soaps, and skincare lotions
- aluminium casing with stainless steel pump, flexible plastic suction pipe
- short arm lever
- suitable for all European Standard bottles



A114725

A114726

REF	material	colour	dimensions
A114725	aluminium	silver / white	capacity: 500 ml European Standard bottles
A114726		silver / white	capacity: 1000 ml European Standard bottles
A108363	brushed stainless steel	silver	capacity: 500 ml wellisaclean bottles
A108364	plastic	white	capacity: 500 ml wellisaclean bottles



Instrument Bath

- capacity of 4 litres use-solution
- transparent lid for improved overview
- contains an inner removable drain tank and lid
- very sturdy
- packaging: 1 piece

Nail Brush Nylon

- scrub brush for hygienic cleaning of hands and nails
- plastic handle
- nylon brush
- assorted colours: white-grey, red, blue
- packaging: single brush



REF	material	dimensions	packaging
A104753	plastic	34 x 21 x 12 cm	1 piece

REF	material	colour	packaging
WNB 010	plastic/nylon	assorted colours	1 piece



Dosing Aid Bottle

- material: plastic



Dosing Aid Canister

- material: plastic

REF	description	dimensions	packaging
A104710	spray head with adapter	for 1 litre bottle	1 piece
A104711	spray head foam	for 1 litre bottle	1 piece
A104707	dosing pump	500 ml/1 litre bottle	1 piece
A104708	dosage bottle, empty	1 litre bottle with 60 ml metering chamber	1 piece
A102494	canister tap universal	for 5 litres/10 litres canister	1 piece
A104709	transfer pump	for 5 litres/10 litres canister	1 piece
A104712	canister key	for 5 litres/10 litres canister	1 piece


Universal Filter for Autoclave

- suitable for different autoclaves
- dimensions: Ø 53 mm
- packaging: 1 filter



Helix Test

- helix tube can be used up to 100 times
- single-use indicator strips with adhesive reverse side for documentation
- indicators acc. to ISO 11140-1 (Class 6)
- PCD conform EN 867-5
- 3.5 min 134 °C/15 min 121 °C
- packaging: 100 indicator strips, 1 helix tube



REF	material	colour	dimensions	packaging
A105040	plastic/filter paper	white	Ø 53 mm	1 piece

REF	packaging
WHT 100	1 helix tube + 100 indicator strips

Irrigation Tubing Set Standard

- sterile, disposable cooling system for surgical units
- WIR 101 suitable for: Aseptico AEU-7000, Kyocera Implantor Neo, NSK Osseocision, NSK Surgic XT/XT Plus, NSK VarioSurg
- WIR 102 suitable for: Nobel Biocare 80, W&H Drilltech, W&H DU900, W&H Frios Unit S//Frios Unit S, W&H Implantmed New Design 2005, W&H Implantmed
- WIR 103 suitable for: GeruDent Futura 1, Nouvag 3000, Nouvag 7000, Nouvag CBM, Nouvag MD10, Nouvag MD20 without tap
- packaging: 1 piece

NEW



WIR 102



WIR 103



WIR 101

REF	material	colour	packaging
WIR 101	plastic	white	1 piece
WIR 102		white	1 piece
WIR 103		white	1 piece

wellsaspray Universal Oil Spray

- for interior maintenance of dental turbines, handpcs., and contra-angles
- routine lubrication extends durability
- easy, safe, and clean application
- scentless, physiologically compatible
- can equipped with 360° valve, can be used overhead
- packaging: aerosol can of 500 ml



WOS 101

wellsamed Universal Oil Spray Adapter

- universally attachable to handpcs. and contra-angles
- high quality metal design with seal ring
- packaging: 1 adapter



WOS 992



WOS 993



WOS 994



WOS 991

REF	description	packaging
WOS 101	Oil Spray	500 ml
REF	description	
WOS 991	Adapter A	for handpcs. and contra-angles
WOS 992	Adapter B	for Roto Quick Turbine Sirona and W+H
WOS 993	Adapter C	for KaVo-turbines
WOS 994	Adapter D	for contra-angle heads



wellsamed Sterilisation Reel, 200 m

- high quality sterilisation reel constructed from paper and film
- three independent seal lines
- medical grade paper, enhanced grammage 60 g/m²
- accurate and clear colour change after sterilisation process by steam and gas
- acc. to DIN 58953-7 and EN 868
- packaging: roll of 200 m

REF	material	width	packaging
WSR 201	paper/ film	50 mm	200 m
WSR 202		75 mm	200 m
WSR 203		100 mm	200 m
WSR 204		150 mm	200 m
WSR 205		200 mm	200 m
WSR 206		250 mm	200 m
WSR 207		300 mm	200 m



wellsamed Sterilisation Reel, 100 m

- high quality sterilisation reel constructed from paper and film
- for packaging of small norm-trays and sterilisation cassettes
- three independent seal lines
- medical grade paper, enhanced grammage 60 g/m²
- double-layer polypropylene/polyester film
- accurate and clear colour change after sterilisation process by steam and gas
- acc. to DIN 58953-7 and EN 868
- packaging: roll of 100 m

REF	material	dimensions (W x L)	packaging
WSR 304	paper/ film	55 x 150 mm	100 m
WSR 305		50 x 200 mm	100 m
WSR 306		65 x 250 mm	100 m
WSR 307		80 x 300 mm	100 m

wellsamed Sterilisation Pouch

- top layer: film PET/PP, bottom layer: medical paper
- at the end an additional, fourth seal line
- medical grade paper, enhanced grammage 60 g/m²
- accurate and clear colour change after sterilisation process by steam and gas
- packaging: pack of 2,000 or 1,000 pcs.



wellsamed Self-Seal Sterilisation Pouch

- top layer: film PET/PP, bottom layer: medical paper, self-sealing
- accurate and fast seal without the requirement of a heat sealing machine
- accurate and clear colour change after sterilisation process
- packaging: pack of 500 or 200 pcs.



REF	material	dimensions	packaging
WSR 102	paper/film	75 x 250 mm	2000 pcs.
WSR 103		100 x 270 mm	2000 pcs.
WSR 106		200 x 300 mm	1000 pcs.
WSR 401		55 x 250 mm	500 pcs.
WSR 402		75 x 250 mm	200 pcs.
WSR 403		100 x 250 mm	200 pcs.

Steri-Tape with Indicators

- autoclave tape with indicator
- packaging: roll of 50 m



REF	width	packaging
WST 052	12 mm	50 m
WST 059	19 mm	50 m
WST 055	25 mm	50 m

wellsamed Non-Woven Sterilisation Wrap

- suitable for packaging cassettes
- after opening, the non-woven sterilisation wrap can be used as a sterile underlay
- suitable for sterilisation by steam, gas or plasma
- water-repellent
- soft, tear-resistant
- after sterilisation fit for storage up to 6 months
- WSR 511: 500 pcs.
- WSR 512: 400 pcs.



REF	material	colour	dimensions (W x L)	packaging
WSR 511	non-woven	blue	60 x 60 cm, 47 g	500 pcs.
WSR 512		green	60 x 60 cm, 53 g	400 pcs.

NEW



DEVICES & PRACTICE EQUIPMENT

LED Curing Light Standard

- three working modes: full, ramp or pulse
- radiant intensity: < 1200 mW/cm²
- range of wavelength: 420-480 nm
- modes of working time range from 5 sec to 30 sec
- handpiece dimension: L 24 cm x Ø 2.7 cm, 160 g
- autoclavable light guide: 8 mm
- incl. 100 protective light sleeves
- packaging: 1 curing light incl. accessories

REF	colour	dimensions
A129398#C01- Blue	blue	24 cm x 2.7 cm
A129404#C01-Lavender	lavender	24 cm x 2.7 cm
A129406#C01-Green	green	24 cm x 2.7 cm
A129407#C01-Lime	hellgreen	24 cm x 2.7 cm
A129411#C01-Yellow	yellow	24 cm x 2.7 cm
A129412#C01-White	white	24 cm x 2.7 cm
A129413#C01-Pink	pink	24 cm x 2.7 cm

NEW

ARTEC Surgical Suction Unit VarioSurg 50

- high suction output
- quiet in operation
- minimum amount of vibration
- continuous suction power adjustment
- suction power: max. 50 l/min
- no maintenance
- 6 year warranty
- also available with case tray
- packaging: 1 suction unit



132398

REF	description	dimensions (W x D x H)
132377	ARTEC Surgical Suction Unit	480 x 460 x 900 mm
132398	Case Tray	-/-



132377

Curing Light Ecolux H-100

- microprocessor-controlled, powerful, and with precise adaption
- radiant intensity: 500 mW/m²
- output wavelength range: 380 nm - 540 nm
- time set: 1 - 99 seconds
- protection fuse to prevent overheating
- halogen reflector bulb: 12 V/75 W
- power input: 120 W
- autoclavable light guide: Ø 8 mm
- dimensions point: W 16 x D 13 x H 14 cm
- dimensions handpiece: W 14 x D 12 x H 4 cm
- length power cord: 225 cm
- length cable connection: 217 cm
- packaging: 1 curing light

REF	colour	dimensions (W x D x H)
A101998	white	16 x 13 x 14 cm



A101998

Curing Light Blue Life LED-320

- cordless, powerful curing light
- working modes: full power, soft-up and impulse
- detailed LED display of time, mode, light intensity, and battery charge
- 2 lithium-ion batteries for nonstop operation
- fast curing (within 10 seconds)
- light intensity: max. 1.600 nW/cm² +/- 100
- wavelength: 450 - 490 nm
- lithium-ion battery (1350 mAh)
- light guide: Ø 8 mm, turbo light guide
- weight handpiece: 190 g
- diameter handpiece: 36 - 54 mm
- length handpiece: 190 mm
- power supply charger: 90 - 240 V/50 - 60 Hz
- power input charger: 15 W
- dimensions charger: W 64 x D 130 x H 55 mm
- packaging: handpiece, light guide, charger, 2 lithium-ion batteries, power cord, and light shield

REF	colour	dimensions (W x D x H)
A102003	blue / white	64 x 130 x 55 mm
A102007	grey / white	64 x 130 x 55 mm



A102003




A102007

Magnifying Lamp

- magnifying lens
- shadow-free illumination
- with ring lightening
- 22 W fluorescent lamp
- base: 6.2 x 6.2 cm (black)
- swinging arm: 80 cm (extended)
- Ø head: 22 cm
- magnifying display, top: 12 cm
- mounting rack: Ø 4.5 cm (inside)
- cable length: 1.80 m
- packaging: 1 lamp

REF	magnification
A104019	75 % (3 dpt.)
A104020	125 % (5 dpt.)

A104019
A104020



A108218

Straight Handpiece

- 1:1 transmission
- internal coolant tube
- maximum speed: 40,000 rpm
- packaging: 1 handpiece

REF	material	packaging
A108218	Stainless steel	1 piece

ARTEC Quick Coupling

- packaging: 1 piece



A121335

REF	compatible with	packaging
A121335	W&H turbine	1 piece



A108217

ARTEC LED Multiflex Coupling with Water Regulation

- with LED light
- packaging: 1 piece

REF	compatible with	packaging
A108217	KaVo	1 piece

ARTEC LED Quick Coupling with Water Regulation

- with LED light
- packaging: 1 piece



A121327

REF	compatible with	packaging
A121327	Sirona	1 piece

ARTEC LED-Bulb for Multiflex/Quick Coupling

- packaging: 1 piece

REF	compatible with	packaging
A108425	ARTEC Multiflex Coupling without water regulation (A107742)	1 piece
A108219	ARTEC Multiflex Coupling with water regulation (A108217)	1 piece
A121333	ARTEC Quick Coupling with water regulation (A121327)	1 piece



A108425



A108219



A121333



V106776

V106777

V106778

V106779

V106780

V106781



V106782

V106783

V106784

V106785

V106786

V106787

LED-Light HPSL-Spot Daylight

- LED light spectrum: ca. 4,500 K
- sterilisable, insensitive to impact and vibration
- much longer durability, resilience, light yield and reliability
- for all standard turbines, motors, and couplings
- low power consumption – very low self-heating
- packaging: 1 light

REF	material	packaging
A106776	KaVo	1 piece
A106777	Sirona	1 piece
A106778	BienAir turbine	1 piece
A106779	BienAir motor	1 piece
A106780	NSK turbine	1 piece
A106781	NSK-coupling	1 piece

LED-Light HPSL-Spot New-White

- for all standard turbines, motors, and couplings
- packaging: 1 light

REF	material	packaging
A106782	KaVo	1 piece
A106783	Sirona	1 piece
A106784	BienAir turbine	1 piece
A106785	BienAir motor	1 piece
A106786	NSK turbine	1 piece
A106787	NSK-coupling	1 piece



Water Distiller Aquaboy

- purifies water by eliminating water contaminants
- with carbon filter
- powerful (1 litre per minute)
- automatic switch-off (overheat protection)
- protects autoclave against damage
- removes 90 % of contaminants
- power: 700 W
- volume water tank: 4 litres
- output: 1 litre per hour
- working voltage: 230 V, 50/60 Hz
- packaging: 1 Aquaboy + accessories

REF	packaging
A103750	1 device + accessories



A104051

Carbon Filter for Aquaboy

- compatible with Water Distiller Aquaboy (A103750)
- packaging: 10 pcs.

REF	packaging
A104051	10 pcs.



A104052

Cleaning Agent for Aquaboy

- compatible with Water Distiller Aquaboy (A103750)
- packaging: pot of 140 g

REF	packaging
A104052	140 g



A104966

Water Tank for Aquaboy

- water tank for Water Distiller Aquaboy (A103750)
- packaging: 1 piece

REF	capacity
A104966	4 litres



122615

Genesys Complete Kit

- cordless, motorised device for root canal filling
- high capacity lithium-ion battery
- 360° rotation at both gutta percha needle and pack tip
- gutta percha needle refillable with gutta percha bars
- simple, easy, and quick to use
- packaging: 1 Kit: 1 x Pack Set (plugger F & FM), 1 Fill Set (gutta percha needle 23 G & 25 G), cleaning and bending tools, 3 x batteries, accessory gutta percha bars

Accessories		
REF	description	packaging
122616	Genesys Plugger F & FM	2 x Plugger (F & FM)
122617	Genesys Plugger M & ML	2 x Plugger (M & ML)
122618	Genesys Plugger XF	1 piece
122626	Genesys gutta percha needle	2 x GP needle (23G & 25G)
122627	Genesys gutta percha needle	2 x GP needle 20G
122628	Genesys gutta percha needle	2 x GP needle 23G
122629	Genesys gutta percha needle	2 x GP needle 25G
122630	Genesys gutta percha bars	100 bars



A101999

REF	dimensions (W x D x H)
A101999	20.5 x 25 x 28.5 cm

Magicmix F2 Alginate Mixer

- power: 220 W power
- mixing frequency: 3.400 rpm
- weight: 16 kg
- dimensions: (W x D x H) 20.5 x 25 x 28.5 cm
- short mixing time: 8 - 12 sec.
- 3 pre-programmed timer settings: 8, 10 and 12 sec.
- low noise level and safe due to automatic locking mechanism
- smooth and bubble-free mixture
- packaging: 1 alginate mixer + accessories



ADV 1011



ADV 1012



ADV 1014

Magicmix Accessories

- packaging: 1 piece

REF	description	colour	dimensions
ADV 1011	Mixing Cup	white	Ø 75 x H 65 mm
REF	description	colour	packaging
ADV 1012	Dosing Bottle	blue-transparent	1 piece
ADV 1013	Spatula	white	1 piece
ADV 1014	Basket	grey	1 piece

Capsule Mixer Mixboy 200

- low-noise and energy-saving capsule mixer for amalgams, cements, and other encapsulated dental materials
- controlled, consistent mix due to micro-processor control
- consumption: 200 W power
- mixing frequency: 4.350 oscillations per minute
- mixing time: adjustable from 1 to 99 seconds
- sturdy, ergonomic design, and easy-to-clean body
- better mix due to 8-shaped rotation
- weight: 3,5 kg
- packaging: 1 device

REF	colour	dimensions (W x D x H)	packaging
A101992	white	25 x 22 x 14 cm	1 piece



A101992



A101993

Capsule Mixer Mixboy 300

- low-noise and energy-saving capsule mixer for amalgams, cements, and other encapsulated dental materials
- controlled, consistent mix due to microprocessor control
- consumption: 200 W power
- mixing frequency: 4,250 oscillations per minute
- mixing time: adjustable from 1 to 99 seconds
- with memory function for time keeping
- better mix due to 8-shaped rotation
- weight: 2.3 kg
- packaging: 1 device

REF	colour	dimensions (W x D x H)	packaging
A101993	silver	24 x 20.5 x 15 cm	1 piece



A107116

A127491RAL1018
Example illustration RS4
with coloured handle



A107117

A129577RAL6027
Example illustration RS5
with coloured front

Cart RS4

- lockable, smooth-running, and low-noise castors
- 4 drawers fitted with pull-out trays and non-slip PVC-bottoms
- top shelf (W 500 x L 430 mm) with rolled edge and 8 mm tempered glass
- lower shelf (W 480 x L 260 mm) usable height of 31 cm
- dimensions: W 500 x D 440 x H 820 mm
- weight: 45 kg
- material: powder-coated stainless steel
- front and/or handle colours can be selected (body: white (RAL 9002))
- handles for front RAL 9002 white will be RAL 9011 black
- packaging: 1 cart

Cart RS5

- lockable, smooth-running, and low-noise castors
- 5 drawers fitted with pull-out trays and non-slip, sterilisable PVC-bottoms
- top shelf (W 500 x L 430 mm) with rolled edge and 8 mm demountable tempered glass
- dimensions: W 500 x D 440 x H 820 mm
- weight: 45 kg
- material: powder-coated stainless steel
- front and/or handle colours can be selected (body: white (RAL 9002))
- handles for front RAL 9002 white will be RAL 9011 black
- packaging: 1 cart

Please replace the X for your colour choice (RAL).

Cart RS4			
Coloured Front		Coloured Handle	
A120473 X	RAL 1018	A127491 X	RAL 1018
	RAL 1021		RAL 1021
	RAL 2003		RAL 2003
	RAL 3000		RAL 3000
	RAL 3003		RAL 3003
	RAL 5002		RAL 5002
	RAL 5010		RAL 5010
	RAL 5015		RAL 5015
	RAL 5022		RAL 5022
	RAL 6000		RAL 6000
	RAL 6018		RAL 6018
	RAL 6019		RAL 6019
	RAL 6027		RAL 6027
	RAL 7040		RAL 7040
	RAL 9005		RAL 9005
	RAL 9006		RAL 9006
A107116	RAL 9010		RAL 9010
	RAL 9016		RAL 9016
	RAL 9002		

Please replace the X for your colour choice (RAL).

Cart RS5			
Coloured Front		Coloured Handle	
A129577 X	RAL 1018	A129597 X	RAL 1018
	RAL 1021		RAL 1021
	RAL 2003		RAL 2003
	RAL 3000		RAL 3000
	RAL 3003		RAL 3003
	RAL 5002		RAL 5002
	RAL 5010		RAL 5010
	RAL 5015		RAL 5015
	RAL 5022		RAL 5022
	RAL 6000		RAL 6000
	RAL 6018		RAL 6018
	RAL 6019		RAL 6019
	RAL 6027		RAL 6027
	RAL 7040		RAL 7040
	RAL 9005		RAL 9005
	RAL 9006		RAL 9006
A107117	RAL 9010		RAL 9010
	RAL 9016		RAL 9016
	RAL 9002		

Additional haulage costs will incur. For further information please contact sales@wellsamed.de



A100897



A100895

Trolley Comfort with Drawer Extension

- ideal tool for different uses both in the operating theatre and in the dentist's office
- lacquered in white
- functional 3-shelf trolley
- with one drawer for accessories
- base frame with 4 easy-running wheels/only 2 front wheels lockable
- 2 universal sockets protected by 2 fuses
- dimensions drawer inside: W 268 x D 251 x H 64 mm
- dimensions drawer outside: W 340 x D 327 x H 80 mm
- working dimensions of top: W 385 x L 360 mm
- gap shelf 1 to 2: 265 mm
- gap shelf 2 to 3: 345 mm
- trolley dimensions: W 540 x D 400 x H 930 mm
- weight: ca. 25 kg
- packaging: 1 trolley

Instrument & Service Cart

- lacquered in white
- 4 drawers
- 2 sheet-steel shelves varnished with epoxy powder
- dimensions (W x D x H): 500 x 480 x 775 mm
- packaging: 1 cart

**TOP
seller**

REF	material	colour	dimensions (W x D x H)
A100897	powder-coated stainless steel	white	540 x 400 x 930 mm

REF	material	colour	dimensions (W x D x H)
A100895	powder-coated stainless steel	white	500 x 480 x 775 mm



A100901

Trolley Standard

- lacquered white
- 3-shelf trolley (W 400 x L 300 mm) with rounded corners
- power supply, 2 universal sockets
- base frame with 4 easy-running wheels
- dimensions (W x D x H): 450 x 400 x 800 mm
- packaging: 1 trolley



A100902

Trolley Standard Plus

- lacquered white
- internal water and air supply
- 3-shelf trolley (W 400 x L 300 mm) with rounded corners
- power supply (2 universal sockets), length power cord: 2 m
- base frame with 4 easy-running wheels
- incl. component hoses for water & air
- dimensions (W x D x H): 450 x 400 x 800 mm
- packaging: 1 trolley

REF	material	colour	dimensions (W x D x H)
A100901	powder-coated stainless steel	white	450 x 400 x 800 mm

REF	material	colour	dimensions (W x D x H)
A100902	powder-coated stainless steel	white	450 x 400 x 800 mm



A100898

**TOP
seller**

Trolley Eco

- trolley for trays and small equipment
- lacquered white
- 3 shelves with raised trim and rounded corners (W 400 x L 300 mm)
- base frame with 4 easy-running wheels
- dimensions (W x D x H): 490 x 370 x 770 mm
- packaging: 1 trolley



A105346

Trolley Implanto Light M

- 3 tempered glass shelves
- 4 lockable, smooth-running und low-noise castors
- with power supply (2 universal sockets)
- shelf dimensions: W 420 x D 361 x H 8 mm
- suitable for Piezotome 1, Intramax, and Implanleo
- packaging: 1 trolley

REF	material	colour	dimensions (W x D x H)
A100898	powder-coated stainless steel	white	490 x 370 x 770 mm

REF	material	colour	dimensions (W x D x H)
A105346	powder-coated stainless steel	white	480 x 400 x 790 mm

Additional haulage costs will incur. For further information please contact sales@wellsamed.de



A101509

**TOP
seller**

Trolley Universal with Drawer MAXI

- stainless steel (powder-coated)
- lacquered white
- sturdy chassis
- practical and handy design
- with 4 smooth-running castors
- 2 front drawers
- 2 sheet-steel shelves (W 550 x L 430 mm)
- dimensions (W x D x H): 620 x 460 x 820 mm
- packaging: 1 trolley

REF	material	colour	dimensions (W x D x H)
A101509	powder-coated stainless steel	white	620 x 460 x 820 mm

Trolley Universal with Drawer

- lacquered white
- sturdy chassis
- practical and handy design
- with 4 smooth-running castors
- 1 front drawer for small instruments
- 2 sheet-steel shelves (W 550 x L 430 mm)
- dimensions (W x D x H): 560 x 440 x 820 mm
- packaging: 1 trolley

REF	material	colour	dimensions (W x D x H)
A100904	powder-coated stainless steel	white	560 x 440 x 820 mm



A100904



A101775

Trolley F3S

- 3 height adjustable, powder-coated metal shelves (W 400 x L 270 mm)
- bearing strength: 10 kg each shelf
- lockable, low noise castors
- weight: 18 kg
- extras: 2 universal sockets attached to the bottom of the body, power cord length: 3 m
- packaging: 1 trolley



A101774

Trolley F3V

- 3 height adjustable, 8 mm tempered glass shelves (W 400 x L 270 mm)
- bearing strength: 10 kg each shelf
- lockable, low noise castors
- weight: 20 kg
- extras: 2 universal sockets attached to the bottom of the body, power cord length: 3 m
- packaging: 1 trolley

REF	material	colour	dimensions (W x D x H)
A101775	chrome plated steel	white	400 x 400 x 820 mm
REF	material		dimensions (W x D x H)
A101774	chrome plated steel/glass		400 x 400 x 820 mm

Additional haulage costs will incur. For further information please contact sales@wellsamed.de



A106174

Trolley Implanto Light L

- 3 tempered glass shelves (8 mm)
- 4 lockable, smooth-running und low-noise castors
- with air and water supply
- with power supply (2 universal sockets)
- shelf dimensions: W 600 x L 500 mm
- trolley dimensions: W 600 x D 540 x H 790 mm
- suitable for all (large) surgical units
- packaging: 1 trolley

REF	material	colour	dimensions (W x D x H)
A106174	stainless steel (powder-coated)/glass	white	600 x 540 x 790 mm

Trolley Implanto Max

- multifunctional cart for surgery/implantology
- 3 pull-out shelves
- practical cable routing
- top and front door: tempered glass and scratch resistant
- shock-resistant rear door with multiple socket outlet
- front door can be opened at a 270° angle
- sterilisable mounts
- antistatic, lockable castors
- dimensions: W 540 x D 510 x H 920 mm
- packaging: 1 cart

REF	material	dimensions (W x D x H)
A108933	powder-coated stainless steel/glass	540 x 510 x 920 mm

Additional haulage costs will incur. For further information please contact sales@wellsamed.de



A108933



A108934

Trolley Implanto Max S

- multifunctional cart for surgery/implantology
- 1 drawer and 2 pull-out shelves
- practical cable routing
- top and front door: tempered glass and scratch resistant
- shock-resistant rear door with multiple socket outlet
- front door can be opened at an 270° angle
- sterilisable mounts
- antistatic, lockable castors
- dimensions: W 540 x D 510 x H 920 mm
- packaging: 1 cart

REF	material	dimensions (W x D x H)
A108934	powder-coated stainless steel/glass	540 x 510 x 920 mm

Additional haulage costs will incur. For further information please contact sales@wellsamed.de



A104056

Attachment Tray Eco for Mounting Rod

- mountable on mounting rods with diameters of 35 to 55 mm
- without swivel joint
- packaging: 1 tray without decorative items



A104057

Attachment Tray Deluxe for Mounting Rod

- mountable on mounting rods with diameters of 35 to 55 mm
- with swivel joint
- packaging: 1 tray without decorative items

REF	material	load capacity	dimensions (W x D)
A104056	plastic	4-5 kg	32 x 23 cm
A104057	plastic	4-5 kg	35 x 25 cm

Dentist Chair Comfort

- circular shaped swivel seat
- especially upholstered seat and back for even more comfort
- gas spring to regulate height from 50 to 65 cm by sideward lever
- bearing strength: max. 150 kg
- 5 smooth-running castors
- colour of upholstery: light yellow
- packaging: 1 chair



A102723

REF	colour
A102723	light yellow

Dentist Chair Eco

- circular shaped swivel seat
- upholstered seat and back
- gas spring to regulate height from 50 to 65 cm by sideward lever
- bearing strength: max. 150 kg
- 5 smooth-running castors
- colour of upholstery: grey
- packaging: 1 chair



A102724

REF	colour
A102724	light grey

Dentist Chair Premium

- includes upholstery and chair frame
- ergonomic and solid design
- gummed castors
- height can be regulated from 48 to 60 cm



A107872

	A107873		A107869
	A107875		A107876
	A109326		A107877
	A109327		

Dentist Chair

- packaging: 1 chair including upholstery and chair frame



A109324

	A109325
	A109324

REF	material	colour
A107872	smooth synthetic leather	medium blue
A107873		light blue
A107875		yellow
A109326	rough synthetic leather	dark blue
A109327		turquoise
A107869		light brown
A107876		russet
A107877		grey

REF	material	colour
A109325	smooth synthetic leather	dark blue
A109324		grey



NEW COLOURS



146192



148191

Saddle Stool

- synthetic leather, soft castors
- lifting range: from 56 to 75 cm
- weight: 10 kg
- bearing strength: approx. 135 kg
- on demand available with foot pedal
- packaging: 1 stool



125439

125440

125441

125442

125443

140168

129555



REF	material	colour
125439	synthetic leather	mint green
125440		yellow
125441		light grey
125442		dark blue
125443		white

REF	material	colour
140168	synthetic leather	black
129555		poppy-red
148191		raspberry
146192		apple-green

INFORMATION ON HAULAGE COSTS

For deliveries within Germany (mainland) additional haulage costs of 50 € will accrue. For deliveries to all other areas (domestic and/or abroad) varying haulage costs will accrue. Please make a request to learn about your individual haulage costs. For further information please contact our service team at sales@wellsamed.de



LABORATORY

Medication Dispenser

- facilitates an accurate and hygienic dosage of medicaments
- removable dappen glasses
- closing lids can easily be opened by hand or by instrument
- silicone inserts make for secure closing of dappen glasses
- secure and non-slip purchase due to silicone feet
- packaging WMD 102: 1 set complete with 2 dappen glasses (clear)
- packaging WMD 103: 1 set complete with 3 dappen glasses (clear)



WMD 103



WMD 102

REF	material	description	dimensions (W x D x H)
WMD 102	stainless steel	two system	70 x 40 x 48 mm
WMD 103	casing	three system	100 x 40 x 48 mm

Microtorch Torch

- handy torch for laboratory and practice
- single-handed / one-handed use
- stand-alone torch
- automatic and reliable lighting /ignition
- flame size continuously adjustable and lockable
- max. burning period: approx. 35 min.
- packaging: 1 torch

REF	dimensions
WLB 020	H 15 cm



WLB 020

Stand Burner, refillable

- steady and non-slip
- with handle
- refillable
- max. temperature: approx. 1,300 °C
- max. burning period: approx. 60 min.
- packaging: 1 torch

REF	dimensions
WLB 030	11.5 x 8 x 8 cm



WLB 030



WLB 101

Universal Gas Lighter Refill

- butane gas
- for laboratory burners of every description (e.g. Blue-Flame, Kager, Kisag etc.)
- depending on burner, flame temperature reaches up to approx. 1,200 °C
- packaging: can of 300 ml

REF	packaging
WLB 101	300 ml

Antkoviak Wax Knife

- knife with wooden handle
- packaging: 1 knife



E 700

Antkoviak Spare Blade

- suitable for Antkoviak wax knife
- packaging: 1 pack of 2 blades



E 700.01

REF	dimensions	packaging
E 700	12.5 cm	1 piece
REF	packaging	
E 700.01	2 pcs.	



Work Pan with Plug-In Board

- made of plastic
- dimensions: 140 x 190 x 41 cm
- packaging: 1 work pan, 1 coloured plug-in board, 1 clip-spring

REF	colour
WMT 180	red
WMT 120	yellow
WMT 130	green
WMT 160	blue
WMT 110	white
WMT 100	black

Dental Pressure Lock Bags

- with insert pocket (kangaroo pouch)
- 3 label fields
- see-through
- with reinforced weld seams
- packaging: 1 pack of 100 pcs.



WVI 110

REF	material	dimensions (W x D x H)	packaging
WVI 110	polyethylene	180 x 250 x 0.05 mm	100 pcs.

PATIENT BROCHURES FOR YOUR DENTAL PRACTICE

- the right brochure for every relevant issue
- important information for your patients
- memo field on the back to note questions, references etc.
- 6-page leaflet

packaging:
1 pack of 25
brochures

All patient
brochures are
only available in
German!




**Angst –
Unterstützung
bei der Angst-
bewältigung**

REF: 121699




**Bleaching –
Professionelle
Zahnauf-
hellung**

REF: 121700




**Brücken –
Eine Brücke
beim Zahnarzt**

REF: 121701



**Gesunde
Ernährung
für gesunde
Zähne**

REF: 121703




**Funktions-
analysen –
Bei Verspan-
nungen in
Kopf, Kiefer &
Nacken**

REF: 121704




**Implantate –
Fast die
eigenen Zähne**

REF: 121705




**Karies –
Prophylaxe
gegen
Zahnverfall**

REF: 121706




**Kinderpro-
phylaxe –
Früherken-
nung und
Erziehung**

REF: 121707




**Die Krone –
Zum Schutz
des Zahns**

REF: 121708




**Narkose und
Sedierung –
In Ruhe zum
Zahnarzt**

REF: 121709




**Parodontitis –
Zahnfleisch-
erkrankung**

REF: 121710




**Schwanger-
schaft –
Hormon-
veränderun-
gen im Mund**

REF: 121712




**Veneers –
Ihr Hollywood
Lächeln**

REF: 121713




**Wurzel-
behandlung
(Endodontie)**

REF: 121702




**Ratschläge
zur Zahnt-
fernung**

REF: 121714




**Zahnerhaltung
durch Zahn-
füllungen**

REF: 121715



**Herausnehm-
barer Zahn-
ersatz – Die
Modellguss-
prothese**

REF: 121716



**Zahnreinigung,
-pflege und
professionelle
Zahnreinigung**

REF: 121711

ADD TO YOUR ORDER

single pocket
rack



REF: 100040

2-pocket
rack



REF: 100039

4-pocket
rack



REF: 126650



INSTRUMENTS

wellsamed Mouth Mirror

- made of plastic and real glass
- slightly bent mirror head
- ergonomic handle
- used for retraction of lips and cheeks
- packaging: 1 mirror

NEW

WDM 030




REF	material	colour	packaging
WDM 030	plastic/glass	assorted colours	1 piece
WDM 020	plastic/glass	assorted colours	1 piece

wellsamed Mouth Mirror

- made of plastic and real glass
- slightly bent mirror head
- individually wrapped in cellophane
- packaging: 1 mirror

NEW

WDM 020



NEW

Mirror Handle, hollow

- narrow form, angular
- polished/brushed stainless steel
- stainless, light
- packaging: 1 handle

E 400.02



NEW

Mirror Handle Finger Shape, hollow

- narrow form, angular
- polished/brushed stainless steel
- stainless, light
- packaging: 1 handle

E 400.07



NEW

Plane Surface Mouth Mirror Head

- with standardised, stainless thread
- suitable for all current disinfection and sterilisation procedures
- material: stainless steel
- packaging: 1 pack of 12 pcs.



E 1003

Mouth Mirror Head Hollow

- with standardised, stainless thread
- suitable for all current disinfection and sterilisation procedures
- material: stainless steel
- packaging: 1 pack of 12 pcs.



E 1013

Front Surface Mirror Rhodium-Coated

- with standardised, stainless thread
- suitable for all current disinfection and sterilisation procedures
- material: rhodium
- packaging: 1 pack of 12 pcs.



E 1023

REF	description		packaging
E 400.02	Mirror Handle, hollow		1 handle
E 400.07	Mirror Handle Finger Shape, hollow		1 handle
E 1003	Plane Surface Mouth Mirror Head	Nr. 3/Ø 20 mm	12 pcs.
E 1004		Nr. 4/Ø 22 mm	12 pcs.
E 1005		Nr. 5/Ø 24 mm	12 pcs.
E 1013	Mouth Mirror Head Hollow	Nr. 3/Ø 20 mm	12 pcs.
E 1014		Nr. 4/Ø 22 mm	12 pcs.
E 1015		Nr. 5/Ø 24 mm	12 pcs.
E 1023	Front Surface Mirror Rhodium-Coated	Nr. 3/Ø 20 mm	12 pcs.
E 1024		Nr. 4/Ø 22 mm	12 pcs.
E 1025		Nr. 5/Ø 24 mm	12 pcs.

Explorer, octagonal handle

- break-proof, thin tip
- firm and flexible
- stainless, sterilisable
- packaging: 1 piece



REF	material	colour	style
E 520.17	stainless steel	silver	Fig. 17
E 520.23		silver	Fig. 23
E 520.32		silver	Fig. 32
E 520.33		silver	Fig. 33

Explorer, round handle

- break-proof, thin tip
- firm and flexible
- stainless, sterilisable
- packaging: 1 piece



REF	material	colour	style
E 511.17	stainless steel	silver	Fig. 17
E 511.23		silver	Fig. 23
E 511.42		silver	Fig. 42
E 511.43		silver	Fig. 43


Pocket Probe CP-11

- break-proof, thin tip
- firm and flexible
- stainless, sterilisable
- packaging: 1 piece



WHO Pocket Probe

- break-proof, thin tip
- firm and flexible
- stainless, sterilisable
- packaging: 1 piece



REF	material	colour	packaging
E 510.09	stainless steel	silver	1 piece
E 510.07		silver	1 piece

Tweezer College

- grooved
- sterilisable
- packaging: 1 piece



REF	material	colour	dimensions
E 407.31	stainless steel	silver	15 cm


Surgical Aspirator Tip, Stainless Steel

- for precise and safe suction during all surgical procedures
- stainless
- slightly bent neck und extra fine tip
- autoclavable up to 134 °C
- inner diameter at the end: 1 mm
- packaging: 1 piece



Surgical Aspirator Tips, Plastic

- slightly bent neck and fine tip
- autoclavable up to 134 °C
- packaging: 1 pack of 20 pcs.



REF	material	colour	dimensions	packaging
WSA 101	plastic	white	L 153.5 mm/Ø 4 mm	20 pcs.
WSA 102		white	L 155 mm/Ø 5.4 mm	20 pcs.
WSA 103		white	L 157 mm/Ø 6.3 mm	20 pcs.
WSA 201	stainless steel	silver	L 155 mm/Ø 2.0 mm	1 piece
WSA 202		silver	L 155 mm/Ø 3.5 mm	1 piece




FIG. 10 FIG. 11 FIG. 12 FIG. 12D FIG. 15

Scalpel Blades, Stainless Steel

- rustproof
- sterile packaging
- packaging: 1 box of 100 pcs.

REF	material	blade size	packaging
E 330110	stainless steel	Fig. 10	100 pcs.
E 330111		Fig. 11	100 pcs.
E 330112		Fig. 12	100 pcs.
E 330113		Fig. 12D	100 pcs.
E 330115		Fig. 15	100 pcs.




FIG. 10 FIG. 11 FIG. 12 FIG. 12D FIG. 15 FIG. 15C


Scalpel Blades, Carbon Steel

- very sharp-edged
- uniform cutting edges, consistent sharpness
- sterile packaging
- packaging: 1 box of 100 pcs.

REF	material	blade size	packaging
E 330210	Carbon Steel	Fig. 10	100 pcs.
E 330211		Fig. 11	100 pcs.
E 330212		Fig. 12	100 pcs.
E 330213		Fig. 12D	100 pcs.
E 330215		Fig. 15	100 pcs.
E 330216		Fig. 15C	100 pcs.

Scalpel Handles, Stainless Steel

- sterilisable, solid handles to hold scalpel blades
- E330.01R with metric ruler
- packaging: 1 handle



REF	material	handle size	packaging
E 330.01	stainless steel	Fig. 03	1 piece
E 330.01R		Fig. 3G/3R	1 piece
E 330.02		Fig. 04	1 piece



Disposable Scalpels

- stainless
- consistent sharpness for clean and smooth incisions
- with blade protector
- plastic handle with metric ruler from 1 to 6 cm
- individual, sterile packaging
- packaging: 1 box of 10 pcs.

Safety Scalpels

- stainless
- consistent sharpness for clean and smooth incisions
- with plastic protection cap to avoid accidental cuts
- protection cap can be locked by pressing against the side
- ready to use by retracting protection cap
- ergonomic handle with scale from 1 to 5 cm
- individual, sterile packaging
- packaging: 1 box of 10 pcs.



REF	material	colour	description	packaging
E 330410	stainless steel	brown	Form 10	10 pcs.
E 330411		red	Form 11	10 pcs.
E 330412		yellow	Form 12	10 pcs.
E 330415		violet	Form 15	10 pcs.

REF	material	colour	description	packaging
E 330510	stainless steel	brown	Form 10	10 pcs.
E 330511		red	Form 11	10 pcs.
E 330515		violet	Form 15	10 pcs.

wellsaview Cheek Retractors, long handle

- autoclavable up to 121 °C
- long handle to ensure optimum visibility and stability
- comfortable for the patient
- excellent for intraoral photography
- for children (WLW 001) or adults (WLW 002)
- material: plastic
- packaging: sold in pairs

REF	material	colour	dimensions	packaging
WLW 001	plastic	transparent	Small	2 pcs.
WLW 002		transparent	Standard	2 pcs.



wellsaview Cheek Retractors, double-ended

- autoclavable up to 121 °C
- for children (WLW 011) and adults (WLW 012)
- material: plastic
- packaging: sold in pairs

REF	material	colour	dimensions	packaging
WLW 011	plastic	transparent	small	2 pcs.
WLW 012		transparent	medium	2 pcs.



wellsaview Cheek Retractors

- autoclavable up to 121 °C
- self-span cheek retractor
- slips into place easily and comfortably
- for children (WLR 001) and adults (WLR 002)
- material: plastic
- packaging: sold in pairs



REF	material	colour	dimensions	packaging
WLR 001	plastic	transparent	small	2 pcs.
WLR 002		transparent	standard	2 pcs.

wellsaview Lip and Cheek Retractor Set

- anatomic design
- maximum visibility and patient comfort
- material: polycarbonate plastic
- packaging: 1 set with 1 x small and 1 x standard retractor



REF	material	colour	dimensions	packaging
WLR 010	plastic	transparent	small + standard	1 set

wellsaview Cheek Retractor ECO

- autoclavable up to 121 °C
- slips into place easily and comfortably
- for children (WLR 011) and adults (WLR 012)
- packaging: 1 retractor



REF	material	colour	dimensions	packaging
WLR 011	plastic	transparent	small	1 piece
WLR 012		transparent	standard	1 piece

wellsaview Cheek Retractors ECO, long handle

- autoclavable up to 121 °C
- packaging: sold in pairs



REF	material	colour	dimensions (W x L)	packaging
WLW 015	plastic	transparent	70 x 175 mm	2 pcs.

Burs Holder

- compatible with EMS
- suitable specifically for the anterior region of the mouth
- compatible with various bur-attachments
- packaging: 1 bur holder



AE8


REF	material	size	packaging
AE8	stainless steel	AE8	1 piece

Bur Cleaning Brush

- designed for rough cleaning of instruments
- packaging: 1 brush



REF	material	packaging
WDC 012	brass	1 piece
WDC 011	nylon	1 piece




WIC 101

Cassette 01

- perforated stainless steel cassette for 3 to 5 instruments
- sterilisable, autoclavable and thermoisinfectable
- with ergonomic twist lock
- detachable lid
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

REF	material	colour	dimensions (W x D x H)
WIC 101	stainless steel	silver	205 x 84 x 32 mm



WIC 102

Cassette 02

- perforated stainless steel cassette for 5 instruments
- sterilisable, autoclavable and thermoisinfectable
- with ergonomic twist lock and separate latch
- detachable lid
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

REF	material	colour	dimensions (W x D x H)
WIC 102	stainless steel	silver	205 x 84 x 38 mm



WIC 103

Cassette 03

- perforated stainless steel cassette for 6 instruments
- sterilisable, autoclavable and thermoisinfectable
- with ergonomic twist lock and separate chamber for large instruments
- detachable lid
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

REF	material	colour	dimensions (W x D x H)
WIC 103	stainless steel	silver	205 x 140 x 33 mm



WIC 104

Cassette 04

- perforated stainless steel cassette for 20 instruments
- sterilisable, autoclavable and thermodisinfectable
- with ergonomic twist lock
- detachable lid
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

REF	material	colour	dimensions (W x D x H)
WIC 104	stainless steel	silver	280 x 203 x 32 mm



WIC 107

Cassette 07

- perforated stainless steel cassette for 6 instruments
- sterilisable, autoclavable and thermodisinfectable
- detachable fixture
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

REF	material	colour	dimensions (W x D x H)
WIC 107	stainless steel	silver	185 x 110 x 20 mm



WIC 108

Cassette 08

- perforated stainless steel cassette for 6 instruments
- sterilisable, autoclavable and thermodisinfectable
- detachable fixture
- bottom part with enlarged gaps
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

REF	material	colour	dimensions (W x D x H)
WIC 108	stainless steel	silver	185 x 110 x 20 mm



WIC 109

Cassette 09

- perforated stainless steel cassette for 6 instruments
- sterilisable, autoclavable and thermodisinfectable
- detachable lid
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

REF	material	colour	dimensions (W x D x H)
WIC 109	stainless steel	silver	185 x 110 x 20 mm



WIC 110


Cassette 10

- perforated stainless steel cassette for 6 instruments
- sterilisable, autoclavable and thermoisinfectable
- detachable lid
- bottom part with enlarged gaps
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

REF	material	colour	dimensions (W x D x H)
WIC 110	stainless steel	silver	185 x 110 x 20 mm

Cassette 11


- perforated stainless steel cassette for 3 instruments
- sterilisable, autoclavable, and thermoisinfectable
- detachable fixture
- bottom part with enlarged gaps
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette



WIC 111


Cassette 12

- perforated stainless steel cassette for 3 instruments
- sterilisable, autoclavable, and thermoisinfectable
- detachable fixture
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette




WIC 112

REF	material	colour	dimensions (W x D x H)
WIC 111	stainless steel	silver	185 x 55 x 25 mm
WIC 112		silver	185 x 55 x 25 mm



WIC 113



WIC 114

Cassette 13

- perforated stainless steel cassette for 3 instruments
- sterilisable, autoclavable, and thermoisinfectable
- detachable lid
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

Cassette 14

- perforated stainless steel cassette for 3 instruments
- sterilisable, autoclavable, and thermoisinfectable
- detachable lid
- bottom part with enlarged gaps
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette

REF	material	colour	dimensions (W x D x H)
WIC 113	stainless steel	silver	185 x 55 x 25 mm
WIC 114		silver	185 x 55 x 25 mm

Cassette 17

- perforated, stainless steel cassette for 10 instruments
- sterilisable, autoclavable and thermoisinfectable
- with ergonomic twist lock
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette




WIC 117

REF	material	colour	dimensions (W x D x H)
WIC 117	stainless steel	silver	205 x 140 x 32 mm


Cassette 18

- perforated, stainless steel cassette for 20 instruments
- sterilisable, autoclavable and thermoisinfectable
- with ergonomic twist lock
- detachable lid
- individually insertable silicone rails
- non-slip
- packaging: 1 cassette



WIC 118

REF	material	colour	dimensions (W x D x H)
WIC 118	stainless steel	silver	280 x 202 x 33 mm




WAT 109

Standard Perforated Aluminium Tray, Lid

- not suitable for thermal disinfection
- suitable for sterilisation reel WSR 305
- instrument inserts sold separately
- packaging: 1 tray

REF	material	colour	dimensions (W x D x H)
WAT 109	aluminium	silver	288 x 187 x 29 mm



WAT 160

Standard Perforated Aluminium Tray, Bottom Part

- not suitable for thermal disinfection
- suitable for sterilisation reel WSR 305
- instrument inserts sold separately
- packaging: 1 tray

REF	material	colour	dimensions (W x D x H)
WAT 110	aluminium	silver	284 x 183 x 17 mm
WAT 160		blue	284 x 183 x 17 mm
WAT 130		green	284 x 183 x 17 mm
WAT 120		gold	284 x 183 x 17 mm
WAT 180		red	284 x 183 x 17 mm



WAT 209

Standard Unperforated Aluminium Tray, Lid

- not suitable for thermal disinfection
- suitable for sterilisation reel WSR 305
- instrument inserts sold separately
- packaging: 1 tray

REF	material	colour	dimensions (W x D x H)
WAT 209	aluminium	silver	288 x 187 x 29 mm



WAT 280

Standard Unperforated Aluminium Tray, Bottom Part


- not suitable for thermal disinfection
- suitable for sterilisation reel WSR 305
- instrument inserts sold separately
- packaging: 1 tray

REF	material	colour	dimensions (W x D x H)
WAT 210	aluminium	silver	284 x 183 x 17 mm
WAT 260		blue	284 x 183 x 17 mm
WAT 230		green	284 x 183 x 17 mm
WAT 220		gold	284 x 183 x 17 mm
WAT 280		red	284 x 183 x 17 mm



Standard Perforated Stainless Steel Tray, Bottom Part

- suitable for sterilisation reel WSR 305
- instrument inserts sold separately
- packaging: 1 tray




WST 100

REF	material	colour	dimensions (W x D x H)
WST 100	stainless steel	silver	284 x 183 x 17 mm

Standard Perforated Stainless Steel Tray, Lid

- suitable for sterilisation reel WSR 305
- instrument inserts sold separately
- packaging: 1 tray



WST 109

REF	material	colour	dimensions (W x D x H)
WST 109	stainless steel	silver	288 x 187 x 29 mm



WST 209

Standard Unperforated Stainless Steel Tray, Lid

- not suitable for thermal disinfection
- suitable for sterilisation reel WSR 305
- instrument inserts sold separately
- packaging: 1 tray

REF	material	colour	dimensions (W x D x H)
WST 209	stainless steel	silver	288 x 187 x 29 mm




WST 200

Standard Unperforated Stainless Steel Tray, Bottom Part

- instrument inserts sold separately
- packaging: 1 tray

REF	material	colour	dimensions (W x D x H)
WST 200	stainless steel	silver	284 x 183 x 17 mm



Instrument Tray Insert – Aluminium

- not suitable for thermal disinfection
- suitable for standard aluminium trays
- for 11 or 17 instruments
- packaging: 1 insert

REF	material	colour	dimensions
WAT 011	aluminium	silver	for 11 Instruments
WAT 017		silver	for 17 Instruments

Instrument Tray Insert – Stainless Steel

- suitable for standard stainless steel trays
- for 16 or 23 instruments
- packaging: 1 insert

REF	material	colour	dimensions
WST 016	stainless steel	silver	for 16 Instruments
WST 023		silver	for 23 Instruments



WCB 107



WCB 101

Sharps Container

- convenient and safe way to dispose needles
- safe transportation of used medical needles and other sharp and contaminated medical instruments
- with extractor tool and lid
- packaging: 1 piece

REF	material	colour	dimensions
WCB 101	plastic	yellow/red	1.5 litres
WCB 107		yellow/red	7 litres

Automatic Needle Remover

- easy and safe disposal of medical needles, burs or files
- collecting tank can be emptied by pulling it out sideways
- battery-powered
- no risk of injury or infection
- packaging: 1 device

clean + safe



A103753



A103749

REF	material	colour	dimensions (W x D x H)
A103749	plastic	dark blue	188 x 108 x 160 mm
A103753		beige	188 x 108 x 160 mm

HANDLE SELECTION

HANDLE SELECTION

You can adapt our Teca Dent instruments to your personal needs. The handles are available in five different designs:

The premium handle **TECA Light™** is super light and still lies firmly in the hand. The instrument is made of high-alloy super steel (440er) and thereby remains consistently sharp.

The slimmer **Artec™** handle is manufactured from smooth, hard plastic and intended for simple clinical procedures. The instrument consists of high-strength stainless steel AISI 419/420 and thermo-resistant plastic.

The slightly broader **Dentec™** handle is made of soft, elastic silicone and thus a bit heavier. To ensure our instruments' long-lasting sharpness, we also use high-strength stainless steel AISI 410/420 for Dentec handles.

Additionally we offer the **ErgoTec™** handle for our single-ended instruments. It consists of silicone and offers a different grip due to its slim and even design.

The wavelike grooves and cross section texture of our new **Teca Colour™** handle provides for better instrument handling and an easier, safer grip.

Every instrument is characterised by perfect form and fine design. When selecting your handle, the most important aspects to consider are its grip and handiness.

We recommend to change between different handles to relieve your hands.

THE CORRECT MAINTENANCE

The correct maintenance of Teca Dent instruments is important to maintain

long durability and optimum efficiency. To maintain our instruments' high quality for a long time, we recommend steam sterilisation which should never exceed a temperature of 134°C. For more detailed information about the correct cleaning and sterilisation procedures, please contact our website www.tecadent.de. For further questions please contact sales@wellsamed.de

*WARRANTY

We are committed to the quality of our products. Therefore we include a 5-year warranty on all defects in material or workmanship. To raise a warranty claim you are required to present correct maintenance of your instruments according to our maintenance instructions which can be found on our website (www.tecadent.de). wellsamed denies liability for any deviations.

Teca Colour™

C for Teca Colour Handles



TECA Light™

L for Teca Light Handles



Artec™

A for Artec Handles



Dentec™

D for Dentec Handles



ErgoTec™

E for Ergo Tec Handles



Mirror Handle ARTEC

- mirror handle thread in accordance with international standard M 2.5 - the mirror shank
- material: hard plastic
- packaging: 1 handle



REF	description	colour	packaging
A 7010 R	Mirror Handle, Artec	red	1 piece
A 7010 Y		yellow	1 piece
A 7010 O		orange	1 piece
A 7010 V		violet	1 piece
A 7010 B		blue	1 piece
A 7010 G		green	1 piece
A 7010 P		pink	1 piece

Mirror Handle Ergo Tec

- mirror handle thread in accordance with international standard M 2.5 - the mirror shank
- material: soft plastic
- packaging: 1 handle



REF	description	colour	packaging
E 7011 Y	Mirror Handle, Ergo Tec	yellow	1 piece
E 7011 O		orange	1 piece
E 7011 P		pink	1 piece
E 7011 B		blue	1 piece
E 7011 V		violet	1 piece
E 7011 R		red	1 piece
E 7011 LB		green	1 piece
E 7011 GR		grey	1 piece

Mirror Handle

- well balanced, better accessibility
- packaging: 1 handle



REF	description	colour	packaging
A 1025 R	Mirror Handle, Artec	red	1 piece
A 1025 B		blue	1 piece
A 1025 V		violet	1 piece
A 1025 G		green	1 piece
D 1025 R	Mirror Handle, Dentec	red	1 piece
D 1025 B		blue	1 piece
D 1025 V		violet	1 piece
D 1025 G		green	1 piece
L 1025	Mirror Handle, Teca Light	silver	1 piece

Mirror Handle Ortho

- millimetre scale/total length: 25 mm
- to measure sagittal and vertical overbite
- packaging: 1 handle



REF	description	colour	packaging
A 1028 R	Mirror Handle ortho, Artec	red	1 piece
A 1028 B		blue	1 piece
A 1028 V		violet	1 piece
A 1028 G		green	1 piece
D 1028 R	Mirror Handle ortho, Dentec	red	1 piece
D 1028 B		blue	1 piece
D 1028 V		violet	1 piece
D 1028 G		green	1 piece
L 1028	Mirror Handle ortho, Teca Light	silver	1 piece

Curettes

- to remove supragingival calculus
- design: triangular cross section, face perpendicular to the lower shank, two cutting edges
- Langer Curette combines the Universal Curette blade and the Gracey Curette shank

Goldman-Fox 3 Curette

A 2031	Artec	1 piece
D 2031	Dentec	1 piece
L 2031	Teca Light	1 piece

**Goldman-Fox 4 Curette**

A 2032	Artec	1 piece
D 2032	Dentec	1 piece
L 2032	Teca Light	1 piece

**Langer 1/2 Curette**

A 2025	Artec	1 piece
D 2025	Dentec	1 piece
L 2025	Teca Light	1 piece

**Langer 3/4 Curette**

A 2026	Artec	1 piece
D 2026	Dentec	1 piece
L 2026	Teca Light	1 piece

**Langer 5/6 Curette**

A 2027	Artec	1 piece
D 2027	Dentec	1 piece
L 2027	Teca Light	1 piece

**Mini Langer 1/2 Curette**

A 2028	Artec	1 piece
D 2028	Dentec	1 piece
L 2028	Teca Light	1 piece

**Mini Langer 3/4 Curette**

A 2029	Artec	1 piece
D 2029	Dentec	1 piece
L 2029	Teca Light	1 piece

**Mini Langer 5/6 Curette**

A 2030	Artec	1 piece
D 2030	Dentec	1 piece
L 2030	Teca Light	1 piece

**Universal Curette McCall 13S-14S**

A 2020	Artec	1 piece
D 2020	Dentec	1 piece
L 2020	Teca Light	1 piece

**Mini Universal Curette McCall 13S-14S**

A 2021	Artec	1 piece
D 2021	Dentec	1 piece
L 2021	Teca Light	1 piece

**Universal Curette McCall 15-16**

A 2022	Artec	1 piece
D 2022	Dentec	1 piece
L 2022	Teca Light	1 piece
C 2022	Teca Colour	1 piece

**Universal Curette McCall 17-18**

A 2023	Artec	1 piece
D 2023	Dentec	1 piece
L 2023	Teca Light	1 piece
C 2023	Teca Colour	1 piece

**Universal Curette McCall 17S-18S**

A 2024	Artec	1 piece
D 2024	Dentec	1 piece
L 2024	Teca Light	1 piece

**Columbia 2L-2R Curette**

A 2033	Artec	1 piece
L 2033	Teca Light	1 piece

**Columbia 4L-4R Curette**

A 2034	Artec	1 piece
L 2034	Teca Light	1 piece

**Columbia 13-14 Curette**

A 2035	Artec	1 piece
L 2035	Teca Light	1 piece

**Interproximal Curette SV 1-3**

A 2039	Artec	1 piece
D 2039	Dentec	1 piece
L 2039	Teca Light	1 piece

**Interproximal Curette SV 2-4**

A 2040	Artec	1 piece
D 2040	Dentec	1 piece
L 2040	Teca Light	1 piece

**Gracey 1/2 Curette**

A 2041	Artec	1 piece
L 2041	Teca Light	1 piece
C 2041	Teca Colour	1 piece

**Gracey 3/4 Curette**

A 2042	Artec	1 piece
L 2042	Teca Light	1 piece



Gracey 5/6 Curette

A 2043	Artec	1 piece
L 2043	Teca Light	1 piece
C 2043	Teca Colour	1 piece

**Gracey 7/8 Curette**

A 2044	Artec	1 piece
L 2044	Teca Light	1 piece

**Gracey 9/10 Curette**

A 2045	Artec	1 piece
L 2045	Teca Light	1 piece
C 2045	Teca Colour	1 piece

**Gracey 11/12 Curette**

A 2046	Artec	1 piece
L 2046	Teca Light	1 piece
C 2046	Teca Colour	1 piece

**Gracey 13/14 Curette**

A 2047	Artec	1 piece
L 2047	Teca Light	1 piece
C 2047	Teca Colour	1 piece

**Gracey 15/16 Curette**

A 2048	Artec	1 piece
L 2048	Teca Light	1 piece

**Gracey 17/18 Curette**

A 2049	Artec	1 piece
L 2049	Teca Light	1 piece

**Macro Gracey 5/6 Curette**

A 2050	Artec	1 piece
D 2050	Dentec	1 piece
L 2050	Teca Light	1 piece

**Macro Gracey 7/8 Curette**

A 2051	Artec	1 piece
D 2051	Dentec	1 piece
L 2051	Teca Light	1 piece

**Macro Gracey 11/12 Curette**

A 2052	Artec	1 piece
D 2052	Dentec	1 piece
L 2052	Teca Light	1 piece

**Macro Gracey 13/14 Curette**

A 2053	Artec	1 piece
D 2053	Dentec	1 piece
L 2053	Teca Light	1 piece

**Macro Gracey 11/14 Curette**

A 2054	Artec	1 piece
D 2054	Dentec	1 piece
L 2054	Teca Light	1 piece

**Macro Gracey 12/13 Curette**

A 2055	Artec	1 piece
D 2055	Dentec	1 piece
L 2055	Teca Light	1 piece

**Mini Gracey 1/2 Curette**

A 2056	Artec	1 piece
D 2056	Dentec	1 piece
L 2056	Teca Light	1 piece
C 2056	Teca Colour	1 piece

**Mini Gracey 7/8 Curette**

A 2057	Artec	1 piece
D 2057	Dentec	1 piece
L 2057	Teca Light	1 piece
C 2057	Teca Colour	1 piece

**Mini Gracey 11/12 Curette**

A 2058	Artec	1 piece
D 2058	Dentec	1 piece
L 2058	Teca Light	1 piece
C 2058	Teca Colour	1 piece

**Mini Gracey 13/14 Curette**

A 2059	Artec	1 piece
D 2059	Dentec	1 piece
L 2059	Teca Light	1 piece
C 2059	Teca Colour	1 piece

**Mini Gracey 15/16 Curette**

A 2060	Artec	1 piece
D 2060	Dentec	1 piece
L 2060	Teca Light	1 piece
C 2060	Teca Colour	1 piece

**Mini Gracey 17/18 Curette**

A 2061	Artec	1 piece
D 2061	Dentec	1 piece
L 2061	Teca Light	1 piece
C 2061	Teca Colour	1 piece



NEW

NEW

Sinus Lift Curette #11

A 9060	Artec	1 piece
D 9060	Dentec	1 piece
L 9060	Teca Light	1 piece



Sinus Lift Curette #12

A 9061	Artec	1 piece
D 9061	Dentec	1 piece
L 9061	Teca Light	1 piece



Sinus Lift Curette #13

A 9062	Artec	1 piece
D 9062	Dentec	1 piece
L 9062	Teca Light	1 piece



Sinus Lift Curette #14

A 9063	Artec	1 piece
D 9063	Dentec	1 piece
L 9063	Teca Light	1 piece



Sinus Lift Curette #15

A 9064	Artec	1 piece
D 9064	Dentec	1 piece
L 9064	Teca Light	1 piece



Sinus Lift Curette #16

A 9065	Artec	1 piece
D 9065	Dentec	1 piece
L 9065	Teca Light	1 piece



Lucas Curette (Spoon 2.4 mm)

A 9070	Artec	1 piece
D 9070	Dentec	1 piece
L 9070	Teca Light	1 piece



Lucas Curette (Spoon 2.7 mm)

A 9071	Artec	1 piece
D 9071	Dentec	1 piece
L 9071	Teca Light	1 piece
C 9071	Teca Colour	1 piece



Lucas Curette (Spoon 3.0 mm)

A 9072	Artec	1 piece
D 9072	Dentec	1 piece
L 9072	Teca Light	1 piece
C 9072	Teca Colour	1 piece



Lucas Curette (Spoon 3.8 mm)

A 9073	Artec	1 piece
D 9073	Dentec	1 piece
L 9073	Teca Light	1 piece



NEW

Miller Curette (Spoon 2.8 mm)

A 9074	Artec	1 piece
D 9074	Dentec	1 piece
L 9074	Teca Light	1 piece



Miller Curette (Spoon 3.5 mm)

A 9075	Artec	1 piece
D 9075	Dentec	1 piece
L 9075	Teca Light	1 piece



Miller Curette (Spoon 4.5 mm)

A 9076	Artec	1 piece
D 9076	Dentec	1 piece
L 9076	Teca Light	1 piece



Hirschfeld #3-7 (Blade 1.4 mm)

A 9077	Artec	1 piece
D 9077	Dentec	1 piece
L 9077	Teca Light	1 piece



Hirschfeld #5-11 (Blade 1.4 mm)

A 9078	Artec	1 piece
D 9078	Dentec	1 piece
L 9078	Teca Light	1 piece



Hirschfeld #9-10 (Blade 1.4 mm)

A 9079	Artec	1 piece
D 9079	Dentec	1 piece
L 9079	Teca Light	1 piece



Orban #10-11 (Blade 2.0 mm)

A 9080	Artec	1 piece
D 9080	Dentec	1 piece
L 9080	Teca Light	1 piece



Orban #12-13 (Blade 2.0 mm)

A 9081	Artec	1 piece
D 9081	Dentec	1 piece
L 9081	Teca Light	1 piece



Sugarman #1-2

A 9082	Artec	1 piece
D 9082	Dentec	1 piece
L 9082	Teca Light	1 piece



Sugarman #3-4

A 9083	Artec	1 piece
D 9083	Dentec	1 piece
L 9083	Teca Light	1 piece



Richard #1-2 (Blade 1.7 mm)

A 9084	Artec	1 piece
D 9084	Dentec	1 piece
L 9084	Teca Light	1 piece



Molt #2-4

A 9085	Artec	1 piece
D 9085	Dentec	1 piece
L 9085	Teca Light	1 piece



Heidbrink #2-3

A 9086	Artec	1 piece
D 9086	Dentec	1 piece
L 9086	Teca Light	1 piece



Davis #11R

A 9087	Artec	1 piece
D 9087	Dentec	1 piece
L 9087	Teca Light	1 piece



Davis #11P

A 9088	Artec	1 piece
D 9088	Dentec	1 piece
L 9088	Teca Light	1 piece



Ochsenbein #4 (Blade 3.5 mm)

A 9089	Artec	1 piece
D 9089	Dentec	1 piece
L 9089	Teca Light	1 piece



Rhodes #36-37

A 9090	Artec	1 piece
D 9090	Dentec	1 piece
L 9090	Teca Light	1 piece



Perio Set

- Gracey 13/14
- Gracey 11/12
- Gracey 1/2
- Mini Sickle
- Curette McCall 17-18
- cassette
- Dentec handle



Prophylaxis Basic-Set

- Explorer 23-52 B
- Micro Sickle Scaler
- Mini Gracy 1/2
- Mini Gracy 11/12
- Mini Gracy 13/14
- cassette
- Dentec handle



REF	description	packaging
D 6010	Perio Set	1 set
D 6070	Prophylaxis Basic-Set	1 set

PERFECT COMPATIBILITY


EMS®

The operating frequency range of the EMS®-scaler is <32.000 Hz. ARTEC tips perfectly adapt to these operating conditions.

SATELEC®

The operating frequency range of the SATELEC®-scaler is 28.000-32.000 Hz. Perfect compatibility of ARTEC tips is ensured.

HIGHEST QUALITY




ARTEC TIPS
The magnifications clearly display the quality and show the difference to competitive products.

TIPS OF OTHER MANUFACTURERS

The illustrations display tip surface 160 times magnified.



Tip of well-known Chinese brand A




Tip of well-known European brand A




Tip of well-known European brand B

MAXIMUM DURABILITY




- › use of selected biomaterial
- › special polishing techniques on surface
- › strictest quality controls
- › finest and rustproof tip surfaces
- › no risk of injury for teeth



- › increased durability due to extra hard titanium nitride coating
- › better check for wear* and corrosion protection
- › low friction and enhanced tissue compatibility

* If the golden coating is no longer visible, changing the tip is recommended.

ORDER FREE OF CHARGE



Even more tips can be found in our Tip Book.
Order a free copy now!
Only available in German!



wellsamed GmbH
Tel.: +49 341 27139-120
E-Mail: sales@wellsamed.de

AG1

Removal of supragingival deposits in all quadrants.



	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AG1	ATG1	EMS® MECTRON®
REF	AGD1	ATGD1	SATELEC® NSK® G8
REF	AGS1		SIRONA®



Mode
G

Power



Irrigation





AG2



Removal of heavy supragingival deposits. Apply flat end to surface of teeth.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AG2	ATG2	EMS® P MECTRON®
REF	AGD2	ATGD2	SATELEC® NSK® G2
REF	AGS2		SIRONA®



Mode
G

Power



Irrigation





AG3



Removal of supragingival deposits in all quadrants, including the interproximal and sulcus areas.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AG3	ATG3	EMS® P MECTRON®
REF	AGD3	ATGD3	SATELEC® NSK®
REF	AGS3		SIRONA®



Mode
G

Power



Irrigation




AG4



Recommended for treatment of interproximal spaces and for supragingival scaling.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AG4	ATG4	EMS® A MECTRON®
REF	AGD4	ATGD4	SATELEC® NSK®
REF	AGS4		SIRONA®



Mode
G

Power



Irrigation





AG5



Recommended for treating simple cases and gross supragingival scaling.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AG5	ATG5	EMS [®] MECTRON [®] S1
REF	AGD5	ATGD5	SATELEC [®] 1 NSK [®]
REF	AGS5		SIRONA [®]



Mode
G



Power	(+) ☒ ☒ ☒ ((•))	Irrigation	☒ ☒ ☒ ☒
-------	---	------------	---

AG6



Recommended for removing voluminous supragingival deposits. Apply flat end to surface of teeth.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AG6	ATG6	EMS [®] MECTRON [®] S3
REF	AGD6	ATGD6	SATELEC [®] 2 NSK [®]
REF	AGS6		SIRONA [®]



Mode
G



Power	(+) ☒ ☒ ☒ ((•))	Irrigation	☒ ☒ ☒ ☒
-------	---	------------	---

AG7



Removal of crowns. Apply the instrument onto surface and activate. Increase pressure until vibrations can no longer be heard and maintain for few seconds.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AG7	ATG7	EMS [®] MECTRON [®]
REF	AGD7	ATGD7	SATELEC [®] NSK [®] G25
REF	AGS7		SIRONA [®]



Mode
G



Power	(+) ☒ ☒ ☒ ((•))	Irrigation	☒ ☒ ☒ ☒
-------	---	------------	---

AGK1



Removal of supragingival deposits in all quadrants.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AGK1	ATGK1	KAVO® 5



Mode
G



Power	(*) ☒ ☒ ☒ ((*))	Irrigation	☒ ☒ ☒ ☒
-------	---	------------	---

AGK2



Removal of severe supragingival calculus and dental plaque.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AGK2	ATGK2	KAVO® 6



Mode
P



Power	(*) ☒ ☒ ☒ ((*))	Irrigation	☒ ☒ ☒ ☒
-------	---	------------	---

AGK3



Removal of supragingival calculus, interdental calculus, and calculus at the neck of the tooth.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AGK3	ATGK3	KAVO® 7



Mode
P

Power	(*) ☒ ☒ ☒ ((*))	Irrigation	☒ ☒ ☒ ☒
-------	---	------------	---

AGK4



Removal of supragingival deposits, including the interproximal and sulcus areas.

STAINLESS STEEL		TITANIUM	COMPATIBLE WITH	
REF	AGK4	ATGK4	KAVO®	8
Power				Irrigation

Mode
G

AGK5



Removal of subgingival calculus.

STAINLESS STEEL		TITANIUM	COMPATIBLE WITH	
REF	AGK5	ATGK5	KAVO®	60
Power				Irrigation

Mode
G

AGK6



Left-angled instrument to remove subgingival calculus.

STAINLESS STEEL		TITANIUM	COMPATIBLE WITH	
REF	AGK6	ATGK6	KAVO®	61
Power				Irrigation



Mode
G

AGK7



Right-angled instrument to remove subgingival calculus.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AGK7	ATGK7	KAVO® 62



Mode
G

Power



Irrigation





AP1



Slim and sharp instrument recommended for cleaning and irrigation of subgingival deposits.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AP1	ATP1	EMS® MECTRON® S2
REF	APD1	ATPD1	SATELEC® NSK®
REF	APS1		SIRONA®



Mode
P

Power



Irrigation




AP2L



Left-angled instrument used to remove calculus from very narrow inter-root spaces and furcation.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AP2L	ATP2L	EMS® MECTRON®
REF	APD2L	ATPD2L	SATELEC® P2R NSK®



Mode
P

Power



Irrigation



AP2LD



Left-angled, diamond-coated (40 µm) instrument for root planing.



STAINLESS STEEL

TITANIUM


COMPATIBLE WITH

REF AP2LD

EMS[®] MECTRON[®]

REF APD2LD

SATELEC[®] H2R NSK[®]



Mode
P

Power



Irrigation



AP2R

Right-angled instrument used to remove calculus from very narrow inter-root spaces and furcation.



STAINLESS STEEL

TITANIUM

COMPATIBLE WITH

REF AP2R



ATP2R

EMS[®] MECTRON[®]

REF APD2R

ATPD2R

SATELEC[®] P2L NSK[®]



Mode
P

Power



Irrigation



AP2RD



Right-angled, diamond-coated (40 µm) instrument for root planing.



STAINLESS STEEL

TITANIUM



COMPATIBLE WITH

REF AP2RD

EMS[®] MECTRON[®]

REF APD2RD

SATELEC[®] H2L NSK[®]



Mode
P

Power



Irrigation



AP3



Recommended for cleaning and irrigating periodontal deep pockets.

STAINLESS STEEL		TITANIUM	COMPATIBLE WITH	
REF	AP3	ATP3	EMS [®]	PL3
			MECTRON [®]	PP10
REF	APD3	ATPD3	SATELEC [®]	NSK [®]
REF	APS3		SIRONA [®]	PE1

Power	(*) ☒ ☒ ☒ ☒ (•••)	Irrigation	☒ ☒ ☒ ☒

Mode
P

AP3D



Diamond-coated (40 µm) instrument for root planing after the periodontal flap surgery. Also suitable for smoothing restoration margins and for expanding root furcations.

STAINLESS STEEL		TITANIUM	COMPATIBLE WITH	
REF	AP3D		EMS [®]	DPL3
			MECTRON [®]	
REF	APD3D		SATELEC [®]	NSK [®]
REF	APS3D		SIRONA [®]	

Power	(*) ☒ ☒ ☒ ☒ (•••)	Irrigation	☒ ☒ ☒ ☒

Mode
P

AP4



Recommended for cleaning and irrigation of periodontal shallow pockets.

STAINLESS STEEL		TITANIUM	COMPATIBLE WITH	
REF	AP4	ATP4	EMS [®]	MECTRON [®]
REF	APD4	ATPD4	SATELEC [®]	10P
			NSK [®]	G4
REF	APS4		SIRONA [®]	

Power	(*) ☒ ☒ ☒ ☒ (•••)	Irrigation	☒ ☒ ☒ ☒

Mode
P

AP4D



Diamond-coated (40µm) instrument for root canal orientation and removal of calcifications in the coronal third of the root canal system. Do not exert too much pressure to avoid hurting the root canal.

STAINLESS STEEL

TITANIUM

COMPATIBLE WITH

REF AP4D


EMS® RT1 MECTRON®

REF APD4D

SATELEC® PM4 NSK®

REF APS4D

SIRONA®



Mode
E

Power



Irrigation



AUP1




Recommended for periodontal gentle curetting scaling.

TITANIUM

COMPATIBLE WITH

REF AUP1

MECTRON® PS1



Mode
ROOT

Power

PERIO

Irrigation



AUP2




Inflammatory tissue removal and fractured root apex extraction.

TITANIUM

COMPATIBLE WITH

REF AUP2

MECTRON® PS2



Mode
ROOT

Power

PERIO

Irrigation



AUP3





Recommended for angled periodontal gentle curetting scaling.

TITANIUM

COMPATIBLE WITH

REF AUP3

MECTRON® PS6



Mode
ROOT

Power

PERIO

Irrigation



AUP4





Diamond-coated (40µm) instrument for root debridement and root planing.

TITANIUM

COMPATIBLE WITH

REF AUP4

MECTRON® OP5



Mode
ROOT

Power

PERIO

Irrigation



AUP5




Recommended for root surface micro-smoothing.

TITANIUM

COMPATIBLE WITH

REF AUP5

MECTRON® PP1



Mode
ROOT

Power

PERIO

Irrigation



AUP6



Mode
ROOT





Diamond-coated (40µm) instrument for micro-osteoplasty: interproximal osteoplasty and root planing.

TITANIUM

COMPATIBLE WITH

REF AUP6

MECTRON® OP6



Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AE1



Mode
E

120° holder used to hold the file (Ø 0,7 or Ø 0,8 mm). Recommended for cleaning and irrigation of the root canal, usually used for the anterior teeth.

STAINLESS STEEL

TITANIUM

COMPATIBLE WITH

REF AE1



ATE1

EMS® 120° FILE HOLDER
MECTRON®

REF AED1

ATED1

SATELEC® NSK® E11



Power



Irrigation



AE2



Mode
E

95° holder used to hold the file (Ø 0,7 or Ø 0,8 mm). Recommended for cleaning and irrigation of the root canal, usually used for the molar teeth.

STAINLESS STEEL

TITANIUM

COMPATIBLE WITH

REF AE2



ATE2

EMS® 90° FILE HOLDER
MECTRON®

REF AED2

ATED2

SATELEC® NSK® E12



Power



Irrigation





AE3



Recommended for melting gutta percha. Immerse the instrument in eugenol before inserting it in the root canal. Slightly press the instrument against the gutta percha point and activate the scaler to condense until the canal is filled.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AE3	ATE3	EMS [®] H MECTRON [®] D3
REF	AED3	ATED3	SATELEC [®] S07 NSK [®] E6
REF	AES3		SIRONA [®]



Mode
E

Power



Irrigation





AE3D



Diamond-coated (40µm) instrument for enlarging the root canal, e.g. to access broken files in the root canal. Do not touch the broken file when enlarging the root canal to avoid pushing it deeper in the canal.

	STAINLESS STEEL	COMPATIBLE WITH
REF	AE3D	EMS [®] MECTRON [®]
REF	AED3D	SATELEC [®] ET18D NSK [®]
REF	AES3D	SIRONA [®]



Mode
E

Power



Irrigation




AE4



Used to remove root filling obstructions during the root canal retreatment. The length of tip slender part is 22 mm.

	STAINLESS STEEL	TITANIUM	COMPATIBLE WITH
REF	AE4	ATE4	EMS [®] MECTRON [®]
REF	AED4	ATED4	SATELEC [®] ET40 NSK [®] E8



Mode
E

Power



Irrigation



AE4D



Diamond-coated (40µm) instrument for removing calcification and bad filling material in the root canal during the root canal retreatment. The length of the diamond-coating on tip is 15 mm.

STAINLESS STEEL




COMPATIBLE WITH

REF AE4D

EMS[®] MECTRON[®]

REF AED4D

SATELEC[®] NSK[®]



Mode
E

Power



Irrigation



AE5

Used to remove root filling obstructions during root canal retreatment. The length of tip slender part is 16 mm.

STAINLESS STEEL

TITANIUM

COMPATIBLE WITH

REF AE5


ATE5

EMS[®] MECTRON[®]

REF AED5

ATED5

SATELEC[®] ET20 NSK[®] E7



Mode
E

Power



Irrigation



AE5D



Diamond-coated (40µm) instrument for removing calcification and bad filling material in the root canal during root canal retreatment. The length of the diamond-coating on tip is 10 mm.

STAINLESS STEEL



COMPATIBLE WITH

REF AE5D

EMS[®] MECTRON[®]

REF AED5D

SATELEC[®] NSK[®]



Mode
E

Power





Irrigation





AUE1



			
TITANIUM		COMPATIBLE WITH	
REF	AUE1	MECTRON® EN1	
			
Mode ROOT	Power	Endo (Maximum power permitted)	Irrigation 



AUE2



TITANIUM		COMPATIBLE WITH	
REF	AUE2	MECTRON® EN2	
			
Mode ROOT	Power	Endo (Maximum power permitted)	Irrigation 

AUE3



			
TITANIUM		COMPATIBLE WITH	
REF	AUE3	MECTRON® EN3	
			
Mode ROOT	Power	Endo (Maximum power permitted)	Irrigation 

AUE4





Used for gentle canal cleaning (soft apical root debridement). The length of tip slender smooth part is 3,5 mm.

TITANIUM

COMPATIBLE WITH

REF AUE4

MECTRON® EN4



Mode
ROOT

Power

Endo (Maximum
power permitted)

Irrigation



ASB1



Diamond-coated (85 µm) instrument for removing caries on occlusal surface of teeth and neck of teeth. Place the instrument on the pit and move slowly into it with light pressure.

STAINLESS STEEL


COMPATIBLE WITH

REF ASB1

EMS® SB MECTRON®

REF ASBD1

SATELEC® EX1 NSK® G62D



Mode
P

Power



Irrigation



ASB2



Diamond-coated (85 µm) instrument for removing caries on mesial surface of the adjacent teeth.

STAINLESS STEEL



COMPATIBLE WITH

REF ASB2

EMS® SBM MECTRON®

REF ASBD2

SATELEC® EX2 NSK® G65D



Mode
P

Power



Irrigation



ASB3



Diamond-coated (85 µm) instrument for removing caries on distal surface of the adjacent teeth.

STAINLESS STEEL



COMPATIBLE WITH

REF ASB3

EMS® SBD MECTRON®

REF ASBD3

SATELEC® EX3 NSK® G66D



Mode
P

Power



Irrigation



ASBL



Set 45° to the left, diamond-coated (85 µm) instrument to remove caries without damaging the adjacent teeth. Half ball diamond tip.

STAINLESS STEEL



COMPATIBLE WITH

REF ASBL

EMS® MECTRON®

REF ASBDL

SATELEC® EXL NSK®



Mode
E

Power



Irrigation



ASBR



Set 45° to the right, diamond-coated (85 µm) instrument to remove caries without damaging the adjacent teeth. Half ball diamond tip.

STAINLESS STEEL



COMPATIBLE WITH

REF ASBR

EMS® MECTRON®

REF ASBDR

SATELEC® EXR NSK®



Mode
P

Power



Irrigation



AUS1





High effectiveness osteotomy of large bone sections during maxillofacial surgery. Can also be used for exodontia/exodontics.

TITANIUM

COMPATIBLE WITH

REF AUS1

MECTRON® OT6



Mode
BONE

Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AUS1R





Right-angled instrument (90°) for horizontal osteotomy technique in maxilla and mandible.

TITANIUM

COMPATIBLE WITH

REF AUS1R

MECTRON® OT8R



Mode
BONE

Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AUS1L





Left-angled instrument (90°) for horizontal osteotomy technique in maxilla and mandible.

TITANIUM

COMPATIBLE WITH

REF AUS1L

MECTRON® OT8L



Mode
BONE

Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AUS2





High effectiveness osteotomy technique in maxilla and mandible (e.g. ridge expansion, corticotomy technique, bone block grafting). With deep marking 7 - 8,5 - 10 mm on the surface.

TITANIUM

COMPATIBLE WITH

REF AUS2

MECTRON® OT7



Mode
BONE

Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AUS3





Micrometric osteotomy: osteotome of great precision in anatomically thin structures.

TITANIUM

COMPATIBLE WITH

REF AUS3

MECTRON® OT2



Mode
BONE

Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AUS4




Universal osteoplasty: periodontal ostectomy, bone chips harvesting, inflammatory tissue removal (e.g. cysts).

TITANIUM

COMPATIBLE WITH

REF AUS4

MECTRON® OP3



Mode
BONE

Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AUS5




High efficiency bone osteoplasty: bone remodelling and harvesting of bone chips.

TITANIUM

COMPATIBLE WITH

REF AUS5

MECTRON® OP1



Mode
BONE

Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AUS6




Osteotomy: osteotome of great precision in anatomically thin structures (ridge expansion, corticotomy technique, Le-Fort-Osteotomies, root splitting during extraction)

TITANIUM

COMPATIBLE WITH

REF AUS6

MECTRON® OT7S-3



Mode
BONE

Power

Special

Irrigation



AUL1




Tip angle 120°, separation of the sinus membrane in internal zones, non-cutting elevator of the sinus membrane.

TITANIUM

COMPATIBLE WITH

REF AUL1

MECTRON® EL3



Mode
ROOT

Power

Endo (Maximum
power permitted)

Irrigation



AUL2




Schneiderian membrane separation from bony walls: separation of the sinus membrane.

TITANIUM

COMPATIBLE WITH

REF AUL2

MECTRON® EL1

Mode
ROOT

Power

Endo (Maximum
power permitted)

Irrigation



AUL3





Diamond-coated (100µm) instrument for micrometric osteotomy or osteoplasty: non-traumatic, to finalise the osteotomy or osteoplasty on thin bone and/or near delicate anatomic structures (e.g. sinus membrane, blood vessels, nerves).

TITANIUM

COMPATIBLE WITH

REF AUL3

MECTRON® OT5

Mode
BONE

Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AUL4





Diamond-coated (100µm) instrument for micrometric osteotomy: to finalise the osteotomy in proximity of soft tissue (e.g. sinus membrane, blood vessels, alveolar nerve).

TITANIUM

COMPATIBLE WITH

REF AUL4

MECTRON® OT1

Mode
BONE

Power

Quality1, Quality2, Quality3
(Cortical/Spongious)

Irrigation



AUL5





Tip angle 95°, separation of the sinus membrane in internal zones, non-cutting elevator of the sinus membrane.

TITANIUM

COMPATIBLE WITH

REF AUL5

MECTRON® EL2



Mode
ROOT

Power

Endo (Maximum
power permitted)

Irrigation



AUC1




Used to cut off the ankylosis and root fraction techniques.

TITANIUM

COMPATIBLE WITH

REF AUC1

MECTRON® EX1



Mode
BONE

Power

Quality1,Quality2,Quality3
(Cortical/Spongious)

Irrigation



AUI1





Ø 1,6 mm implantation site preparation insert. The working length of the tip is 9 mm.

TITANIUM

COMPATIBLE WITH

REF AUI1

MECTRON®



Mode
BONE

Power

Quality1,Quality2,Quality3
(Cortical/Spongious)

Irrigation



AUI2





Ø 2 mm implantation site preparation insert. Irrigation from tip point to prevent overheating.

TITANIUM

COMPATIBLE WITH

REF AUI2

MECTRON®



Mode
BONE

Power

Quality1,Quality2,Quality3
(Cortical/Spongious)

Irrigation



AUI7





Diamond-coated (85 µm) instrument for finalizing the implantation site preparation close to the alveolar nerve.

TITANIUM

COMPATIBLE WITH

REF AUI7

MECTRON® OT4



Mode
BONE

Power

Quality1,Quality2,Quality3
(Cortical/Spongious)

Irrigation



ZEG-Torque Wrench

- secure fastening of scaler-tips
- can also be used to pick up instruments
- packaging: 1 piece



ATW02S
ATW02E



ATU-1



ATW01S
ATW01E



ATK-4I
ATS-4L



ATW-A
ATD-A



REF	compatible with	packaging
ATW02S	Satelec-Tips	1 piece
ATW02E	EMS-Tips	1 piece
ATU-1	surgery	1 piece

REF	compatible with	packaging
ATW01S	Satelec-Tips	1 piece
ATW01E	EMS-Tips	1 piece
ATK-4I	Kavo	1 piece
ATS-4L	Sirona	1 piece
ATW-A	EMS-Tips	1 piece
ATD-A	Satelec-Tips	1 piece

A lginates	6-7	F ace Masks	12, 14	Patient Towels	16
Amalgam Gun	63	Facial Tissues	20	Periodontic Instruments	104-105
Applicator Brushes/Holder	60-62	Files	47-48	Periodontic Tips	113-115
Applicator Tips	62	Filter	74	Pipettes, Disposable	64
Apron	18	Foam Inserts	46	Pocket Probes	93
Apron Dispenser	19	Foam Pellets	20	Prophy Angles	35
Aspirator Tips, Cleaning Brush	24	G as Lighter Refill	89	Prophy Brushes	34
Aspirator Tips, Disposable	24	Gloves	23	Prophy Cups	34
Aspirator Tips, Surgical	93	H and Disinfectant	70	Prophylactic Paste	35
Attachment Tray	86	Handpiece	77	R apid Disinfectant	70
B ib	16	Headrest Cover	15	Retainer Boxes	29-30
Bib Clips	17	Helix Test	74	Rinse Beakers, Disposable	13
Bite Registration	9-10	I ce Spray	45	Rinse Beakers	30
Bone Surgery Tips	124-126	Implantation Tips	128-129	Root Canal Sealer/Filling Materials	50
Brushing Timer	30, 43	Impression Disinfectant	72	S aliva Ejectors	24
Bur Bath	72	Injection Needles	23	Scaling Tips	108-113
Bur Cleaning Brush	96	Instrument Bath	74	Scalpels	94
Bur Stand	53	Instrument Carts and Trolleys	82-86	Sharps Container	101
Bur Holder	96	Instrument Disinfectant	72	Silicone Mixing Bowls	64
C apsule Mixers/Accessories	81	Instrument Tray Inserts	101	Silicones	8-10
Cassettes	96-99	Interdental Brush Holder	27	Sinus-Lifting Tips	126-128
Cavity Preparation Tips	122-123	Interdental Brushes	25-27	Spatulas, Disposable	63
Cellulose Pads	20	Intraoral Tips	10	Stain Remover	33
Cheek Retractors	95	Irrigation Tubing Set	75	Stand Burner	89
Colour Sets	22	L ED-Lights/-Bulbs	79	Sterilisation Reel/Pouch	76
Composite Brushes/Holder	61-62	Lining Material	65	Sterilisation Wrap, Non-Woven	76
Composites	65-67	M agnifying Lamps	78	Steri-Tape	76
Cooling Compress	45	Medication Dispenser	88	Surface Disinfectant	72
Cotton Rolls/Cotton Roll Dispenser	20-21	Meta	65-68	T ooth Whitening Strips	31
Couplings	79	Micro Applicators	60-61	Tongue Cleaner	33
Cup Dispenser	13	Milk Tooth Boxes	39	Tooth Jewellery	54-58
Cups	13	Mixing Bowl	8	Toothbrush Cup	30
Curettes	104-107	Mixing Gun	10	Toothbrush Cup, Kids	41
Curing Lights	77-78	Mixing Pad	9	Toothbrush, Kids, Electronic	40
D appen Dish, Disposable	63	Mixing Spatula	8	Toothbrushes	30-33
Dappen Dish, Glass	64	Mixing Tips	11	Toothbrushes, Kids/Baby	39-41
Demonstration Model	33	Mixing Tray, Disposable	63	Toothbrushes, Disposable	33
Dental Floss	28	Motivational Aids for Kids	38-43	Torch	89
Dental Pressure Lock Bags	89	Mouth Mirrors, Mirror Handles	92, 103	Tray Cover/Paper	16
Dentist Chairs, Saddle Stool	86-87	N ail Brush	74	Tray Material	7
Device for Root Canal Filling	80	Napkins	17	Trays	99-100
Disinfection Wipes	69-71	Needle Remover	101	Tweezer	93
Disinfection Wipes Empty Dispenser	71	O il Spray/Adapters	75	W ash Lotion	73
Dispenser, Refillable	69-71	Orange Clean	71	Water Distiller/Accessories	80
Displays	36-37	Oshibori	44	Wax Knife	89
Disposable Dispensing Syringe	10	P aper Towel Dispenser	19	wellsakids	40-43
Dosing Aids	74	Paper Towels	19	Wipes	70
E ndo Boxes	51-53	Paper, Gutta Percha Points	48-49	Work Pan	89
Endo-Cleaner	53	Patient Aprons/Bibs	18	Z EG-Torque Wrench	129
Endodontic Instruments	47-49	Patient Brochures	90		
Endodontic Tips	119-122				
Etching Gel	59				
Exelcymosis Tip	128				
Explorers	93				

**Germany**

Brahestraße 16
D-04347 Leipzig
sales@wellsamed.de

Phone +49 341 27139 142
Fax +49 341 27139 140

Switzerland

Betpurstrasse 14
CH-8910 Affoltern am Albis

Phone +41 44 419 20 20
Fax +41 44 419 20 19

Bulgaria

Milin kamak Str. 50A
BG-1421 Sofia

Phone +359 2 936 06 44

www.wellsamed.de

SUPERSOFT 3



GERMAN PREMIUM ORAL CARE

wellsamed SUPERSOFT 3

- with over 2,000 soft Tynex®-filaments
- particularly gentle on gums
- compact brush head
- plane bristles
- gummed and ergonomic handle

compact
head

0,15 mm
tynex
filaments

ergonomic
handle



NEW

